MiniLoc[®]

Safety Infusion Set

R_x Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.

Bard and MiniLoc are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.

©2008 C. R. Bard, Inc. All rights reserved.

Bard Access Systems, Inc.
Salt Lake City, UT 84116
www.bardaccess.com

BARD

EC REP

Bard Limited
Forest House, Brighton Road
Crawley, West Sussex
RH11 9BP UK

CE
0344

STERILE EO	Sterile	Use By	Priming Capacity	Patent Information
Single Use	LOT Lot Number	Needleless Y-Injection Site	Registered Trademark	
Attention, See Instructions For Use	Contents	Manufactured By	Trademark	
REF Reorder Number	PVC Tubing Does Not Contain DEHP	Y-Site	Latex Free	

Indications for Use:

This product is a safety intravascular administration set and is indicated for the administration of fluids and drugs, or blood sampling through surgically implanted vascular ports. It is designed to help protect against exposure to bloodborne pathogens caused by accidental needlestick injuries. This product utilizes a non-coring right angle Huber needle and is available with and without an adaptable Y-site which will accommodate the addition of a needleless access or standard injection port. The product is a sterile, non-pyrogenic, single use device. It is recommended that this product be changed in accordance with U.S. Centers for Disease Control (CDC) guidelines for administration sets, local or country specific guidelines, professional standards of practice, and/or according to your institution's policy for Huber needle IV administration sets. Do not use this product if it appears damaged or if the package has been previously opened or damaged.

Warnings:

- Verify needle length is correct based on port reservoir depth, tissue thickness, and the thickness of any dressing beneath the bend of the needle; if too long, needle and/or portal may be damaged at insertion; if too short, needle may not completely pierce portal septum, and medication may be delivered into surrounding tissue and/or needle may be blocked.
- Fully tighten all connections, Y-site end caps, or needleless connectors before use. Failure to attach an end cap or appropriate needleless device after removing a male luer locking end cap or needleless connector can result in an embolism or bleeding.
- Needleless Y-injection site valve may leak slightly at certain pressures resulting in a small bead of fluid on the valve.
- Failure to use the device correctly (securely hold finger tabs down; firmly pull wings up until you hear or feel a "click" and visually observe the orange dot 6) when removing needle from port site could result in the needle tip re-emerging from the base, resulting in an accidental needlestick with a contaminated needle. A needlestick with a contaminated needle may cause infectious disease.
- Handle and discard in accordance with accepted medical practice and all applicable regulations. After use, this product is a potential biohazard.
- This is a single-use product. **DO NOT REUSE. DO NOT RESTERILIZE.**
- Do not alter this device.


Cautions:

- Care must be taken to avoid accidental needlesticks.
- Only qualified healthcare practitioners should insert, manipulate and remove these devices.
- Carefully read and follow all instructions prior to use.
- Follow all instructions, contraindications, warnings, cautions and precautions for all infusates, ports, IV pumps, IV sets and needleless systems, as specified by the respective manufacturer.
- Follow standard infection control precautions as specified by the CDC (USA) or local equivalent.
- High pressure or use with power injectors may cause leakage or damage. It is necessary to know the maximum pressure rating for the implanted port which you are accessing.
- Confirm correct needle placement in the port reservoir by aspiration of blood before infusion of any substance. If there is doubt regarding proper needle placement, perform a radiographic dye procedure to confirm placement per institutional protocol.
- Do not remove and reinsert the needle into the port.
- Avoid excessive manipulation once the needle is in the port.
- Do not access the needleless Y-injection site valve with a needle. Valve puncture may result in an embolism.

Instructions For Use:

Use aseptic technique. Refer to the numbered illustrations when using this device.

- Prepare port site for sterile needle insertion. **1**
- Attach a 10 ml syringe containing normal saline to the proximal Luer lock connection of the safety Huber needle infusion set. Grasp the colored wings and remove the needle cover.
- Prime and flush the infusion set. **2** Insert needle perpendicular into the port septum. **3** Verify patency. Dress and secure site per institutional protocol. **4** Begin injection, infusion, or blood aspiration per institutional procedure.
- PORT DEACCESS:** After therapy completion, flush port per institutional protocol. Close the clamp while injecting the last 0.5 ml of solution. Stabilize the port by **Securely Holding** the finger tabs down; **5A Firmly Pull** the wings up until you hear or feel a "click" **5B** and visually observe the orange dot. **6**
- Dispose of set in a sharps container.

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