



## STRETCHING GUIDE AND INDICATOR



## 1 HOW TO REMOVE PAPER BACKING - END




**Tear**  
Twist and pull to tear paper 2" from the end.




**Peel**  
Remove the 2" of paper from end.

## HOW TO REMOVE PAPER BACKING - MIDDLE




**Tear**  
Twist and pull to tear paper in the center.



**Peel**  
Peel paper leaving 2" of paper on both sides for handling.

## 2 HOW TO APPLY TAPE - END




**Apply**  
Apply as shown.

Apply as needed. Can be worn up to 8 days.

- ANCHOR first 2 inches without stretch.
- Apply stretch according to application
- ANCHOR last 2 inches without stretch.
- After applying, rub tape for best adhesion.

## HOW TO APPLY TAPE - MIDDLE



**Apply**  
Apply as shown.

Read all cautions on box.

- Apply stretch according to application
- ANCHOR first 2 inches without stretch.
- ANCHOR last 2 inches without stretch.
- After applying, rub tape for best adhesion.

## 3 HOW TO REMOVE TAPE



Begin by peeling tape up at a corner. Hold skin down where tape is being pulled off with one hand while gently pulling tape off with the other hand. Remove slowly.

See instructions for strong adhesive tape removal on our website

## WATCH OUR STEP-BY-STEP VIDEOS



[KT TAPE.COM/INSTRUCTIONS](http://KT TAPE.COM/INSTRUCTIONS)

[KT TAPE.COM](http://KT TAPE.COM)

## TIPS FOR BEST RESULTS



### APPLY BEFORE ACTIVITY

Apply tape one hour before beginning activity for best adhesion.



### TRIM HAIR

Trim excessive hair for better adhesion.



### CLEAN SKIN

Clean dirt, oils and lotions from area with rubbing alcohol.



### DON'T RIP OFF!

Gently remove tape from skin, use baby oil if necessary and remove slowly.



### DON'T OVERSTRETCH

Overstretching tape may result in skin irritation and/or abrasions.



### END ON SKIN, NOT ON TAPE

Apply directly to skin when possible, especially the ends.

## HOW TO CHOOSE KT TAPE FOR EVERY MOMENT

KEY BENEFIT	PRODUCT NAME	LIGHT ACTIVITY	MODERATE ACTIVITY	INTENSE ACTIVITY
EASY REMOVAL	GENTLE	█	█	█
ALL-NATURAL COTTON	COTTON	█	█	█
BREATHABILITY	PRO®	█	█	█

KEY BENEFIT	PRODUCT NAME	LIGHT ACTIVITY	MODERATE ACTIVITY	INTENSE ACTIVITY
LARGE AREAS	PRO® WIDE	█	█	█
FOR HUMID CONDITIONS	PRO® EXTREME®	█	█	█
INCREASED CELLULAR OXYGENATION	PRO® OXYGEN	█	█	█



## PREPARE

Use KT products while training for your long-term goals, warming up for workouts, and for day-to-day support to keep your muscles, tendons, and ligaments ready to go.

## PERFORM

Quality products from KT provide the flexible, lightweight support your body needs to perform without holding you back.

## RECOVER

Recovery is a priority in any fitness routine. KT products amplify your body's recovery process so you can focus on the next session.

**CAUTION:** Improper application, overstretching, or applying on areas of sensitive skin, and improper or quick removal may result in skin irritation, blisters or skin abrasion. In rare cases, a user may experience a burning or stinging sensation which may indicate a mild allergic skin reaction. Do not place KT TAPE® directly over cuts or open wounds. One or more of the components of this product is a dyed material. There is a possibility that some color migration may occur during use. Instructions provided are for educational use only. KT TAPE® is not a replacement for professional medical care. Warranties and remedies limited to product replacement cost.

## VIEW OUR OTHER PRODUCTS AT KTTAPE.COM

PREPARE	PERFORM	RECOVER	PRODUCT NAME
✓	✓	✓	PRO®
	✓	✓	BLISTER PREVENTION TAPE
		✓	PAIN RELIEF GEL

SCAN FOR OTHER LANGUAGES



NEW USER GUIDE

# KT TAPE™