How To Use

step2bed

4 STEPS TO GET IN BED



Approach step2bed, turn sideways towards bed and grab tall grab bar.



With hands on tall bar, step onto step and turn backside to bed.



With hands on tall grab bar sit down.



Lay down to resting position and lift legs over small rail into bed.

5 STEPS TO GET OUT OF BED



Reach closest arm over to grab tall grab bar.



Use rail to lift shoulders off the bed while swinging legs over short rail place feet on step.



Take 15-30 seconds while sitting up to gain your composure.



(When ready) Use one hand on tall grab rail or two hands (one on tall rail, one on short rail) and stand up.



While keeping one hand on tall grab bar slowly step off step2bed onto floor.