

Cosmopor® Sterile

General Product Description/Intended Purpose

Cosmopor® Sterile is a non-active, single-use, self-adhesive, sterile device consisting of an adhesive coated backing material (carrier material) and a woundpad, protected with a silicone paper.

The main constituents of this products is:

- The backing / carrier material: a white non-woven.
- The medical grade adhesive: a transparent synthetic acrylic adhesive.
- The wound pad has two different layers: a wound contact layer made of a polyethylene net to prevent the wound pad to stick to the wound and an absorbent layer made of viscose for the absorption of wound exudate.
- Silicone paper: protective paper in contact with the adhesive side of the device.


The product does not contain resorbable components, medicinal products / drug substances (derived from human blood or plasma or not), tissues or cells of human or animal origin, radioactive substances, nanomaterials or other hazardous components.

The dressings are sealed individually in peel-pouches and packed in cardboard folding boxes.

Application/Indication

Cosmopor® Sterile is a non-active, non-invasive, single-use, sterile self-adhesive dressing consisting of an adhesive coated backing material (carrier material) and a woundpad, protected with a silicone paper, for the protection of postsurgical / primary wounds and absorption of exudate (wounds with a normal exudate production).

The product is intended for humans, without restrictions (in age group, weight range, health status or condition), for use by healthcare professionals and lay-users, in contact with intact skin (healthy skin) and breached skin, and short-term (60 minutes ≤ Duration ≤ 30 days, acc. to MDR) / prolonged use (24 hours < Duration ≤ 30 days, acc. to ISO 10993-1:2018) duration of use.



Cosmopor® Sterile

Reference Numbers

Packaging	Dimensions	Reference
P5	7,2 x 5 cm	900 700
	10 x 6 cm	900 701
	10 x 8 cm	900 707
	15 x 8 cm	900 703
	20 x 10 cm	900 709
P10	7,2 x 5 cm	900 833
P25	10 x 6 cm	900 802
	10 x 8 cm	900 806
	10 x 10 cm	900 820
	15 x 6 cm	900 804
	15 x 8 cm	900 808
	15 x 15 cm	900 823
	20 x 8 cm	900 810
	20 x 10 cm	900 812
	25 x 10 cm	900 814
P50	35 x 10 cm	900 816
	7,2 x 5 cm	900 800

Contraindications

Not indicated to be applied on non-treated infected wounds, ulcers or chronic wounds.

Sterile Device

Prevention symbols added:

Do not use if package is damaged.



Do not resterilize.



Single Use Device

Reusing a single-use medical device is dangerous. Reprocessing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request.

Cosmopor® Sterile

Product Disposal

To minimize the risk of potential infection hazards, or environmental pollution, disposable components of Cosmopor® Sterile should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards.

Incident Reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

Product Performance Characteristics

Cosmopor® Sterile			
Property	Test method	Unit	Specification *
Peel adhesion – Non-woven with acrylic coating adhesion.	Evaluation of the peel test based on DIN ISO, 6133:2004-05, procedure B.	N/25mm	3 ± 1
Absorption capacity	Method acc. to (DIN) EN 13726-1:2002:06.	g/100cm ²	>17

(*) Specification obtained through internal test results. The tests are performed periodically to ensure that the product meets the specifications

Labelling

Lot-No. with 9-Digit Code

	Eg.: 0	XXX	Eg.: 03	XXX
	Year	For internal purposes only	Manufacturing week	For internal purposes only

Manufacturing Date

E.g.: 	2025	01	14
	Year	Month	Day

Use-by-Date

	Eg.: 2025	Eg.: 01	01
	Year	Month	Day

Shelf Life: 5 years

Cosmopor® Sterile

Medical Device



Unique Device Identification (UDI)



Single Sterile Barrier System



Single sterile barrier system

Legal Manufacturer



PAUL HARTMANN AG
Paul-Hartmann-Straße 12
89522 HEIDENHEIM, GERMANY
www.hartmann.info

Latest Date of Revision: 2020.03.25