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Adult

Mölnlycke® Z-Flo™ Fluidized positioner

Mold the positioner to desired shape

The positioner maintains its shape

Please read all of the instructions before use. Use as directed by appropriate clinician. Always follow your facility protocol for turning and positioning.

Mölnlycke® Z-Flo™ Fluidized positioner Adult

Positioner product for Adult and Adolescent hospital population.

Product description

The positioner consists of a fluidized medium composed of viscous fluid, microspheres and polyethylene beads which is encapsulated in a polyurethane film exterior. The positioner can be molded to support multiple anatomical shapes and clinical needs.

Indications for use

Mölnlycke® Z-Flo™ Fluidized positioner is designed to provide patient positioning, pressure redistribution and/or to offload anatomical sites, which may be used as part of a prophylactic therapy to help prevent pressure ulcers.

Indications for use

- Adult and Adolescent
- Patient positioning, turning and repositioning
- Off-loading of heel, sacrum, occiput and other bony prominences
- Contracture management

Lifetime

Product properties will deteriorate over time in use. The rate of deterioration is dependent on use environment and application. The product can be used with one patient until the product has deteriorated and does no longer fulfil its intended use.

Other information:

If the product is used after expiry date product properties cannot be ensured. The product expiry date is printed on the packaging labels.

Warning(s)

- Patients must be properly supervised and monitored in health care facilities when using positioners.
- Do not allow patients to lie face down in positioner unless closely monitored due to risk of suffocation.
- To avoid puncturing, be careful when handling sharp objects in proximity of the product.
- When handling the product, especially during molding, do not use fingernails and avoid wearing rings with raised or sharp edges.

Precaution(s)

- Replace promptly if punctured or damaged.
- If punctured, the positioner should be placed in a sealed plastic bag to avoid spillage of the content, which may increase risk of slips and falls.
- Do not autoclave or microwave.
- Do not X-ray through the positioner

Patent numbers: US 8.171.585 B2, US 9 120 666 B2
Toll free number: USA 1-800-882-4582, Canada 1-800-494-5134
Australian Sponsor Address: Mölnlycke Health Care Pty. Ltd., 12 Narabang Way, Belrose, NSW 2085, Australia
UK Responsible Person: Mölnlycke Health Care Ltd., Medlock Street, Oldham, Lancashire, OL1 3HS, UK

Mölnlycke Health Care AG
Brandstrasse 24, 8952 Schlieren

Manufacturer / Fabricant

Mölnlycke Health Care AB
Entreprenörstråket 21, SE-431 53 Mölndal, Sweden

www.molnlycke.com

The Mölnlycke trademark, name and logo are registered globally to one or more of the Mölnlycke Health Care group of companies

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