

COMFORT COVERUP™

COMPONENTS OF COMFORT COVERUP™: Nylon, Spandex.

COMPONENTS OF HOOK: Polyethylene.

COMPONENTS OF LOOP: Nylon, Urethane Foam.

THIS PRODUCT DOES NOT CONTAIN LATEX.

The Comfort CoverUp™ can be worn over any compression garment. Its sleek design prevents your compression garments from irritating other parts of your body or snagging on surrounding fabrics. The Comfort CoverUp™ also fits well under clothing.

The Comfort CoverUp™ does not provide therapeutic compression. CircAid Medical Products, Inc. assumes no liability for injury or damages resulting from contraindicated use of its products. CircAid is a registered trademark and Comfort CoverUp is a trademark of CircAid Medical Products, Inc. All Rights Reserved. ©2012 CircAid Medical Products, Inc.

 CircAid Medical Products, Inc.
9323 Chesapeake Drive, Suite B2
San Diego, CA 92123 U.S.A.
01-858-576-3550

EC REP

MDSS
Schiffgraben 41
30175 Hannover, Germany


by medi

C-CU-IN-111512-01

SIZING CHART

Note: Take measurements while wearing your compression garment to ensure the best fit

| | SMALL (cm) | LARGE (cm) |
|--|-------------|--------------|
| ARMSLEEVE COMFORT COVERUP™ | | |
| Largest Arm Circumference w/garment | 26cm - 39cm | 40cm - 60cm |
| LOWER LEG COMFORT COVERUP™ | | |
| Largest Calf Circumference w/garment | 34cm -50cm | 51cm - 80cm |
| WHOLE LEG COMFORT COVERUP™ | | |
| Largest Thigh Circumference w/garment | 50cm - 71cm | 72cm - 108cm |
| A custom Comfort CoverUp™ is available for those who do not fit in the above size ranges. Contact CircAid® Customer Service at (800) 247-2243 or info@circaid.com for assistance. | | |

DIRECTIONS FOR USE:

DONNING: Remove the small piece of loop secured between the hook tabs on the Comfort CoverUp™. Lay the Comfort CoverUp™ out flat and fold over the top edge so that the two hook tabs are visible. This will prevent the hook tabs from catching on the garment while sliding it on. Slide the Comfort CoverUp™ over the compression garment and secure the hook tabs on the Comfort CoverUp™ to the garment once it is in its desired position.

REMOVAL: To remove the Comfort CoverUp™, disengage the two hook tabs from the garment, fold over the top edge so the hook tabs are visible, and slide it off your arm or leg before removing your compression garment.

CARE INSTRUCTIONS:

 MACHINE WASH WARM

 TUMBLE DRY LOW

 DO NOT BLEACH

 DO NOT IRON

The Comfort CoverUp™ can be machine washed using a gentle cycle and dried on a low heat setting. It is recommended that the Comfort CoverUp™ be washed inside of a pillowcase to protect other items from the hook tabs. To extend the life of your Comfort CoverUp™ hand wash and drip or flat dry.