



# Max-Fast™ Forehead Sensor

When Timing is Critical



## Challenging Conditions

Some patients represent a challenge for monitoring due to:

- Intense vasoconstriction
- Hypovolemia
- Hypothermia
- Therapeutic hypothermia
- Low cardiac index
- Septic shock
- Severe peripheral vascular diseases


For these patients, it is even more important to have:

- Accurate data correlated to blood gas
- Rapid detection of SpO<sub>2</sub> changes



**ACCURACY, QUICKLY**

When timing is critical, you can be alerted of changes in SpO<sub>2</sub> about two minutes earlier than with digit sensors.<sup>1</sup>




Because forehead circulation is fed by the supraorbital artery, this area is not as prone to vasoconstriction during low perfusion as other sites. Therefore the pulse amplitude is not affected.

## SpO<sub>2</sub> OBVIOUS

Arterial blood travelling from the heart reaches the head sooner than it reaches distal sites such as fingers, especially when patients have poor pulse perfusion.


### IMPACT OF COLD-INDUCED VASOCONSTRICTION



■ Warm room 74 °F= 23°C  
■ Cold room 58 °F= 14°C



The head and fingers are warm and well perfused upon initial exposure to the cold room.




After 45 minutes, the fingers, ears and nose are cold, indicating vasoconstriction and low regional perfusion. The forehead temperature is well maintained.

## Succeed In Challenging Conditions

- When digit sensors fail to detect pulsatile signal, Max-Fast™ forehead sensors can often obtain SpO<sub>2</sub> readings<sup>1</sup>
- During poor peripheral perfusion, the Max-Fast™ forehead sensor typically detects changes in SpO<sub>2</sub> about two minutes earlier than conventional sensors<sup>3</sup>
- During conditions of poor peripheral circulation, the forehead site is less prone to peripheral vasoconstriction
- Forehead sensor is less prone to interference from natural movement
- For patients subject to critical desaturations, the unique LoSat™ expanded accuracy range offers clinical information validated down to 60% SpO<sub>2</sub>.<sup>4</sup>

## Know Sooner To Act Faster

The Max-Fast™ forehead sensor detects changes in SpO<sub>2</sub> faster than finger sensors, and with an accuracy that correlates closer to arterial blood data.<sup>3</sup>



*A technology overview of the Nellcor™ OxiMax™ pulse oximetry system. Internal testing.<sup>3</sup>*

# Specialty Sensor With Special Design

Max-Fast™ forehead sensors have been designed specifically for the forehead:

- Dedicated LED for reflectance technology
- Specific emitter-receptor distance
- Black contact surface to avoid light interference




A single Max-Fast™ forehead sensor can be used up to two days with appropriate site inspection and changes, thanks to its four layers of adhesive.




For accurate and reliable reading, use the soft, adjustable headband packaged with the Max-Fast™ forehead sensor to help prevent venous pulsation at the sensor site and maintain proper sensor position.

“Despite providing a seemingly satisfactory waveform and stable oxygen saturation measures, the use of a disposable (transmittal) finger oxygen saturation sensor on the forehead provided inaccurate readings in 11 out of 20 hypoxic ED patients.”<sup>5</sup>




Finger sensors are not designed to be placed on the forehead and should never be used on this site.



Studies show that finger sensor usage on the forehead can lead to wrong SpO<sub>2</sub> readings despite a normal pleth curve.<sup>5</sup>

The data obtained gives false reassurance as SpO<sub>2</sub> remains high even when the true value is at 70%.<sup>3</sup>

## Using The Max-Fast™ Forehead Sensor Is Very Easy



# Sensors That Span The Entire Patient Spectrum

The Max-Fast™ forehead sensor is part of Covidien's extensive line of sensors, designed to meet the entire patient spectrum and to provide readings in most challenging perfusion and saturation situations.

## OxiMax™ Single-Patient-Use Sensors



MAX-A-I  
MAX-AL-I



MAX-P-I



MAX-I-I



MAX-N-I



MAX-R-I



SoftCare™  
SC-PR-I



SoftCare™  
SC-NEO-I



SoftCare™  
SC-A-I



D-YS



DS100A



OXI-A/N



OXI-P/I



D-YS  
with D-YSE

## OxiMax™ Reusable Sensors

**YOU CAN BE SURE OF IT**




# Nellcor™ OxiMax™ Technology At Work

Due to the unique operating characteristics and calibration curve inside the sensor, Max-Fast™ forehead sensors operate only with Nellcor™ OxiMax™ technology, which is integrated into multiparameter solutions of more than 70 OEM partners worldwide.

**Nellcor™ OxiMax™ technology is available in most OEM monitoring solutions.**



**ASK IF YOU ALREADY CAN USE MAX-FAST™ FOREHEAD SENSORS WITH YOUR CURRENT EQUIPMENT**



## WHEN TIMING IS CRITICAL

### Max-Fast™ forehead sensors provide,

- **Reliable information**

Max-Fast™ forehead sensors give readings when conventional finger sensors fail<sup>6</sup>

- **Accurate readings**

Max-Fast™ forehead sensors readings are correlated with blood gas analysis and are not affected by vasoconstriction<sup>7-9</sup>

- **Fast response**

Max-Fast™ forehead sensors detect changes in SpO<sub>2</sub> about two minutes earlier than conventional sensors<sup>10</sup>

### USE THE RIGHT SENSOR ON THE RIGHT PATIENT

Use Max-Fast™ SpO<sub>2</sub> forehead sensors to know faster, to secure your reading and to increase your accuracy compared to conventional digit sensors<sup>1, 6-10</sup>

1. Bebout DE, Mannheimer PD, Wun C-C. Site-dependent differences in the time to detect changes in saturation during low perfusion. Crit Care Med. 2001;29(12):A115.  
 2. Bebout DE, Mannheimer PD. Effects of cold-induced peripheral vasoconstriction on pulse amplitude at various pulse oximeter sensor sites. Anesthesiology. 2002;96:A558. [Abstract]  
 3. A technology overview of the Nellcor OxiMax pulse oximetry system. Internal testing.  
 4. FDA 510(k).  
 5. Smithline HA, Rudnitzky N, Macomber S, Blankrn FS. Pulse oximetry using a disposable finger sensor placed on the forehead in hypoxic patients. J Emerg Med. April 28, 2009 Apr 28. [Epub ahead of print]  
 6. Schallom L, Sona C, McSweeney M, Mazuski J. Comparison of forehead and digit oximetry in surgical/trauma patients at risk for decreased peripheral perfusion. Heart Lung. 2007;36(3):188-194.

7. Fernandez M, Burns K, Calhoun B, George S, Martin B, Weaver C. Evaluation of a new pulse oximeter sensor. Am J Crit Care. 2007;16(2):146-152.  
 8. Bebout, DE, Mannheimer PD, Asbagh NA. Detection of hypoxemia during peripheral vasoconstriction at the radial artery and various pulse oximeter sensor sites. Critical Care Med. 2003;31(2):A72 [Abstract].  
 9. Berkenbosch JW, Tobias JD. Comparison of a new forehead reflectance pulse oximeter sensor with a conventional digit sensor in pediatric patients. Respir Care. 2006;51(7):726-731.  
 10. MacLeod DB, Cortinez LI, Keifer JC, et al. The desaturation response time of finger pulse oximeters during mild hypothermia. Anaesthesia. 2005;60(1):65-71.



**IMPORTANT** : Please refer to the package insert for complete instructions, contraindications, warnings and precautions

COVIDIEN, COVIDIEN with logo, Covidien logo and "positive results for life" are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company. © 2009 Covidien. - C-MN-MaxFastSensors/GB - 11/2011

COVIDIEN (UK)  
COMMERCIAL LTD

4500 PARKWAY  
WHITELEY  
FAREHAM  
HAMPSHIRE  
PO15 7NY

+44 (0) 1329 224000 [T]  
+44 (0) 1329 220213 [F]

COVIDIEN IRELAND  
COMMERCIAL LTD

BLOCK G, GROUND FLOOR  
CHERRYWOOD BUSINESS PARK,  
LOUGHLINSTOWN,  
CO. DUBLIN

+353 (0) 1.4073173 [T]  
+353 (0) 1.4073174 [F]

WWW.COVIDIEN.COM