# healthsmart®

## Standard Infrared Ear Thermometer

Item #18-220-000

Please read this guidebook  
completely before operating this unit.

Includes: : 1 Thermometer, 1 Storage Case, 30 Lens Covers,  
2 AAA Batteries, 1 User Manual

## ***Index***

---

Introduction . . . . .	2
Intended for Use . . . . .	3
Safety Precautions . . . . .	3-4
Recommendations . . . . .	5
Body Temperature Chart . . . . .	5
Product Overview / LCD Display . . . . .	6
Temperature Taking Hints . . . . .	7
Before Use . . . . .	7-8
Setting the Date and Time . . . . .	8
Switching Between Fahrenheit and Celsius . . . . .	9
Changing Measurement Modes . . . . .	9
How To Take A Reading - Ear . . . . .	10-11
How To Take A Reading - Object . . . . .	11-12
Memory Recall . . . . .	12
Replacing The Battery . . . . .	13
FCC Statement and EMC Compatibility . . . . .	13-14
Maintenance, Storage and Calibration . . . . .	14
Troubleshooting . . . . .	15-16
Product Specifications . . . . .	17
Normalized Symbols . . . . .	17-18
Important Information Regarding Electromagnetic Compatibility (EMC) . . . . .	18 - 23
Warranty . . . . .	23

## ***Introduction***

---

- Thank you for purchasing the Standard Infrared Ear Thermometer.
- For safe and proper use of this product, be sure to read and fully understand Safety Precautions contained in this instruction manual.
- Keep this Manual with you for a convenient reference.
- The reference body site is core and mode of operation is adjusted mode.
- The user can measure and change the battery under normal circumstances and maintain the device and its accessories according to the user manual.

## ***Intended for Use***

---

The device is intended for the intermittent measurement and monitoring of human body temperature in the ear by consumers in the home. This device is intended for use on people of all ages.

## ***Safety Precautions***

---

Proper technique is critical to obtaining accurate temperatures. Therefore, in order to obtain an accurate measurement, please read this manual carefully before using.

- To avoid inaccurate readings, always attach a new, clean lens cover for each ear temperature reading. Lens covers are for single use only. Only HealthSmart branded lens covers may be used with this thermometer.
- Missing, damaged, perforated or dirty lens covers can lead to inaccurate temperature readings.
- Use of this thermometer is not intended as a substitution for consultation with your physician. Measurement results are for reference only. Contact your physician if you have or suspect any health concerns.
- The infrared sensor must be kept clean, dry and undamaged at all times to ensure accurate measurements. If the probe tip is placed into the ear canal without a probe cover installed, please remove from the ear canal and properly clean the lens.
- Earwax and other debris in the ear canal may cause an inaccurate measurement. Ensure the ear canal is clean when taking a measurement.
- Do not touch or blow on the infrared sensor. A dirty sensor may cause inaccuracy.
- Avoid measuring the ear temperature if the ear is inflamed or infected.
- If the temperature of the storage area differs greatly from that of the measurement area, please wait for the thermometer temperature to equalize to the room temperature for about 30 minutes before use.
- Keep the device out of the reach of children and pets to avoid inhalation or swallowing of small parts. Do not allow children to take their own temperatures unattended. Children may not be able to use the device according to the instructions in this user manual. This is not a toy.

- Do not discard the battery into a fire. The battery is likely to explode.
- Do not store the thermometer in extreme temperatures below -13°F (-25°C) or over 131°F (55°C) or in humidity below 15% RH or over 90% RH. Failing to do so may cause inaccuracy.
- Do not use the thermometer if the main body is damaged (for example the infrared sensor is broken) The continuous use of a damaged unit may result in false readings or injury.
- Do not take the thermometer apart, repair or change any parts on the unit at any time.
- Do not use a mobile phone nearby when the thermometer is working.
- If your thermometer will not be used for a long time, it is advised to remove the battery to prevent possible damage due to chemical leakage. If the battery does leak, remove them carefully. Do not allow bare skin to touch the leaking fluid.
- This product needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided. This unit can be affected by portable and mobile RF communications equipment.
- The materials (e.g. ABS) of expect contact with patient had passed the ISO 10993-5 and ISO 10993-10 standards test, no toxicity, allergy and irritations reaction. However, based on current science and technology, other potential allergic reactions are unknown. If you have allergic reaction to the materials, please stop using immediately and consult your physician.
- Store the thermometer in a dry, clean place. Do not leave the thermometer exposed to any chemical solvent, lint, dust direct sunshine or high temperatures.
- When using, do not touch the battery and the patient simultaneously.
- Do not maintain or service the device while the device is in use.

**CAUTION:** This device should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this machine should be observed to verify normal operation in the configuration in which it will be used.

## ***Recommendations***

---

- Don't use the product for other purposes than what it is intended to be used. Only designed for use on the human body.
- It is not a waterproof device. Clean it with a soft, wet cloth or cotton swab.
- Do not leave the product exposed to any chemical solvent, direct sunshine or high temperature.
- Don't drop or impose any vibration or impact on the product.
- Don't use a mobile phone nearby when the thermometer is working.
- Please act according to the native law in handling and disposing of the batteries.
- Take out the batteries if you are not going to use the unit for a long period of time.

## ***Body Temperature Chart***

---

A person's normal temperature will vary and trends to decrease with age. The following table shows normal temperature ranges by age:


### **Normal temperature ranges by age:**

<b>Age</b>	<b>Temperature (°F)</b>	<b>Temperature (°C)</b>
0 - 2 years	97.5 - 100.4 °F	36.4 - 38.0 °C
3 - 10 years	97.0 - 100.0 °F	36.1 - 37.8 °C
11 - 65 years	96.6 - 99.7°F	35.9 - 37.6 °C
> 65 years	96.4 - 99.5°F	35.8 - 37.5 °C

The range of normal temperature varies from person to person and can be influenced by many factors such as time of day, measure from different sites, level of activity, medications, emotion, etc. We recommend that you practice with the thermometer when you or family members are healthy, so you are familiar with the readings when they are ill.


## Product Overview

---

## LCD Display

---

This package contains: (1) Thermometer, (1) Storage Case, (30) Lens Covers, (2) 1.5v AAA Batteries, (1) User Manual

## Temperature Taking Hints

---

1. A temperature measurement taken from the right ear may differ from measurement taken from the left ear. For this reason, always try to take the temperature in the same ear.
2. External factors may influence ear temperature, particularly when an individual has:
  - been lying on one ear or the other
  - had their ears covered
  - been exposed to very hot or very cold temperatures
  - been recently swimming or bathingIn any of these cases listed above, remove the person from the situation or uncover their head and wait 30 minutes prior to taking a temperature.
3. For people wearing hearing aids or ear plugs, remove the device and wait 30 minutes prior to taking a temperature.
4. Use the untreated ear if ear drops or other medications have been placed in the ear canal.
5. When using the thermometer on infants under age 1, pull the ear up to ensure the probe tip faces the eardrum, Fig. 1.
6. When using the thermometers on individuals over the age of 1, pull the ear back to ensure the probe tip faces the eardrum Fig. 2.

FIG. 1

FIG. 2

## Before Use

---

1. To achieve an accurate measurement, it is important to check the infrared sensor and ensure it is clean before use. To clean the sensor, gently wipe its surface with a cotton swab slightly moistened with alcohol and immediately wipe dry with a clean cotton swab. After cleaning, allow at least 20 minutes drying time before taking temperatures.
2. Earwax in the ear canal may cause an inaccurate measurement, so please ensure the ear canal is clean. If you clean the ear canal, wait 5~10 minutes before measuring.

3. If the temperature of the storage area differs greatly from that of the measuring area, please wait 30 minutes for the thermometer and room temperature to equalize before measuring.
4. Open the battery cover and remove the insulating plastic from the battery compartment. The thermometer should turn on automatically.

**CAUTION:** In Standby Mode, the LCD display will show the environment temperature along with the current date and time.

## Setting The Time And Date

**NOTE:** When using this thermometer for the first time, or after replacing the batteries, you will need to reset the time and date.

1. If installing batteries for the first time, skip to step 2. Otherwise, follow this step if resetting the date/time. With the thermometer OFF, press and hold the POWER button for 3 seconds to enter into the setup mode.
2. First, select between a 12hr or 24hr format. Press the SCAN button to switch between formats. Once the desired format is chosen, press the POWER button to confirm.
3. The hour followed by the minute will begin to flash, Fig 1.
4. Scroll through the numbers by pressing SCAN button. Confirm the desired time by pressing the POWER button.
5. The year followed by the month and day will begin to flash.
6. Repeat step #4 to set the date. After setting the day, the thermometer will automatically turn OFF.


FIG. 1

## Switching Between Fahrenheit and Celsius

---

This thermometer is set to measure in Fahrenheit, but can be switched to measure in Celsius.

### To switch from Fahrenheit to Celsius

1. With the thermometer OFF, press and hold the SCAN button for 3 seconds to enter into the setup mode.
2. Press the SCAN button again to switch between °F to °C.
3. Once the preferred mode is selected, press the POWER button set the mode and turn off the thermometer.

## Changing Measurement Modes


---

**NOTE:** The thermometer is set at the factory to measure in Ear mode. However, you may follow the steps below to switch between Ear and Object mode.


### To Change Measurement Modes:

1. Press the POWER button to turn the thermometer ON.
2. Press and hold the POWER button first, then, while still holding down the POWER button, press and hold the SCAN button until the thermometer beeps and switches modes.
3. LCD will display the following icons to indicate the current measurement mode.

 = Ear Measurement Mode


 = Object Measurement Mode

## How To Take A Reading - Ear

The ear is an ideal measurement site because body temperature is regulated by the hypothalamus, which shares the same blood supply as the tympanic membrane or more commonly called eardrum.

**NOTE:** Make sure the measurement sensor is clean before taking a reading. Attempting to take a measurement anywhere else on the body may result in an inaccurate reading.

1. Firmly place a new clean lens cover on the probe tip. Ensure the probe cover fits snugly onto the probe tip, Fig 1.
2. If the thermometer is OFF, press the POWER button to turn it ON. The screen will run a quick self test. The thermometer is ready to take a measurement.

FIG. 1

3. Insert the probe VERY gently and slowly into the ear canal.
  - Hold the thermometer so the probe faces straight in the direction of the eardrum, Fig 3.
  - It is important to point the probe tip toward the eardrum and is inserted inside the ear canal until slight resistance is felt.
  - Hold the outer ear and gently pull it up and toward the back of the head to straighten the ear canal, Fig 3.
  - Always insert the thermometer into the same ear with the same direction and depth.


FIG. 2


FIG. 3

- Press and hold the SCAN button. The thermometer will take a reading and make a 'beep' sound, indicating the measurement is complete.

**NOTE:** The thermometer will beep twice and display the 'ear' symbol when the thermometer is ready to take another measurement.


- Remove the thermometer from the ear canal. The thermometer will display the result on the backlit LCD screen.

**NOTE:** The thermometer is ready to take another measurement once the EAR temperature symbol appears again on the screen. However, it is recommended to wait at least 20 seconds before taking another measurement to avoid excessive skin cooling which could affect the accuracy.

- Since the lens covers are recommended to only be used 1 time, please discard the used lens cover by squeezing the Lens Cover Ejector.

- After 1 minute, the thermometer will automatically shut OFF.

**NOTE:** This thermometer is equipped with a Fever Indicator showing the symbols below on the LCD.

## How To Take A Reading - Object

**NOTE:** Before taking the measurement of an Object, please ensure the thermometer has been switched to Object Mode - see Changing Measurement Modes. No lens cover is required to measure an object. Eject the lens cover if necessary.

- If the thermometer is OFF, press the POWER button to turn it ON. The screen will run a quick self test then appear as shown in Fig. 1 indicating the thermometer is ready to take a measurement.
- Aim the probe tip of the thermometer at the center of the object you want to measure. Keep the probe tip 0.25" ~ 0.75" (1cm ~ 2cm) away from the object.

FIG. 1

## How To Take A Reading - Object (continued)

---

3. Press and hold the SCAN button. The thermometer will take a reading and make a 'beep' sound, indicating the measurement is complete.

**NOTE:** in Object mode, the temperature shown is the actual surface temperature, which is different from body temperature.

## Memory Recall

---

There are a total of 10 memory readings.

**NOTE:** The memory bank holds memories for both ear and object mode.

The current measurement is always stored as number 1, while the oldest reading is always stored as number 10.

### To recall measurements in memory:

1. If the thermometer is OFF, press the POWER button to turn it ON. The screen will run a quick self test then appear as shown in Fig. 1 indicating the thermometer is ready to take a measurement.
2. Press and hold the POWER button for approximately 3 seconds.
3. The screen will display the memory number (Fig. 2) and then the actual measurement (Fig. 3). Press the POWER button to scroll through additional saved readings.

**NOTE:** All subsequent measurements taken will fill the memory spaces until they reach the number 10 spot. Once all of the spaces are filled, the oldest measurement will automatically be deleted from memory.

FIG. 1


FIG. 2

FIG. 3

## Replacing the Battery


---

When the low battery indicator “” appears on the screen, or the thermometer does not turn ON, it is time to replace the batteries.

### To Replace the batteries:

1. Gently slide the battery cover back.
2. Carefully remove the old batteries and properly discard.
3. Insert 2 new AAA 1.5v batteries. Ensure they are inserted correctly according to the proper polarity,

*Fig 1.*

4. Slide the battery cover back on until it snaps into place.

To protect the environment, please dispose of exhausted batteries according to current federal, state and local regulations. Keep batteries out of the reach of children.

## FCC Statement and EMC Compatibility

---

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

## ***FCC Statement and EMC Compatibility***

---

15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

## ***Maintenance, Storage and Calibration***

---

### **Maintenance:**

1. To ensure accurate measurements, keep the infrared sensor and distance sensor tip clean and free of scratches. Finger prints, ear wax or dirt will affect the accuracy of the thermometer.
2. To clean the infrared sensor gently wipe the sensor with a cotton swab slightly moistened with 75% alcohol and immediately wipe dry with a clean cotton swab (Fig 1). The thermometer can be used after 10 minutes when the alcohol dries.
3. Use a soft dry cloth to clean the body of the thermometer (Fig 2).
4. Never clean the thermometer with an abrasive cleanser, thinner, or submerge the thermometer into water or other liquids.

FIG. 1

FIG. 2

### **Storage:**

1. After each use clean the thermometer.
2. Store the thermometer and lens covers in a clean dry place at room temperature. NEVER leave the thermometer exposed to direct sun, high temperatures, humidity, fire, flames, shocks or impacts.
3. If you do not use the thermometer for a long time, take out the batteries.

### **Calibration:**

1. The thermometer is factory-calibrated. If the thermometer is used according to the user manual, periodic calibration is not required. Do not try to modify or reassemble the thermometer.

## Troubleshooting

Error Message	Problem	Solution
Er 1	A measurement was attempted before the thermometer was ready.	Only attempt to take a measurement when EAR or OBJECT symbols appear on the display.
Er 2	The ambient temperature is not within the proper measurement range between 50°F ~ 104°F (10° ~ 40°C)	Place the thermometer in a room for at least 30 minutes at room temperature between 50°F ~ 104°F (10° ~ 40°C)
Er 3	The thermometer is placed incorrectly or is unsteady.	Refer to "How To Take A Reading - Ear" and take a new temperature reading.
Er 4	The thermometer is detecting rapid ambient temperature change.	Place the thermometer in a room for at least 30 minutes at room temperature between 50°F ~ 104°F (10° ~ 40°C)
Er 5	The thermometer is not functioning properly.	Remove the batteries and wait 1 minute before replacing the batteries.
Er 6	The thermometer is in OBJECT mode, but is detecting a lens cover.	Remove the lens cover and attempt a new measurement.

## Troubleshooting

Error Message	Problem	Solution
	In Ear mode: Temperature taken is higher than 109°F (43.0°C) In Object mode: Temperature taken is higher than 212°F (100°C)	Refer to “How To Take A Reading - Ear” and check the integrity of the lens cover. Or refer to “How To Take A Reading - Object” and take a new temperature reading.
	In Ear mode: Temperature taken is lower than 82.4°F (28.0°C) In Object mode: Temperature taken is lower than -4°F (-20°C)	Refer to “How To Take A Reading - Ear” and check the integrity of the lens cover. Or refer to “How To Take A Reading - Object” and take a new temperature reading.
	The battery is low	Refer to “Replacing the Battery” and insert new AAA alkaline batteries

## Specifications

---

- **Product Name:** Standard Infrared Ear Thermometer
- **Power supply:** DC 3V (2 x AAA battery)
- **Ear Measuring Range:** 82.4°F ~ 109.4°F (28.0°C ~ 43.0°C)  
**Ear Measuring Accuracy:**  $\pm 0.4^{\circ}\text{F}$  (or  $\pm 0.2^{\circ}\text{C}$ )  
from 95.9°F ~ 107.6°F (or 35.5°C ~ 42.0°C)  $\pm 0.5^{\circ}\text{F}$  ( $\pm 0.3^{\circ}\text{C}$ ) for other temperature range
- **Object Measuring Range:** -4°F ~ 212°F (-20°C ~ 100°C)
- **Object Measuring Accuracy:**  $\pm 4^{\circ}\text{F}$  or  $\pm 4^{\circ}\text{F}(\pm 2^{\circ}\text{C})$  which ever is greater
- **Measurement Time:** 1 Second
- **Dimension:** 1.9 (L) x 1.50 (W) x 6.1 (H) in. [47(L) x 38(W) x 155(H) mm]
- **Weight:** 83g (with batteries)
- **Operation Temperature:** 50.0°F ~ 104.0°F (10°C ~ 40°C)
- **Operation Relative Humidity:** 15% RH – 85%
- **Storage & Transportation Temperature:** -13°F ~ 131°F (-25°C ~ 55°C)
- **Storage & Transportation Relative Humidity:** 15% RH ~ 95% RH
- **Service Life:** 3 years
- **Battery Life:** Approximately 1 year / 6,000 readings

## Normalized Symbols

---

Applied part of type BF

Disposal in accordance with Directive  
2012/19/EU (WEEE)

The name and the address of the manufacturer.

Caution

Refer to Instruction Manual.

## ***Normalized Symbols***

---

**IP22** The first number 2: Protected against solid foreign objects of 12,5 mm  $\Phi$  and greater. The second number: Protected against vertically falling water drops when enclosure tilted up to 15°. Vertically falling drops shall have no harmful effects when the enclosure is tilted at any angle up to 15° on either side of the vertical.

## ***Important Information Regarding Electromagnetic Compatibility (EMC)***

---

With the increased number of electronic devices such as PC's and mobile (cellular) telephones, medical devices in use may be susceptible to electromagnetic interference from other devices. Electromagnetic interference may result in incorrect operation of the medical device and create a potentially unsafe situation. Medical devices should also not interfere with other devices.

In order to regulate the requirements for EMC (Electro Magnetic Compatibility) with the aim to prevent unsafe product situations, the EN60601-1-2 standard has been implemented. This standard defines the levels of immunity to electromagnetic interferences as well as maximum levels of electromagnetic emissions for medical devices. Medical devices manufactured by HealthSmart International conform to this EN60601-1-2:2007 standard for both immunity and emissions.

Nevertheless, special precautions need to be observed: The use of accessories and cables other than those specified by manufacturer, with the exception of cables sold by manufacturer as replacement parts for internal components, may result in increased emission or decreased immunity of the device.

- The medical devices should not be used adjacent to or stacked with other equipment. In case adjacent or stacked use is necessary, the medical device should be observed to verify normal operation in the configuration in which it will be used.
- Refer to further guidance below regarding the EMC environment in which the device should be used.

## ***Guidance and Manufacturer's Declaration***

(Table 1)

<b>Guidance and manufacturer's declaration – electromagnetic emissions</b>		
The thermometer are intended for use in the electromagnetic environment specified below. The customer or the user of these thermometer should assure that it is used in such environment.		
<b>Emissions test</b>	<b>Compliance</b>	<b>Electromagnetic environment - guidance</b>
RF emissions CISPR 11	Group 1	The thermometer use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The thermometer are suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC61000-3-3	Not applicable	

Table 2

Guidance and declaration of manufacturer-electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment guidance
Electrostatic discharge (ESD) IEC 61000-4-2	$\pm 6$ kV contact $\pm 8$ kV air	$\pm 6$ kV contact $\pm 8$ kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrostatic transient/burst IEC 61000-4-4	$\pm 2$ kV for power supply lines $\pm 1$ kV for input/output lines	N/A	
Surge IEC 61000-4-5	$\pm 1$ kV differential mode $\pm 2$ kV common mode	N/A	
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycle 70% UT (30% dip in UT) for 25 cycle <5% UT (>95% dip in UT) for 5 sec	N/A	
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Table 3

Guidance and declaration of manufacturer-electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environmentguidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 Mhz	N/A	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 Ghz	3 V/m	<p>Recommended seperation distance</p> $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[ \frac{7}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ Ghz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

Table 4

Recommended separation distances between portable and mobile RF communications equipment and the device
The device is intended for use in an electromagnetic environment in which radiated therefore disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter  W	Separation distance according to frequency of transmitter m	
	80 MHz to 800 MHz $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12	0.23
0.1	0.38	0.73
1	1.2	2.3
10	3.8	7.3
100	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## Warranty

---

Three-year warranty is available from purchasing date, excluding user-caused failures listed below:

- (1) Failure resulted in unauthorized disassembly and modification
- (2) Failure resulted in unexpected drop during application or transportation
- (3) Failure resulted in operation away from proper instruction in user's manual

## Declaration of Conformity

---

EN12470-5, ISO 80601-2-56, EN 60601-1-11, EN 60601-1-2(EMC),  
IEC/EN 60601-1, ISO 13485

### NOTES:

- 1) Under the environment with electrostatic discharge, the unit may malfunction and may require user to reset the unit.
- 2) Dispose with care, please consult the retailer for details.

**hs**  
**Satisfaction  
Guaranteed**

Manufactured for:

**HealthSmart® International**

1931 Norman Drive  
Waukegan, IL 60085

Register this product at [LiveHealthSmart.com](http://LiveHealthSmart.com)

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©2019 HealthSmart® International