medi GmbH & Co. KG
Medicusstraße 1
D-95448 Bayreuth
Germany
T +49 921 912-0
www.medi.de

medi Australia Pty Ltd
83 Fennell Street
North Parramatta NSW 2151
Australia
T +61-2 9890 8696
F +61-2 9890 8439
sales@mediaustralia.com.au
www.mediaustralia.com.au

medi Austria GmbH
Adamgasse 16/7
6020 Innsbruck
Austria
T +43 512 57 95 15
F +43 512 57 95 15 45
vertrieb@medi-austria.at
www.medi-austria.at

medi Belgium bvba
Posthoornstraat 13/1
3582 Koersel
Belgium
T: + 32-11 24 25 60
F: +32-11 24 25 64
info@medibelgium.be
www.medi-belgium.be

medi Brasil
Rua Cristovan de Vita, n.º 260, Un. 30,
Centro Logístico Raposo Tavares,
Vila Camargo
Vargem Grande Paulista - SP
CEP 06730-000
Brasil
T: +55-11-3500 8005
sac@medibrasil.com
www.medi-brasil.com

medi Canada Inc / médi Canada Inc
597, Rue Duvernay, Verchères
QC Canada J0L 2R0
T: +1 450-583-3317 / +1 800 361 3153
F: +1 888-583-6827
service@medicanada.ca
www.medicana.ca

MAXIS a.s.,
medi group company
Slezská 2127/13
120 00 Prague 2
Czech Republic
T: +420 571 633 510
F: +420 571 616 271
info@maxis-medica.com
www.maxis-medica.com

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medi Trading (Shanghai) Co. Ltd.
Room 105, Building No. 4,
Lane 500 ShengXia Road
Pudong New District
201210 Shanghai
People's Republic of China
T: +86-21 61761988
F: +86-21 61769607
infocn@medi.cn

Medi Turk Ortopedi Medikal İc ve Dis Tic. AS
Mustafa Kemal Mah. 2156 Sokak No: 16 /4-5
06510 Çankaya Ankara
Turkey
T: +90 312 435 20 26
F: +90 312 434 17 67
info@medi-turk.com
www.medi-turk.com

medi UK Ltd.
Plough Lane
Hereford HR4 OEL
Great Britain
T +44-1432 37 35 00
F +44-1432 37 35 10
enquiries@mediuk.co.uk
www.mediuk.co.uk

medi

circaid® reduction kit™ lower leg inelastic compression system

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


circaid® reduction kit lower leg

Purpose

The compression system is designed to provide compression to the leg for patients with venous and lymphatic disorders.

The device consists of a series of juxtapositioned, inelastic bands extending from the central portion of the garment. The wrap contours the leg and is designed to cover from just above the ankle bone to just below the knee.

Indications

- Lymphedema
- other forms of edema
- Lipedema
- CVI (C3-C6)
- Varicose veins (C3)
- DVT/Thrombosis Prevention
- Post Thrombotic Syndrome
- Venous Stasis Ulcer (C6)
- Post Sclerotherapy
- Lipodermatosclerosis (C4b)

Contraindication

- Contraindications
- Severe Peripheral Arterial Disease
- Uncontrolled Congestive Heart Failure
- Septic Phlebitis
- Phlegmasia Cerulea Dolens
- Untreated infection
- Any circumstance where increased venous and lymphatic return is undesirable
- Suspected or known untreated acute DVT

Cautions

- Mild to Moderate Peripheral Arterial Disease
- Impaired Sensation - must be able to sense applied pressure

- Intolerance to material(s)
- Pediatric & assisted use - must be able to sense & communicate applied pressure

Washing instructions

- Machine wash warm - gentle cycle
- Do not bleach
- Tumble dry, low heat
- Do not iron
- Do not dry clean
- Do not wring

When washing it is recommended that the circaid reduction kit component be placed in a pillow case to keep the hook material lint-free and to prevent possible damage to other garments. Should lint collect on the tabs, it can be removed with tweezers.

Storage instructions

- Store in a dry place
- Keep out of sunlight

Lifespan

circaid reduction kit has a lifespan of 2 months.

Due to material wear-and-tear the medical efficacy can only be expected for a defined lifespan. This assumes correct handling of the garment (e.g. proper care, application, and removal).

Disposal

Please dispose of the garment in household waste. No special criteria for disposal apply for medical compression garments.

medi Ukraine LLC
Evhena Sverstiuka str, 11.
Kiev 02002
Ukraine
T: +380 44 591 11 63
F: +380 44 392 73 73
info@medi.ua
http://medi.ua

medi USA L.P.
6481 Franz Warner Parkway
Whitsett, N.C. 27377-3000
USA
T +1-336 4 49 44 40
F +1-888 5 70 45 54
info@mediusa.com
www.mediusa.com

This Package contains

- 1 circaid reduction kit lower leg
- 6 hook stays
- 1 pair of circaid undersleeves/undersocks
- 1 circaid reduction kit shelf strap
- 1 circaid BTS guide card
- 1 paper measuring tape
- 1 booklet instructions for use

Material composition

circaid compression garment

Polyurethane 54%

Nylon 42%

Elastane 4%

For all other components, please refer to the inserts included with each item.

This garment does not contain latex.

Individualisation**Step 1:**

Separate the body and spine from one another, unroll the bottom band and locate the number markings along the bottom and top edge of the body piece.

Step 2:

Locate the patient's ankle circumference measurement amongst the printed numbers along the bottom edge. Then, locate the patient's calf circumference measurement amongst the printed numbers along the top edge of the garment.

Step 3:

Take the spine and place the bottom left edge on the appropriate number. Take the top top left edge of the spine and place it on the appropriate number. Adjust the center spine tab to create a smooth spine curve.

Note:

There is no one correct spine shape. The garment may not lay flat once the spine has been set. This is to be expected when dealing with more extreme spine angles. If there is a large difference between the two circumferences, the spine may not reach both the bottom and top edges of the body piece. In this case, center the spine.

When dealing with more extreme spine angles, excess material sticking out beyond the top of the spine piece can be trimmed.

Step 4:

With the spine in place, fold over the spine bands exposing the body piece. Cut the excess material from the body piece leaving approximately 5 cm of overlap. Be sure not to cut any other part of the garment.

Note:

If desired, attach the separate hook tab stays along the spine for additional adhesion.

Tip:

You may also trace the spine to easily reattach it in its original position in the unlikely event of the spine becoming detached.

Step 5:

If the lower leg component is too long for the patient, it can be trimmed by 3 cm. Trim along the indicated cut line on the black side of the material. Trim both the body and the spine portion.

Donning instructions**Step 1:**

If the foot or toes require any bandaging, apply prior to donning the undersleeve and legging.

Step 2:

Slide the enclosed circaid undersleeve or undersock onto the leg up to the knee (for lower leggings) or to the groin (for upper leg and whole leg).

The undersleeve may be trimmed in length if desired.

If a knee and/or upper leg component will be applied it is not recommended to trim the undersleeve until these components have been fitted over the undersleeve to determine the excess length.

Step 3:

Align the garment just above the ankle bone and below the knee crease. The black side of the material should face the skin.

Unroll the top two bands and loosely secure them to the outside of the garment to anchor it in place. This enables you to easily raise or lower the garment for correct positioning.

Step 4:

Unroll the bottom two bands and secure the bottom most band to the outside of the garment to a firm and comfortable compression level.

Step 5:

While holding the second band, detach the next band. Secure the second band and continue this process with the other bands, going up the leg.

Ensure that the bands are placed in an alternating order.

Step 6:

Once all bands are secured, inspect the garment for any gaps or creases. Adjust the bands as necessary. If parts of the leg require additional support or a gap exists, use the provided shelf strap as needed.

Note:

If no reduction kit knee component or upper leg component are being worn, you may fold the top and bottom part of the circaid undersleeve over the top and bottom edge of the circaid compression garment.

Step 7:

Apply an appropriate foot compression garment as needed. If using another circaid garment, please refer to the specific product instructions included in its packaging.

Doffing instructions

Detach all bands of the legging starting with the top band. Fold each band onto itself to keep the hook and loop tabs lint-free.

Remove the legging and then remove the circaid undersleeve or sock.

Using the Built-In-Tension System (BTS)

Note:

Bands should be pulled to a firm but comfortable tension to provide therapeutic compression within the tolerance level of the patient.

Step 1:

Starting with the bottom band, line up the black triangle on the card with one of the BTS lines on the bottom band.

Step 2:

Note where the second BTS line lines up with the card's compression ranges.

Step 3:

If the BTS line on the garment either falls short or goes beyond the desired compression range, readjust the band as necessary so that the second BTS line on the garment is aligned with the desired compression range. Completely loosen

the band before reapplying it if too much tension has been applied.

Note:

Record tension setting for each set of markings so the patient can easily repeat the process and ensure it is donned to the same or higher tension.

Step 4:

Repeat the steps 1 - 3 for each band going up the leg.
It is not required to have all bands set to the same tension level. If more than one level of tension is appropriate for treatment, use the BTS card to identify and record the appropriate tensions levels.

Bands may need adjusting throughout the day to maintain the desired compression range as edema is reduced.

circaid reduction kit shelf strap

Material composition

55% Nylon
38% Polyurethane
4% Polyethylene
3% Elastane

This garment does not contain latex.

Donning the circaid reduction kit shelf strap:

Step 1:

Apply the appropriate circaid compression system as directed.

Step 2:

If using the reduction kit shelf strap, trim the strap to the desired lengths or keep as one long band, and secure hook tabs to the strap ends.

Step 3:

If there is an overhang of oedematous skin and additional support or coverage is desired, use the shelf strap to wrap up underneath the area in need, pulling up and around the limb until firm and comfortable compression is achieved and secure to the circaid compression system of choice.

Step 4:

If a gap between the bands occurs and cannot be corrected by adjusting the bands then cover the gap with the shelf strap pulling to firm and comfortable compression and securing the two hook tabs to the circaid compression system of choice.

Removing the circaid reduction kit shelf straps:

Step 1: Detach the two hook tabs and remove the shelf strap.

Note:

Do not tightly roll the hook tabs. This will cause the tabs to curl and no longer lay flat while in use.

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