
Instructions for Use

---

**COLD RUSH**<sup>®</sup>  
Cold Therapy System

3

EN | Instructions for Use





5

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

17

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30

## DESCRIPTION OF SYMBOLS

Symbol	Description
	Caution Symbol Indicates that important cautionary information, such as warnings and precautions, cannot be presented on the medical device itself and are presented in the instructions for use.
	Type BF Applied Part
	Manufacturer address
	Consult instructions for use
	This marking on the product, packaging, accessories, or literature indicates that the product contains electronic components and/or batteries that should not be disposed of in regular waste at the end of its usable life. Users should contact Össur or their local government office for information on how these items can be recycled or disposed of in an environmentally sound manner. To protect natural resources and to promote material reuse, please separate batteries and electrical components from other types of waste and recycle them.
<b>IP21</b>	Protected against solid foreign objects of 12.5 mmØ and greater. Protected against vertically falling water drops.
	Negative center connection (ISO. 5926)
	Preferred Output Polarity Symbol
	Class II Equipment
	For indoor use only
	UL's Recognized Component Mark for Canada and USA
	Wah Hing Transformer Mfy Ltd. - Power Supply Manufacturer

## **PRODUCT DESCRIPTION**

Cold Rush is a mechanical circulation cold therapy device. The device is composed of a bucket holding a mixture of water and ice, a circulation pump that forces the coolant through a closed-loop hydraulic circuit and a cooling pad that is used to apply cold therapy to the region that is being treated. The system operates on a 12-volt power which is supplied via an external wall-mounted power supply.

## **LIST OF PARTS (Figure 1)**

- Cold Rush Cold Therapy System
- Hose
- Couplings
- Flow Control Knob
- DC Connector
- 12V Power Supply
- Cold Rush Pad (Sold Separately)

## **INDICATIONS FOR USE**

- Provides cold therapy
- Indications requiring the application of cold therapy to a joint or limb such as post-surgical, trauma pain, and swelling

## **CONTRAINDICATIONS**

The Cold Rush Cold Therapy device should not be used on patients with the following conditions:

- Known hematological dyscrasias that predispose to thrombosis (e. g., paroxysmal cold hemoglobinuria, cryoglobulinemia, sickle cell disease, serum cold agglutinins)
- Compromised local circulation (including arteriosclerosis, ischemia)
- Neurologic impairment (including paralysis or localized compromise due to multiple surgical procedures or diabetes) in the affected region, polyneuropathy, or other nerve damage causing decreased skin sensitivity
- Cognition or communication impairments that prevent the user/patient from giving accurate and timely feedback including incapacitated patients with severe cardiovascular disease, anesthetic skin, hypercoagulation disorders, poor circulation, extremities sensitive to pain, extremely low blood pressure
- Raynaud's disease, phenomenon, or other vasospastic conditions
- Buerger's disease, significant vascular impairment
- Hypersensitivity to cold (cold urticaria) or history of cold related injury (including chilblains or frostbite)
- Severe cardiovascular disease, anesthetic skin, hypercoagulation disorders, poor circulation, extremities sensitive to pain, extremely low blood pressure, incapacitation, decreased skin sensitivity, vein ligation or recent skin grafts, or pheochromocytoma
- Localized unstable skin condition (e. g., dermatitis, vein ligation, gangrene, or recent skin graft) in the affected region or potential impairment of healing in the treatment area, including infection
- Tissues inflamed as a result of recent injury or exacerbation of chronic inflammatory condition

## **INTENDED PURPOSE**

The Cold Rush is a non-invasive device, intended for single patient use and designed to reduce localized pain and swelling following surgical procedures or trauma.

## **SAFETY PRECAUTIONS**

**WARNING:** The Cold Rush® cold therapy system can be cold enough to seriously injure skin. Read the instructions which are located inside the unit, the instructions for the pad, and the labels on the device carefully before operating the Cold Rush cold therapy system.

## **ELECTRIC SHOCK HAZARD**

- Unplug unit before filling it with ice and water. To avoid danger of electric shock be certain that hands are dry before inserting or removing the power supply from wall socket.
- Avoid getting water or ice on the flow control panel or power adapter.
- Always unplug this product immediately after each use.
- Do not use while bathing.
- Do not place or store product where it can fall or be pulled into a tub or sink.
- Do not place in or drop into water or other liquid.
- Do not reach for a product that has fallen into water. Unplug immediately.
- Never operate this product if it has a damaged electrical cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return the product to your healthcare provider for examination and repair.
- Keep cord away from heated surfaces.

**NOTE:** To create the safest operating conditions possible, Össur has designed the Cold Rush Cold Therapy System to operate on a 12-volt power supply, significantly reducing the chance of dangerous electric shock. Great care should always be taken when using electricity near water.

## **WARNING: RX ONLY.**

The device should only be used with a prescription from a physician that includes the following treatment information: (1) the duration of the entire cold therapy protocol, (2) the length of individual cold therapy sessions, as well as the length of breaks from cold therapy treatment, (3) setting of flow control knob, (4) instructions about how to inspect the skin, and (5) frequency of skin checks.

## **WARNING: SKIN INSPECTION AND ADVERSE REACTIONS.**

The patient or his healthcare provider must check the skin at least every 1–2 hours during use of this device, regardless of whether the patient is asleep or awake. The patient must discontinue the use of the Cold Rush Cold Therapy Unit if he experiences any adverse reaction to the treatment area such as: increase in pain, burning, blistering, itching, increased swelling, discoloration of the skin, increased redness, welts, or any other changes in the appearance of the skin.

## **WARNING: NO DIRECT CONTACT WITH SKIN: INSULATION BARRIER REQUIRED (NOT INCLUDED).**

The Cold Rush Pad should never make direct contact with the skin. An insulation barrier between the pad and the skin should be used at all times without exception in accordance with instruction from the physician so as to prevent the skin from becoming too cold. The Cold Rush pad operating temperature is too cold to be placed directly on the skin.

- Unit is not suitable for use in the presence of flammable anesthetic mixture with air, oxygen, or nitrous oxide.

- Follow local ordinances or regulations for proper disposal of device, accessories, and packaging.
- **CAUTION:** Federal law requires this device to be sold by a physician or by the order of a physician.

## WARNING

Read the following treatment warning protocol instructions carefully before operating the Cold Rush cold therapy system.

- The flow control knob should never be set at the maximum settings (number “4” or “MAX”) for more than 30 minutes.
- The Cold Rush pad should never be applied directly to the patient’s skin. Do not use the Cold Rush cold therapy system without first applying an insulation barrier between the Cold Rush pad and the patient’s skin.
- Patients should never utilize cold therapy while they are asleep or unable to check their skin every 1–2 hours unless specifically instructed to do so by the prescribing health care provider.
- Patients must check their skin every 1–2 hours once steady flow of water through the pad is achieved. Turn the flow control knob counterclockwise to increase the temperature (warmer), and clockwise to decrease the temperature (colder).
- When turning the unit on, turn the flow control knob clockwise to the number “3” setting wait 10 minutes to allow the temperature to stabilize before making any additional temperature adjustments.
- The following patients should ALWAYS turn the flow control knob counterclockwise to the number “7” setting when using cold therapy
  - Patients whose treatment area has been desensitized due to local anesthesia or regional nerve blockers
  - Patients taking prescribed medications that may negatively affect peripheral vascular circulation, including beta adrenergic blockers and local epinephrine use (such as in local anesthetics)
  - Patients taking prescribed medications that negatively affect their mental capacity or judgment
  - Patients suffering from diabetes
  - Patients who have recently undergone hand/wrist or foot/ankle surgery
  - Patients suffering from cognitive disabilities
  - Young children and the elderly
- **CAUTION:** Contact your healthcare provider immediately and discontinue the use of your Cold Rush Cold Therapy Unit if you experience any adverse reaction to the treatment area such as: increase in pain, burning, blistering, itching, increased swelling, discoloration of the skin, increased redness, welts or any other changes in the appearance of the skin.
- Make sure that the insulation barrier is dry.
- Care should be taken during prolonged use for children, diabetics, incapacitated patients and those with decreased skin sensitivity or poor circulation.
- A licensed healthcare practitioner must determine the correct frequency and duration of treatment for each patient.
- If prescribed by a licensed healthcare practitioner for home use, the healthcare practitioner must provide each patient with training for appropriate use and proper application of the Cold Rush Cold Therapy System.

- The health care practitioner must monitor the patient's use of the device to assure appropriate use and proper application of cold therapy, including but not limited to application of an insulation barrier between the patient's skin and the Cold Rush Pad.
- The product should never be left unattended when plugged in.
- Never drop or insert any object into any opening or hose.
- Do not operate where aerosol (spray) products are being used or where oxygen is being administered.
- Do not attempt to drag or reposition the Cold Rush Cold Therapy unit by pulling or tugging on the insulated hoses, as this may cause damage to the unit and/or cause the hoses to disconnect. Use handle when moving the Cold Rush Cold Therapy unit.
- If you experience disturbance in any other electrical equipment near the Cold Rush unit, please try to unplug it from the wall socket to verify that it is not the source of disturbance.

## **SAFETY STANDARDS**

Cold Rush™ Cold Therapy System is tested and certified to comply with the IEC60601-1 standard of electrical safety of medical devices and IEC/EN60601-1-2, electromagnetic compatibility for medical electrical devices. The company fulfills the requirements of ISO13485.

## **SETUP AND OPERATION**

1. Ensure cold therapy unit is turned off, turn flow control knob counterclockwise to the off position. Do not plug in to wall socket **(Figure 2)**
2. Fill reservoir with water to the indicated fill line **(Figure 3)**
3. Add chopped or cubed ice to the reservoir up to the indicated fill line **(Figure 4)**
4. Attach the provided 12V power supply to DC power port, do not plug into wall circuit **(Figure 5)**
5. Apply insulation barrier over the patient's skin in the area to be treated **(Figure 6a)**
6. Apply Cold Rush Pad over insulation barrier (not included). The Cold Rush Pad should never be applied directly over the patient's skin **(Figure 6b)**
7. Connect Cold Rush pad to hose extending from the unit. To ensure proper connection, push couplings together until they snap into place **(Figure 7)**
8. Plug the provided 12V power supply to wall socket (120V US, 240V Europe) **(Figure 8)**
9. Turn flow control knob clockwise to the number "3" setting, this will turn unit on **(Figure 9)**
10. Wait 10 minutes to allow temperature to stabilize before adjusting the temperature for comfort **(Figure 10)**
11. Turn flow control knob counterclockwise to increase temperature.
12. Turn flow control knob clockwise to decrease temperature.
13. To turn the cold therapy unit off, turn flow control knob counterclockwise until it stops.
14. Always turn the cold therapy unit off before disconnecting the Cold Rush pad or any hoses. To disconnect the hoses, press the release tabs on the coupling and pull the connectors apart.
15. Drain all of the water out of the cold therapy unit after each use.
16. Follow directions 2–10 before each use.

## **MAINTENANCE**

Do not try to repair your Cold Rush unit as this may create a user hazard and set basic safety of the user at risk in regard to electromagnetic disturbances. An attempt of repairing the device will automatically void warranty. After talking with your healthcare provider, to ensure you are following the operating instructions properly, the unit shall be sent to Össur for repairs or maintenance. No other company or individual is allowed to service or repair the unit in any way.

## **CLEANING AND PROLONGED STORAGE**

- When you have finished all of your cold therapy treatments, prior to prolonged storage, fill the unit with one gallon of water (or enough so the water level is at the ice fill line.
- Add 2 tablespoons of liquid bleach or 3 ounces of 3% hydrogenated peroxide to the water and mix.
- Connect hoses to the Cold Rush Pad, turn unit on to MAX setting.
- Allow water mixture to circulate for 5 minutes.
- Turn unit off and unplug the power supply.
- Carefully drain the water mixture out of the unit and store in a cool dark place leaving the lid slightly open. (Exposure to sunlight and extreme heat may damage hoses and unit.)
- When reusing the unit after storage, fill the reservoir with clean water and turn the pump on to circulate the water through the hoses and flush the remaining water mixture out of the system (use old pad for this procedure, discard the pad after cleaning is completed).
- Dispose water and refill the reservoir with fresh ice and water combination to begin cold therapy treatment. Follow operating instructions 2–10 before each use.
- Before you use the unit, inspect hoses and new pad for any tears or breaks. Contact your healthcare provider to obtain new Cold Rush pad.
- Clean or replace insulation barrier as recommended by your healthcare practitioner.

## **ELECTRICAL SPECIFICATIONS – POWER SUPPLY**

Do not use any other power supply other than the one provided with the unit.

- Model no: SA9505-120040US
- Input: 100-240V
- Input frequency: 50/60 Hz
- Output power: 4.8W
- Output Voltage: 12V
- Output Current: 0.4A

Or

- Model no: UES06WU-120050SPA
- Input: 120V US, 240V EU
- Input frequency: 60Hz US, 50Hz EU
- Output power: 6W
- Output Voltage: 12V
- Output Current: 0.5A

## ENVIRONMENTAL CONDITIONS

- Operating Temperature: +10°C to 40°C
- Operating Humidity: 30% - 75% Relative Humidity
- Atmosphere Pressure: 700h Pa – 1060 h Pa
- Shipping and Storage Temperature: -20°C to 60°C
- Shipping and Storage Humidity: 10% - 90% Relative Humidity, non- condensing

## ELECTROMAGNETIC COMPATIBILITY PRECAUTIONS AND WARNINGS

The Cold Rush Therapy System needs special precautions regarding electromagnetic compatibility (EMC). Specifically, it needs to be installed and put into service according to the EMC information provided as follows:

- The Cold Rush Therapy System should not be used adjacent to or stacked with other equipment. In case adjacent or stacked use is necessary, the Cold Rush Therapy System should be observed to verify normal operation in the configuration in which it will be used.
- The Cold Rush Therapy System may be susceptible to electromagnetic interference from portable and mobile RF communications devices such as mobile (cellular) telephones.
- The Cold Rush Therapy System may be interfered with by other equipment, even if that other equipment complies with CISPR EMISSION requirements.

Cold Rush accessories, various pads, will not negatively affect EMC performance of the Cold Rush Therapy System. To maintain basic safety, do not use any other power supply other than the one provided with the unit. Cold Rush Therapy System is intended for use in EM environment tabulated below.

<b>Guidance and Manufacturer's Declaration – Electromagnetic Emissions</b>		
The device is intended for use in the electromagnetic environment specified below. The customer or user of device should assure that it is used in such an environment.		
<b>Emissions Test</b>	<b>Compliance</b>	<b>Electromagnetic Environment – Guidance</b>
RF Emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonics Emissions IEC 61000-3-2	Class A	
Voltage Fluctuations/ Flicker Emissions IEC 61000-3-3	Complies	

## Guidance and Manufacturer's Declaration – Electromagnetic Emissions

The device is intended for use in the electromagnetic environment specified below. The customer or user of device should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV	±2 kV	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	±1 kV line(s) to line(s) ±2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, Short interruptions and voltage variations on power supply lines IEC 61000-4-11	0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°.  0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0°	0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°.  0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0°	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note:  $U_T$  is the a.c. mains voltage prior to application of the test level.

**Recommended separation distances between portable and mobile RF communications equipment and the device.**

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz d=1,2√P	80 MHz to 800 MHz d=1,2√P	800 MHz to 2,7 GHz d=2,3√P
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1	At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
NOTE 2	These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## Guidance and Manufacturer's Declaration – Electromagnetic Emissions

The device is intended for use in the electromagnetic environment specified below. The customer or user of the device should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:
Conducted RF IEC 61000-4-6	3 Vrms, 150 kHz to 80 MHz	3 Vrms	$d=1.2\sqrt{P}$
Radiated RF IEC 61000-4-3	10 V/m, 80 MHz to 2.7 GHz	10 V/m	$d=1.2\sqrt{P}$ 80 MHz to 800 MHz $d=2.3\sqrt{P}$ 800 MHz to 2.7 GHz
			where “P” is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).  Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey*, should be less than the compliance level in each frequency range**.
			Interference may occur in the vicinity of equipment marked with the following symbol: 
Note 1	At 80 MHz and 800 MHz, the higher frequency range applies.		
Note 2	These guidelines may not apply in situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.		
*	Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.		
**	Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.		

RF Wireless Communication Compliance Information					
Test Frequency (MHz)	Band (MHz)	Service	Modulation	Distance (m)	Immunity test level (V/m)
385	380-390	TETRA 400	Pulse modulation 18 Hz	0,3	27
450	430-470	GMRS 460, FRS 460	Pulse modulation 18 Hz	0,3	28
710	704-787	LTE Band 13, 17	Pulse modulation 217 Hz	0,3	9
745					
780					
810	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	0,3	28
870					
930					
1720	1700-1990	GSM 1800, CDMA 1900, GSM 1900, DECT, UMTS LTE Band 1, 3, 4, 25	Pulse modulation 217 Hz	0,3	28
1845					
1970					
d2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	0,3	28
5240	5100-5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0,3	9
5500					
5785					

## DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

## LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

## **LIMITED WARRANTY**

Össur warrants, to the original purchaser, that this product is free from defects in materials and workmanship for six: 6 months from initial purchase. This warranty does not apply if this product: 1 was not purchased from Össur or its authorized dealer,; 2 has been altered in any way, or: 3 has not been used in accordance with Össur's Instructions For Use (IFU). This Limited Warranty does not cover damage due to accidents, neglect, misuse, or operation beyond capacity, parts damaged by improper installation, substitution of parts not approved by Össur, or any other alteration or repair by others that, in Össur's judgment, materially or adversely affects the product or part.

For complete warranty guidelines please visit [www.ossur.com](http://www.ossur.com)

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