

RXCRUSH PILL CRUSHING SYSTEM


The RxCrush® Pill Crushing system is a sealed crushing system that pulverizes pills through its patented crushing and rolling action. Nurses are protected from hazardous drug exposure and clogged feeding tubes are eliminated. Dissolved medication is quickly and easily delivered to feeding tubes from the ENFit® compatible pill pouches.

Directions

Follow these steps to crush, dissolve, and deliver medication to feeding tube.









We're Here to Help

To reorder or for assistance, visit us at www.RxCrush.com or call +01.877.343.2783.



Move the pill pouch under roller and roll back and forth to pulverize pill. Syringe should be held to the side again.



Remove syringe and pouch from the crusher. Inject water from syringe into pill pouch.



Keep syringe attached to pouch and shake pill pouch until medication is completely dissolved. Take the extra time to ensure it's completely mixed.



Draw dissolved solution from pill pouch into syringe. Hold pouch up and syringe down to ensure contents flow into syringe.



Hold syringe up and pouch down so contents do not leak. Detach syringe from pouch. Attach syringe to feeding tube and inject solution from syringe into feeding tube.



Refill syringe with water, inject water into pouch, withdraw contents and inject into feeding tube as many times as necessary to remove all pill particles from pill pouch.

Cleaning

RxCrush Pill Crushing system is made of non-corrosive material and can be cleaned with water and disinfectant cleaners. Pill pouches are intended for single-patient use and are not intended to be cleaned or re-used.

Warranty

The manufacturer warrants the RxCrush Pill Crusher to be free of functional defects for five years. If the RxCrush Pill Crusher is used to crush anything other than pills in the RxCrush Pill Pouch, the warranty is null and void. The warranty is also void if damage is due to accident, misuse, or misapplication. To make a claim, contact us at www.RxCrush.com/contact.