




# MESALT EFFECTIVELY STIMULATES THE CLEANSING OF HEAVILY DISCHARGING AND DISCHARGING INFECTED WOUNDS

- Absorbs exudate, bacteria and necrotic material
- Effectively stimulates the cleansing of wounds
- Easy application

Absorbent viscose/  
polyester nonwoven

- Impregnated with sodium chloride.



## HOW MESALT WORKS


The discharge from the wound releases the sodium chloride from the dressing. Mesalt effectively stimulates the cleansing of wounds in the inflammatory phase by absorbing exudate, bacteria and necrotic material from the wound, thereby facilitating the natural wound healing process.

## AREAS OF USE


Mesalt is intended for the management of heavily discharging and discharging infected wounds in the inflammatory phase and deep cavity wounds such as:

- Pressure sores
- Surgical wounds


## HOW TO USE MESALT



1. Clean the wound and open the package.



2. Apply or loosely pack Mesalt (dry) to the wound.  
Do not soak beforehand.



3. Apply Mepilex Border as cover dressing. Depending on the amount of exudate, Mesalt is normally changed daily, or more often if necessary.

## BENEFITS OF MESALT

- Absorbs exudate, bacteria and necrotic material
- Effectively stimulates the cleansing of wounds
- Easy application

## NOTE

Mesalt is not to be used in dry wounds. The use of Mesalt has been presented in papers and study reports. The data are available from Mölnlycke Health Care on request.

### MESALT ORDERING INFORMATION

Product	Art. No	Size cm	Pieces Per Inner	NHS Code
Mesalt	285580	5 × 5 cm	30	ELS085
	285780	7.5 × 7.5 cm	30	ELS086
	286080	10 × 10 cm	30	ELS087
Mesalt Ribbon	285280	2 × 100 cm	10	ELS088



Mölnlycke Health Care Ltd, Arenson Centre, Arenson Way, Dunstable, Bedfordshire LU5 5UL, UK. Tel. +44 (0) 870 6060766, [www.molnlycke.co.uk](http://www.molnlycke.co.uk)

The Mölnlycke Health Care name and logo, Mesalt® are registered trademarks of Mölnlycke Health Care AB  
©2009 Mölnlycke Health Care. All rights reserved.

UK122550611