DAY DÍA	1	2	3	4	5	6	7
2500ml							
2250ml							
1750ml							
1500ml							
1250ml							
1000ml							
750ml							
500ml							
250ml							
GOAL META							

VOLUME ATTAINED / VOLUMEN OBTENIDO

**Rx ONLY**  

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**Warning:** Single patient use: The use of products intended for single patient use on more than one patient creates a potential risk of serious injury and/or infection which may lead to death. Reprocessing of medical devices intended for single patient use only may result in degraded performance or a loss of functionality.

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Hudson RCI, Teleflex and Voldyne are registered trademarks of Teleflex Incorporated or its affiliates.

All other trademarks are the properties of their respective owners.

© Copyright 2019 Teleflex

Please scan here for instructional video.  
Por favor escanee aquí para ver el video instructivo.

For more information on incentive spirometry, please refer to the recommendations section of the 2011 AARC Clinical Practice Guidelines.

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The Hudson RCI® Voldyne® Volumetric Exerciser is intended as an inspiratory deep breathing positive exerciser for adult and pediatric (above 5 yrs) patients. It is a single patient, multi-use device used in hospital or home care setting.

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**1 Remove components from package, expand the tubing to its full length and attach open end of tubing to stem on front side of exerciser.**

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**2 Slide the green pointer to the prescribed volume level.**

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**3 Hold or stand exerciser in an upright position. Exhale normally. Then, place lips tightly around mouthpiece.**

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**4 Inhale slowly to raise the piston in the chamber. When inhaling, maintain the green flow chip in line with the smiley face flow range. Continue inhaling and try to raise the piston to prescribed volume level (top of piston indicates inspired volume).**

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**5 When inhalation is complete, remove mouthpiece, hold the breath as prescribed, and exhale normally. Allow piston to return to bottom of chamber. Rest and repeat exercises. Frequency of use and recommended inspiratory volumes should be at the direction of physician or therapist.**

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## CLEANING INSTRUCTIONS

- Remove the mouthpiece from the expandable tubing.
- Dilute a ratio of 1 mL of liquid Dawn® dishwashing detergent into 1 L of warm tap water in a clean tub or bowl.
- Immerse the mouthpiece in the soapy solution and soak for approximately 5 minutes.
- After soaking but while the component is immersed, wipe the interior and exterior of the mouthpiece with a

non-linting wipe for a minimum of 30 seconds. Rinse the mouthpiece under warm running tap water for a minimum of 1 minute or until all soap is removed.

- Dry the mouthpiece using non-linting wipes.
- Visually inspect the mouthpiece to ensure it is clean. Repeat cleaning steps if residual soil is present.
- Reconnect the mouthpiece to the expandable tubing.

**CAUTION/WARNING: After cleaning, visually inspect the corrugated tubing to ensure that the tubing is clean.**

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\*. Milliliter - Inspiratory capacity measured in milliliters rounded off to the nearest 50mL.

\*\* Formula used in the above Nomogram published in The American Review of Respiratory Diseases, official journal of the American Thoracic Society, September 1979, Vol. 120, Number 3, by G. Polgar and V. Promadhat.

† Inspiratory capacity measured in milliliters rounded off to the nearest 50 mL.

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