

# CanDo® Folding Pedal Exerciser

item # 10-0718

## For Safety:

- Consult your doctor before starting any exercise program.
- Keep children away from Pedal Exerciser while in use.
- Ensure that the Pedal Exerciser is positioned on a solid flat surface.
- **Do NOT** wear loose clothing that could become caught or trapped during exercise.
- Please exercise at your own pace and **Do NOT** over-exercise which may cause injury.
- **Do NOT** exercise in the standing position.
- Capacity: 150 lbs.

## Maintenance:

- Regularly check to see all the nuts, bolts and fittings are securely tightened. Periodically, check all moving parts for signs of wear or damage.
- Clean with a dry cloth. **Do NOT** use solvent cleaners.

## Usage:



### **unfold to use:**

pull each leg all the way open until the push button pops to lock in place



### **fold to store:**

push and hold the push button on both legs and bend each leg to fold for storage



**Fabrication Enterprises Inc.**  
PO Box 1500 White Plains, NY 10602 (USA)  
tel: 1-914-345-9300 • 800-431-2830  
fax: 1-914-345-9800 • 800-634-5370  
[www.Fab-Ent.com](http://www.Fab-Ent.com)

**CE** Authorized CE representative:  
RMS UK Ltd.  
28 Trinity Road  
Nailsea, Somerset BS48 4NU (UK)