




# Male Catheterization Simulator LF00855U Instruction Manual

To view the complete manual in English, Spanish, German, French, Italian, and Portuguese, please go to eNasco.com/healthcare and click on "Instruction Manuals."



Life/form. Products by Nasco

## About the Simulator

The Life/form® Male Catheterization Simulator is designed to duplicate the human condition as closely as modern plastics technology allows. Care and treatment should be used since abuse or rough treatment will damage the simulator almost as it would injure or cause pain to a patient. With responsible care, the simulator will last for thousands of catheterizations.

**Note:** To avoid the possibility of leakage, make sure to use the 16 French Foley catheter supplied with the simulator.

# **List of Components**

- Life/form<sub>®</sub> Male Catheterization Simulator
- Life/form<sub>®</sub> Simulator Lubricant
- · 16 French Foley Silicone Catheter
- Administration Unit

## Internal Structure

- · Simulated bladder
- Simulated mucosal folds
- Simulated bulbous urethra
- Simulator internal urethra sphincter

## General Instructions for Use

**Note:** Administration unit is a modified enema bag.

# A. To Fill the Administration Unit (Water Reservoir)

- Hang the administration unit 18" above the simulator. Use ring stand or similar device.
- 2. Close clamp on tube.
- Spread opening at top of bag and fill container with DISTILLED WATER.
- 4. Press to close zipper. Start at one end, then continue along length of zipper.

# B. To Prepare Simulator for Catheterization

- Place the Male Catheterization Simulator on the bottom of the legs with the plastic top in an upright position. Place the drain tube in a water-safe container.
- Connect the quick-disconnect fittings of the administration unit. Check to be sure the Dura-Clamp on the drain tube is closed. (See figure 1.)




Figure 1

- Open the clamp on the administration unit and the Dura-Clamp on the simulator drain tube.
- Allow water to flow out of the simulator drain tube to be sure tubing is clear of any obstructions. Close Dura-Clamp. The system is now filled with water and pressurized.
- Place simulator on its back for proper catheterization position.

### C. Catheter Selection

Use only 16 French Foley catheters. The Male Catheterization
 Simulator is precisely engineered to provide the most lifelike experience possible. Use of the 16
 French Foley catheter facilitates catheterization and helps minimize damage to the simulated urethra. Use of larger catheters results in unnecessarily difficult catheterization procedures. Use of smaller sizes will result in leakage around the catheter.

Note: Special care should be taken when using a Foley catheter. Cuff inflation should only be attempted when the catheter is in the proper position inside the bladder. Just as in a real patient, the cuff must be completely deflated before the catheter is removed. Improper use of a Foley catheter can result in damage to the simulator.

## D. Lubrication

- Lubricate the simulator urethra and catheter liberally EVERY TIME a catheterization procedure is attempted. This provides maximum realism and minimizes damage to the urethra wall of the simulator.
- Use only Life/form® Simulator Lubricant (LF00985U) or Ivory® liquid soap. Make NO OTHER substitutions. Even water-soluble lubricants such as K-Y Jelly tend to build up in the urethra over a period of time.
- Initial use after long storage or heavy usage during a long class session will require heavy lubrication. Do this by inserting the lubricated catheter about half way, then remove and relubricate several times before inserting the catheter fully. This procedure facilitates a successful initial catheterization by students.

# Catheterization Procedure

A. Lubricate catheter by coating with <u>Life/form</u>® simulator lubricant or Ivory® liquid soap. (See figure 2.)



Figure 2



Figure 3

B. Slowly insert the catheter about 2" into urethra with the penis in the normal hanging position until the restriction simulating the mucosal fold is encountered. (See figure 3.)



Figure 4

C. Withdraw the catheter a little and stretch the penis slightly to straighten the blockage of the mucosal fold. (See figure 4.)



Figure 5

D. Proceed by inserting the catheter another 2". You have now approached the bulbous urethra. Elevate the penis about 60° to straighten the urethra and proceed with the insertion. (See figure 5.)



# Figure 6

E. The next restriction is the simulated sphincter entering the bladder. With experience, you will feel a "pop" as the catheter passes into the bladder. (See figure 6.) It will require about 12" of catheter. Continue insertion until water flows from the catheter.

# Causes for Failure in Function

- Lack of sufficient lubrication see <u>General Instructions for</u> Use — D. Lubrication.
- Water does not flow even after proper insertion check administration unit water level. If there is still no water, remove catheter and repressurize the system following <u>General Instructions for Use B. To Prepare Simulator for Catheterization</u>.

### Care of Simulator

Normal soil on the surface of the simulator can be removed with mild soap and warm water. Use Nasco Cleaner (LF09919U) to remove stubborn stains from simulator. Simply spray soiled area and wipe clean with cloth or paper towels.

Do not store simulator for extended periods of time with water in the unit. Before returning simulator to the case,



# Figure 7

drain as much liquid as possible from the bladder. To do this, place the simulator upright on the lid. Open the Dura-Clamp, then close the clamp on the administration unit and pull the quick-disconnect apart over a water-safe container. (See figure 7.) Any residual water inside the simulator should drain out through the quick-disconnect fitting. Dry all surfaces to prevent mildew formation while stored in the case.

Before storage, be sure to remove the catheter from the urethra. The lubricant can dry out and act as a mild adhesive if the catheter is left in place. This results in damage to the simulator when removed.

# **Cautions**

Solvents or corrosive materials will damage the simulator. Never place the simulator on any kind of printed paper or plastic. These materials will transfer indelible stains. Ball-point pens will also make indelible stains.

# Supplies/Replacement Parts for Male Catheterization Simulator

**LF00985U Life/form**® Simulator Lubricant Kit

LF09919U Nasco Cleaner

**LF01127U** Foley Urethral Catheters

(pkg. of 1)

**LF01128U** Foley Urethral Catheters

(pkg. of 10)

# Other Available *Life/form* Simulators

Officer F	Available <b>Breporti</b> ® Simi	Jidiors	
LF00698U	Adult Injectable Arm (White)		Venatech IM & Sub Q
LF00855U	Male Catheterization	LF01193U	Special Needs Baby
LF00856U	Female Catheterization	LF03000U	CPARLENE® Series
LF00901U	Prostate Examination	LF03601U	Adult Airway Management
LF00906U	Ostomy Care		Trainer with Stand
LF00929U	Surgical Bandaging	LF03602U	Adult Airway Management
LF00957U	Enema Administration		Manikin
LF00958U	Pediatric Injectable Arm	LF03609U	Child Airway Management
LF00961U	Intramuscular Injection		Trainer with Stand
LF00984U	Breast Examination		Child <b>CRISIS™</b> Manikin
LF00995U	Arterial Puncture Arm	LF03617U	Deluxe Child <b>CRISis™</b>
LF00999U	Pediatric Injectable Head	1500/0011	Manikin with Arrhythmia Tutor
LF01005U	First Aid Arm		PALS Update Kit
LF01008U	Intradermal Injection Arm	LF03623U	Infant Airway Management
LF01012U	Heart Catheterization (TPN)	1502/2011	Trainer with Stand
LF01019U	Ear Examination	LFU3032U	Child Intraosseous Infusion/ Femoral Access Leg on a Stand
LF01027U	Peritoneal Dialysis	I E0343311	Child Airway Management
LF01028U	Suture Practice Arm	LFU3U33U	Trainer Torso
	Suture Practice Leg	LF03693U	Basic Buddy® CPR Manikin
	Spinal Injection		"Airway Larry" Airway
	Hemodialysis Practice Arm	<b>2</b> 1 000770	Management Trainer
	Episiotomy Suturing Set	LF03709U	Infant <b>CRiSis™</b> Manikin
LF01042U			Baby Buddy™ Infant CPR Manikin
	Pelvic, Normal & Abnormal		Fat Old Fred
	Stump Bandaging, Upper		Airway Management/Cricoid
	Stump Bandaging, Lower		Pressure Trainer
	Cervical Effacement	LF03770U	Chest Tube
	Birthing Station	LF03953U	<b>CRISIs™</b> Manikin, Complete
	Cricothyrotomy	LF03955U	Deluxe <i>CRiSis</i> ™ Manikin
	Tracheostomy Care	LF03956U	Deluxe "Plus" <b>CRiSis™</b> Manikin
LF01084U	Sigmoidoscopic	LF03965U	Adult <b>CRISIs™</b> Auscultation
150100711	Examination		Manikin
	Central Venous Cannulation	LF03966U	Adult <b>CRISIs™</b> Auscultation
	Blood Pressure Arm		Manikin with ECG Simulator
	Infant Intraosseous Infusion		<b>GERi™/KERi™</b> Manikin Series
-	Advanced IV Arm	LF04200U	Adult Sternal Intraosseous
	Venipuncture and Injection Arm		Infusion
	Advanced IV Hand Auscultation Trainer	LFU6001U	CPR Prompt® Adult/Child
		I E04012!!	Manikin CPR Prompt® Infant Manikin
	Testicular Exam		CPR Prompt® Keychain
	Male & Female Catheter	LFU0ZUUU	Rescue Aid
	Advanced CPR Dog Venatech IV Trainer	I E0620411	CPR Prompt® Rescue and
	NG Tube & Trach Skills	LIVUZUAU	Practice Aid
LFU11/4U	ING TUDE & ITACH SKIIIS		Tructice Alu


asco Fort Atkinson

901 Janesville Avenue, P.O. Box 901 Fort Atkinson, Wisconsin 53538-0901 1-800-558-9595 eNasco.com • E-mail: lifeform@eNasco.com





# Female Catheterization Simulator LF00856U Instruction Manual



Life form. Products by Nasco

# About the Simulator

The Life/form® Female Catheterization Simulator is designed to duplicate the human condition as closely as modern plastics technology allows. Care and treatment during use should be the same as with a patient, since abuse or rough treatment will damage the simulator almost as it would injure or cause pain to a patient. With reasonable care, the simulator will last for thousands of catheterizations.

**Note:** To avoid the possibility of leakage, make sure you use the 16 French Foley catheter supplied with the simulator.

# **List of Components**

- Administration Unit
- Plastic Pad
- 16 French Foley Catheter
- Life/form® Lubricant Kit

# **General Instructions for Use**

# A. To Fill the Flexible Administration Unit (water reservoir)

**Note:** Administration unit is a modified enema bag.

- Hang the flexible administration unit 18" above the simulator. Use ring stand or similar device.
- 2. Close clamp on tube.
- Spread opening at top of bag, then fill container with DISTILLED WATER
- Press to close zipper. Start at one end, then continue along length of zipper.

# B. To Prepare Simulator for Catheterization

- The Female Catheterization Simulator should be placed without the plastic pad in position. The pad will be used later. The legs should rest on a flat surface.
- Hang the administration unit 18" above the simulator. Do not hang higher because the increased pressure can cause leakage.



Figure 1

- Connect the quick-disconnect fittings of the administration unit to the simulator. Check to be sure the Dura-Clamp<sup>®</sup> (drain) on the simulator is closed. (See figure 1.)
- Open clamp on tube of administration unit.
- Open Dura-Clamp® (drain) on simulator until distilled water flows out of plastic tubing. (See figure 2.) Close Dura-Clamp® (drain). The system is now full, pressurized, and ready to use.



Figure 2

 Position the plastic pad between legs and under buttocks to allow the proper angle of the simulator for catheterization.

# C. Catheter Selection

Use only 16 French Foley catheters. The Female Catheterization Simulator is precisely engineered to provide the most lifelike experience possible. Use of the 16 French Foley catheter facilitates catheterization and helps minimize damage to the urethra of the catheterization simulator.

Use of larger catheters results in unnecessarily difficult catheterization procedures. Use of smaller catheters results in leakage around the catheter.

**Note:** Special care should be taken when using a Foley catheter. Cuff inflation should only be attempted when it is in the proper position inside the bladder. Just as in a real patient, the cuff must be completely deflated before the catheter is removed. Improper use of a Foley catheter can result in damage to the simulator.

# D. Lubrication

- Lubricate the simulator urethra and catheter liberally EACH TIME a catheterization procedure is attempted. This provides maximum realism and minimizes damage to the urethra wall of the simulator.
- Use only Life/form® Lubricant Kit (LF00985U) or Ivory® liquid detergent. Make NO OTHER substitutions. Even water soluble lubricants such as K-Y Jelly tend to build up in the urethra over a period of time.
- 3. Initial use after long storage or heavy usage during a long class session will require heavy lubrication. Do this by inserting the lubricated catheter about halfway, then remove and relubricate several times before inserting the catheter fully. This procedure facilitates a successful initial catheterization by students.

# How to Perform Catheterization on This Simulator

- Lubricate catheter.
- Slowly insert catheter through simulated sphincter. With experience, you will feel a "pop" as the catheter passes into the bladder. The simulator will require about 3½" of catheter. Continue insertion until water flows from catheter.

### Care of Simulator

Normal soil accumulated on the surface of the simulator can be removed with mild soap and lukewarm water. Use Nasco Cleaner (LF09919U) to remove stubborn stains from simulator. Simply spray soiled area and wipe clean with a soft cloth or paper towels.



# Figure 3

Do not store simulator for extended periods of time with water in the unit. Before returning simulator to the case, drain as much liquid as possible from the bladder. To do this, open the Dura-Clamp® (drain), holding it upside down and high in the air. Then close the clamp on the administration unit and pull the quick-disconnect fittings apart over a tray. (See figure 3.) Dry all surfaces to prevent mildew formation in the case.

Remove the catheter from the urethra after using. If catheter is left in, the lubricant will dry out and act as a mild adhesive. This results in damage when the catheter is finally removed.

#### Cautions

Solvents or corrosive materials will damage the simulator. Never place simulator on any kind of printed paper or plastic. These materials will transfer indelible stains. Ball-point pens will also make indelible stains.

# Supplies/Replacement Parts for Female Catheterization Simulator

LF00985U Life/form® Lubricant Kit

LF09919U Nasco Cleaner

LF01127U Foley Urethral Catheter.

Pkg. of 1

**LF01128U** Foley Urethral Catheters. Pkg. of 10

# Other Available Life form, Simulators

Other F	Available <i>Liepolii</i> , Simi	ulators	
LF00698U	Adult Injectable Arm (White)	LF01184U	Venatech IM & Sub Q
LF00855U	Male Catheterization	LF01193U	Special Needs Baby
LF00856U	Female Catheterization		CPARLENE® Series
LF00901U	Prostate Examination	LF03601U	Adult Airway Management
LF00906U	Ostomy Care		Trainer with Stand
	Surgical Bandaging	LF03602U	Adult Airway Management
	Enema Administration		Manikin
LF00958U	Pediatric Injectable Arm	LF03609U	Child Airway Management
LF00961U	Intramuscular Injection		Trainer with Stand
LF00984U	Breast Examination		Child <b>CRISIS™</b> Manikin
LF00995U	Arterial Puncture Arm	LF03617U	Deluxe Child <b>CRISis™</b>
LF00999U	Pediatric Injectable Head		Manikin with Arrhythmia Tutor
LF01005U	First Aid Arm		PALS Update Kit
LF01008U	Intradermal Injection Arm	LF03623U	Infant Airway Management
	Heart Catheterization (TPN)	1 = 00 / 0011	Trainer with Stand
LF01019U	Ear Examination	LF03632U	Child Intraosseous Infusion/
LF01027U	Peritoneal Dialysis	150242211	Femoral Access Leg on a Stand
LF01028U	Suture Practice Arm	LFU3033U	Child Airway Management Trainer Torso
LF01034U	Suture Practice Leg	1 EU340311	Basic Buddy® CPR Manikin
LF01036U	Spinal Injection		"Airway Larry" Airway
LF01037U	Hemodialysis Practice Arm	LF030770	Management Trainer
LF01038U	Episiotomy Suturing Set	LF03709U	Infant <b>CRISis™</b> Manikin
LF01042U	Suture Kit		Baby Buddy™ Infant CPR Manikin
	Pelvic, Normal & Abnormal		Fat Old Fred
	Stump Bandaging, Upper		Airway Management/Cricoid
LF01064U	Stump Bandaging, Lower		Pressure Trainer
LF01069U	Cervical Effacement	LF03770U	
LF01070U	Birthing Station	LF03953U	<b>CRISIS™</b> Manikin, Complete
	Cricothyrotomy		Deluxe <b>CRISIS™</b> Manikin
	Tracheostomy Care	LF03956U	Deluxe "Plus" <b>CRISIS™</b> Manikin
LF01084U	Sigmoidoscopic	LF03965U	Adult <b>CRISIS™</b> Auscultation
	Examination		Manikin
	Central Venous Cannulation	LF03966U	Adult <i>CRiSis</i> ™ Auscultation
	Blood Pressure Arm		Manikin with ECG Simulator
	Infant Intraosseous Infusion		<b>GERi</b> <sup>™</sup> / <b>KERi</b> <sup>™</sup> Manikin Series
	Advanced IV Arm	LF04200U	Adult Sternal Intraosseous
	Venipuncture and Injection Arm		Infusion
	Advanced IV Hand	LF06001U	CPR Prompt® Adult/Child
	Auscultation Trainer		Manikin
	Testicular Exam		CPR Prompt® Infant Manikin
	Male & Female Catheter	LF06200U	CPR Prompt® Keychain
	Advanced CPR Dog	. = 0 / 0 0	Rescue Aid
	Venatech IV Trainer	LF06204U	CPR Prompt® Rescue and
LF01174U	NG Tube & Trach Skills		Practice Aid

asco Fort Atkinson

901 Janesville Avenue, P.O. Box 901 Fort Atkinson, Wisconsin 53538-0901 1-800-558-9595