
# Airflo Mk3

user instructions

English

Deutsch

Français

Nederlands

Italiano

Español

Polski

*Simple solutions for everyday independence*

# English

## Intended Use

The Airflo Mk3 by Mangar International, is a low pressure air compressor, which has been designed for use with Mangar International pneumatic lifting equipment and is controlled by means of the supplied 2-way hand control.

## Applied Parts (Accessories) for the Airflo Mk3:

Bathing Cushion  
Surfer Bather  
Handy Pillow Lift  
Sit-u-Up  
Leg Lifter

English	Pg 2
Deutsch	S 10
Français	P 18
Nederlands	P 26
Italiano	Pag 34
Español	Pg 42

## Contents


Standard Package Contents	Pg 3
Getting Started	Pg 3
Battery Status Indicator	Pg 3
Battery Charging	Pg 4
Connecting the Air Hose	Pg 4
Push-Button Hand Control	Pg 5
Cleaning & Storage	Pg 5
Maintenance	Pg 6
Trouble shooting guide	Pg 7
Specifications	Pg 7
Symbols used on Airflo and Charger	Pg 8
Symbols used on packaging	Pg 8
Warranty	Pg 9

## WARNINGS

- **BE SURE** that you have read and understood these instructions and those provided with any other mangar international product, before operating the equipment.
- **DO NOT** place the battery charger in a bathroom or in a damp environment.
- **ALWAYS** disconnect from the battery charger before cleaning.
- **USE ONLY** the supplied battery charger.
- **NOT SUITABLE** for use in the presence of flammable anaesthetic mixtures with air, oxygen or nitrous oxide.
- Ensure that when in use, the Airflo and its connected cable and/or hoses are laid out in such a way as to avoid any risk of tripping or strangulation.
- Ensure that children and pets do not interfere with the safe operation of the Airflo and inflatable product.
- Keep the Airflo and inflatable product away from sources of heat and any sharp surfaces.
- The Airflo should only be used for it's stated purpose and should not be modified in any way.
- Disconnect the Airflo from the inflatable product before transporting so as to avoid accidental inflation.

*This product and associated packaging does not contain latex.*

### Standard Package Contents - fig. ①

- 1 x Airflo compressor fitted with one 12v- 2.3Ah battery (A)
- 1 x Battery charger (B)
- 1 x 2-Way hand control (C)
- 1 x User instructions

### Getting Started

The Airflo does not have an 'on/off' switch. To operate the compressor or activate the Battery Status Indicator, connect the hand control and momentarily press either the 'Up' or 'Down' button. If the Airflo is not used for approximately 5 minutes it will go into 'sleep' mode and the battery status indicator light will go out.

### Battery Status Indicator - fig. ②

Before use, check that the battery is sufficiently charged (see 'Getting Started'). If the battery is fully charged, the indicator light (A) will show a steady Green light.

**Steady Green** = Fully charged.

**Flashing Green** = Recharge soon.


**Steady Red** = Recharge immediately (will not deflate).

**Flashing Red** = Over temperature cut-out (will not deflate or inflate).

NOTE: If charge is too low the 'down' button on the hand control will not work.

The Airflo needs to be recharged when the indicator light is flashing Green.

**Charge the battery overnight before first using the Airflo.**

**Battery Charging - fig. ③**


- Check the battery charger casing and lead for signs of damage and if damaged DO NOT USE.
- Fit the appropriate mains plug adaptor to the battery charger (A).
- Plug the battery charger into a mains outlet socket and connect the charger output lead (B) to the Airflo charging socket (C).
- The battery status indicator light (D) on the Airflo will flash green while the battery is charging and then turn to a steady (non flashing) light when the Airflo is ready for use.

**NOTE:** The charger is not capable of providing sufficient power to charge the battery while the compressor is running.

To maintain battery condition, It is recommended that the Airflo is left continuously on charge when not in use. The battery cannot be overcharged using this charger.

**WARNING:** DO NOT recharge the battery in a bathroom or in a damp environment.


Unplug the charger from the mains power socket to disconnect the charger from the mains power supply.


**Connecting the Air Hose - fig. ④**

Disconnect the Airflo from the battery charger and place in a safe and convenient position.

Connect the air hose from the Mangar International equipment to the air outlet socket (A). The connector will 'click' into place when firmly pushed onto the outlet. To release the hose, press the paddle down on the connector (B).

### 2-Way Hand Control - fig. ⑤

Connect the hand control (A) to the socket (B). DO NOT connect any other plug to this socket.

- Press the 'Up' button (C) to raise.
- Press the 'Down' button (D) to lower.

NOTE: When either the 'Up' or 'Down' buttons are pressed, the internal microprocessor carries out a battery voltage check and if the voltage becomes too low, the 'Down' button will not work. This feature will always ensure that the seat of a Mangar International bathing product cannot be lowered into a bath if there is insufficient battery power to raise it again.

When inflating Mangar International products which are not fitted with a pressure relief valve, the compressor will stop automatically once the product is fully inflated.

### Cleaning

Unless contamination is suspected, periodic cleaning is all that is recommended. The Airflo should be cleaned/disinfected as follows:

- Unplug the Airflo from the battery charger. Clean the external surfaces using a non abrasive damp cloth with warm water. Wipe dry thoroughly with a dry towel before reconnecting to the battery charger. Do not allow moisture to enter any part of the compressor.
- Do not use solvent-based products.

### Disinfecting

- Pre-clean using a damp cloth and then disinfect using a hard surface alcohol impregnated wipe.

### Storage

- When not in use, store in a dry place together with these instructions. To maintain the battery when being stored for long periods, put the Airflo on charge overnight every 3 months.



**Maintenance - fig. ⑥**

At approximately 12 month intervals, inspect the air filter located in the base of the Airflo (**A**), and if necessary, fit a new filter (a kit of replacement filters is available through your supplier - part number: CA0074).

Battery replacement (part number: CB0405)

- Turn the Airflo upside down and remove the two screws (**B**) from each of the two retaining plates (**C**) and rubber pads. Turn the Airflo the right way up to remove the battery.
- Fit the replacement battery pack ensuring that the two contacts (**D**) align with the two contacts in the battery housing.
- Refit the rubber pads, retaining plates (**C**) and screws (**B**).

**IMPORTANT:** Always replace the battery with the same type.

⑥

When it becomes apparent that the Airflo is producing fewer than expected lifting cycles, the battery is probably reaching the end of its working life and should be replaced.

The Airflo contains no other user serviceable parts. If a fault should develop or if you require assistance in setting up using or maintaining the equipment, please contact our Customer Service Department on **Freephone 0800 2800 485**. A service manual (MI0338) is available for service personnel.

**Disposal**

When the Airflo has reached the end of its working life, remove the batteries for safe disposal in accordance with local regulations, and dispose of the Airflo at your local household waste recycling centre (civic amenity site).

### Trouble shooting guide

<p> Compressor does not run when 'up' is selected on the 2-way hand control.</p> <p> Ensure that the air hose from the product being inflated is connected to the Airflo.</p> <p> Compressor keeps stopping when trying to inflate the lifting aid.</p> <p> Check that there are no kinks in the air hose.</p> <p> Battery State Indicator light flashing 'Red/Green'.</p> <p> This happens when the current overload cut-out has tripped. Leave for 3 seconds to self reset.</p>	<p> Battery State Indicator light flashing 'Red'.</p> <p> This happens when the compressor's internal temperature cut-out has tripped. Leave the Airflo for 60-90 minutes to cool.</p> <p> Compressor causes electrical interference with another piece of electrical equipment.</p> <p> Move either the compressor or other item further away from each other until the interference stops.</p>
---	--

### Specifications

#### Airflo

Weight: 4.6Kg (10lbs)

Output pressure: 0.35 bar (5 psi)

Battery: 1 x 12Vdc- 2.3Ah (Maintenance free Valve Regulated Sealed Lead Acid)

Battery Life: Approx 1000 charging cycles (2- 5yrs Depending on usage)

Lifting performance (*from full battery*): Bathing Cushion- 4 bathing cycles;

Surfer Bather- 10 bathing cycles; Handy Pillow Lift- 11 lifts; Sit-u-Up- 20 lifts; Leglifter- 8 lifts

Environmental protection:

*(Airflo housing)* IP51 (Protection against large dust deposits & vertical falling drops of water)

*(Push-button hand control)* IP67 (Totally protected against dust & immersion in water)

Type BF equipment

This equipment has been tested and meets the requirements of current EMC legislation.



**Charger: Model number** - UE10WCP1-150035SPA

Input: AC 100v- 240v, 50-60Hz, 500mA




Output: DC 15v, 350mA

Double insulated Class II equipment

**Symbols used on Airflo and Charger labels**

	Warning		For indoor use only. Do Not connect to mains power in a bathroom or outdoors.
	Attention consult instructions		Duty cycle I = Operating time O = Non operating time (Max. continuous operating time = 5 minutes)
	Attention consult instructions		Max/Min operating environment temperature
	Type BF equipment (protection against electric shock)		Recycle
	Direct current (DC)		Do not dispose of in household waste
	Alternating current (AC)		Date of manufacture
	Risk of electric shock		'Up' button
	Class II Double Insulated equipment		'Down' button

**Symbols used on packaging**

	Keep packaging dry		Recycle packaging
	Handle with care		Do not dispose of in household waste
	Max/Min storage temperature		Max/Min humidity when transporting or storage
	Net & gross weights		Max/Min pressure when transporting or storage
	Stack no more than five high		Max/Min temperature when transporting or storage

## Warranty

The guarantee is given against defects in material and workmanship and comprises piece, parts and labour. The guarantee is not valid in the case of fair wear and tear, wilful damage, common neglect, misuse or unauthorised alteration or repair by unqualified persons.

The following statements do not affect the consumer's statutory rights. This product is sold on the understanding that in the event of any defect of manufacture or material appearing within the warranty period, the defect will be rectified free of charge providing that:

- (a) Reasonable evidence is provided (e.g. purchase invoice, delivery note) that the claim is being made within the warranty period.
- (b) The defect is not attributable to accidental damage (either in transit or otherwise), misuse or unauthorised repair.

Airflo serial number: .....

Date of purchase: .....

Invoice/Delivery note number: .....

We strive to provide our customers with reliable and durable products of the highest quality. In addition to our standard warranty, and for added peace of mind, our Customer Service Department offers an **Annual Service and Clean package**.

To find out more about this service please call our Customer Service Department on **Freephone 0800 2800 485**.

Mangar International reserves the right to alter product specifications and/or any of the information contained within this document without notice.

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- Bathing Cushion
- Surfer Bather
- Handy Pillow Lifgt
- Sit-u-Up
- Leg Lifter

English	Pg 2
Deutsch	Pg 10
Français	Pg 18
Nederlands	Pg 26
Italiano	Pg 34
Español	Pg 42

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The  
Airflo Mk3 Compressor - CB0025

**complies with:-**

*The requirements for Class 1 devices according to Annex VII of the Medical Devices Directive 93/42/EEC as amended by 2007/47/EC*

*IEC 60601-1  
Medical Electrical Equipment - General requirements for basic safety and essential performance*

*IEC 60601-1-11  
Medical Electrical Equipment - Requirements for use in the home care environment*

*IEC 60601-1-2  
Medical Electrical Equipment - Electromagnetic compatibility*

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**Mangar International Ltd**