# Installation and Operating Instructions

## ***Model FA4-MC***

### *Fire Alarm with Voltage Input 418MHz RF Transmitter*

#### **Transmitter Features:**

FA4-MC transmitters operate using auxiliary power source. FA4-MC is coded at the manufacturer with a unique address ID. The test and low battery LED are not operational on these units. With proper voltage input the RED LED will illuminate and the FA4-MC will transmit a Fire alert to the Silent Call Medallion series receiver.

#### **Transmitter Mounting And Location:**

The FA4-MC Medallion Series Transmitter is designed for indoor use only. The transmitter case has provisions for mounting to a nail or screw on the wall. Velcro strip is also provided for mounting on desktop or wall. Avoid mounting transmitter near large metal objects, other transmitting devices in close proximity or high current appliances that will degrade the transmitter / receiver operation. Test transmitter mounting location before mounting to ensure proper system operation. Adjust antenna for best results. Transmission range with a Silent Call Medallion Series receiver is approximately 500 feet.

#### **Installation:**

**Please contact a licensed fire alarm installer for installation and working with this product. Improper installation can damage your fire alarm and/or prevent proper operation of the fire alarm. Technical support is available at (800) 572-5227.**

**The FA4-MC in conjunction with a Silent Call Medallion Series receiver is considered a *Supplementary Visible Notification Appliance*.**

**The FA4-MC is designed to interface with a filtered supervised fire panel notification appliance circuit. 12VDC to 24VDC Reverse Polarity voltage will power the transmitter for supervision reporting to Silent Call Receiver if applicable.**

1. Identify Fire Panel NAC specifications. **FILTERED PANELS ONLY.** NAC alert voltage 12VDC to 24VDC positive polarity. 12VDC to 24VDC reverse polarity required for Silent Call supervision.
2. Connect the FA4-MC to a switched 12 or 24VDC voltage output from the fire panel NAC-- OR -- to the 24V horn / strobe output.
3. Follow Fire panel wiring requirements. Transmitter is a supervised power limited device.
4. Activate transmitter and check for proper receiver operation before mounting.
5. Transmitters interfaced with Fire Panel with no supervision voltage will require a non-supervised receiver. Test Button will not operate without reverse polarity voltage.
6. Pair FA4-MC intended Medallion Receiver. Follow Programming procedure in Receiver installation manual.
7. Reference Silent Call FA4-MC Installation Drawing 20-300-022-01FA4.

#### **Operation:**

1. FA4-MC will transmit when the input voltage has positive polarity. Fire Alarm active.
2. When the FA4-MC transmitter is activated the Red LED on the face of the transmitter lights up. The transmitter then sends an RF signal to a preprogrammed Silent Call Medallion Series Receiver activating the FIRE Alert on the Receiver.
3. Applicable receivers will supervise fire transmitter operation and alert the user to Transmitter/NAC fault.
4. Weekly testing of Auxiliary alerting system is recommended.

#### **Regulatory Information Notice:**

***THIS DEVICE COMPLIES WITH PART 15 OF THE FCC RULES.***

***This Device Complies With Industry Canada License-Exempt RSS Standard(s).***

***Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation of the device. This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:***

- ***Reorient or relocate the receiving antenna.***
- ***Increase the separation between the equipment and the receiver***
- ***Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.***
- ***Consult the dealer or an experienced radio/television technician for help.***

**Unauthorized changes or modifications could void the user's authority to operate the equipment**

# Installation and Operating Instructions

## ***Model FA4-MC***

*Fire Alarm with Voltage Input 418MHz RF Transmitter*

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## Installation and Operating Instructions

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### **Technical Support**

For technical support on this or any other Silent Call product, please feel free to contact us. You can reach us by phone at 800-572-5227 (voice or TTY) or by Email at [support@silentcall.com](mailto:support@silentcall.com)

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### **Limited Warranty**

Your transmitter is warranted to be free of defects in material and workmanship for five years from the date of initial purchase. During that time, the unit will be repaired or replaced free of charge when shipped prepaid to Silent Call Communications. This warranty is void if the defect is caused by customer abuse or neglect.

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**Return for Serving To:**

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## ***Silent Call Communications***

5095 Williams Lake Road, Waterford Michigan 48329

(800) 572-5227 v/tty (248) 673-7360 fax

Website: [www.silentcall.com](http://www.silentcall.com) Email: [silentcall@silentcall.com](mailto:silentcall@silentcall.com)