# EN CPAP SILICONE FULL FACE MASK

REF:

LOT

Consult instructions for use

Not made with natural rubber latex.

Date of Manufacture:

Use By:

DeVilbiss  
HEALTHCARE

## INTENDED USE

The CPAP/VPAP Mask is intended to be used for patients who have been prescribed CPAP or VPAP therapy and are under the specific direction of a physician. It is used in conjunction with continuous positive airway pressure (CPAP/VPAP) devices in home, hospital and Institutional environments.

## CONTRAINDICATIONS

This mask may not be suitable for patients with the following conditions: impaired cardiac sphincter function, excessive reflux, impaired cough reflex and hiatal hernia. It should not be used if the patient is unresponsive, obtunded, unable to remove the mask, or uncooperative.

## WARNINGS

- This mask should only be used with CPAP systems recommended by your physician or respiratory therapist.
- Do not use this mask unless the CPAP or VPAP system is turned on and operating properly.
- Never block the mask's vent holes or exhalation port.

**NOTE:** CPAP systems are intended to be used with special masks with vent holes in either the connectors or the mask that allow continuous flow out of the mask. When the CPAP device is turned on and functioning properly, new air from the CPAP device flushes the exhaled air out through the mask's exhalation port. However, when the CPAP device is not operating, enough fresh air will not be provided through the mask and exhaled air may be rebreathed. Rebreathing of exhaled air for several minutes can in some circumstances lead to suffocation.

- If oxygen is used with the CPAP system, the oxygen flow must be turned off when the CPAP device is not in operation.

**NOTE:** When the CPAP device is not in operation and the oxygen flow is left on, oxygen delivered into the ventilator tubing may accumulate within the CPAP device's enclosure and create a risk of fire.

**NOTE:** At a fixed flow rate of supplemental oxygen, the inhaled oxygen concentration will vary depending on the pressure setting, patient's breathing pattern and leak rate.

- Do not smoke if oxygen is in use.
- Contact a health care professional in the event chest discomfort, shortness of breath, stomach distention, belching, or severe headache are experienced.
- Avoid eating or drinking three hours before using this mask to minimize the risk of vomiting during sleep.
- Do not use this mask if you are taking medication that may cause vomiting.

**CAUTION:** At low CPAP pressures the flow through the exhalation port may be inadequate to clear all exhaled gas from the tubing and some rebreathing may occur.

- Must only be used by persons who have received adequate training from qualified medical personnel.
- Used in combination with other medical devices might alter the performance of the mask.

## CAUTION

USA federal law restricts this device to sale by or on the order of a physician. The use of bleach, chlorine, alcohol or aromatic-based solutions may cause hardening and reduce the life of the product.


## FEATURES

- One piece molded silicone rubber; small, light and comfortable.
- Clear mask allows a good view of mouth and nasal area.
- Low dead space results in low CO<sub>2</sub> build up.
- Dual swivel 360° rotating elbow and tubing connector.
- Headgear with Velcro fasteners provides maximum versatility and movement.

## CLEANING INSTRUCTIONS

- Detach headgear from mask.
- Handwash gently with mild soap and warm water (30°C), rinse the mask and components thoroughly. **NOTE - Headgear should only be hand washed with mild soap by gently rubbing in warm water at 30°C. No washing machine, no ironing.**
- Allow mask and components to air dry before reusing. Do not place in direct sunlight.
- Replace mask if there are cracks, tears or excessive discoloration.
- The D100 is dishwasher safe if placed in top rack.

<b>Size:</b> Assembled Mask (without Headgear): 161mm (H), 90mm (W), 69mm (D)	
<b>Weight:</b> Mask incl. Headgear: 122g; Headgear: 31g	
<b>Materials used:</b>	
Mask frame: Polycarbonate (PC)	Forehead Pad: Silicone
Headgear: Nylon + Spandex	Cushions: Silicone
The mask system does not contain latex materials. This product does not contain PVC or Phthalates.	

RESISTANCE (cmH <sub>2</sub> O)	50LPM	100LPM
D100F	0.24	0.62


MODEL	VOLUME (ML)
D100F-S	296.6
D100F-M	339.4
D100F-L	367.8

Storage relative humidity: < 85%  
Operation temp: -18°C - 50°C  
Storage temp: -40°C - 60°C


## Expected Service Life

- Mask Frame Components = 6 months
- Headgear = 3 months
- Mask Cushions = 3 months

## Please follow assembly instructions below:

Place headgear with Velcro side down on a flat surface. Detach 2 Velcro tabs on A1 & A2.

Identify headgear fastening positions on the mask.

Pull the upper headgear strap A1 through the slot B1 on the mask, making sure the Velcro faces outward.

Fasten the Velcro tab.

Pull the upper headgear strap A2 through the slot B2 on the mask, making sure the Velcro faces outward. Fasten the Velcro tab.

Place the headgear over your head and adjust both upper headgear straps to a comfortable tension.

Hook the plastic buckle A3 onto B3 and A4 onto B4 on mask frame.

Readjust both upper headgear straps A1 & A2 to a comfortable tension.

Adjust both lower headgear straps A3 & A4 to a comfortable tension.

Ensure the overhead strap is in proper position to avoid slippage.

- Important Note:** When removing the mask for cleaning, there is no need to unfasten the Velcro tabs. Simply unclip the lower buckles and remove the upper headgear straps from the slots in the mask frame. Replace after cleaning.
- To remove the mask and headgear, simply unbuckle both buckles (A3 & A4) from the mask frame and pull off the mask.
- To reuse the mask and headgear, simply place upper headgear over head and hook both buckles.
- This product does not require maintenance.

Manufactured for  
DeVilbiss Healthcare LLC  
100 DeVilbiss Drive  
Somerset, PA 15501-2125 • USA  
800-338-1988 • 814-443-4881

DeVilbiss Healthcare GmbH  
Kamenzer Straße 3  
68309 Mannheim GERMANY  
+49-621-178-98-230

Manufactured by HSINER Co., Ltd  
312, Jhongsan Rd., Shengang Dist  
Taichung City 429, Taiwan  
Tel: 886-4-25152480 • Fax: 886-4-25152482  
Email: sales@hsiner.com • Web Site: http://www.hsiner.com

EU Representative  
mdi Europa GmbH  
Langenhagener Str. 71  
D-30855 Langenhagen, Germany • Tel: +49 511 39089530  
Email: werner.sander@mdi-europa.com • internet: www.mdi-europa.com

CE  
0434

DeVilbiss® is a registered trademark of  
DeVilbiss Healthcare.  
© 2013 DeVilbiss Healthcare LLC. 08.13  
All Rights Reserved. A-D100F Rev. B

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DeVilbiss® is a registered trademark of DeVilbiss Healthcare.