



# Upeasy Power Lifting Seat

## User Guide

### Read all instructions before using the Upeasy Power Seat

#### IMPORTANT SAFETY INSTRUCTIONS

**DANGER** - To reduce the risk of electric shock:

Always unplug the Upeasy Power Lifting Seat from the electrical outlet before cleaning.

**WARNING** - To reduce the risk of burns, fire, electric shock, or injury to persons:

1. Unplug from outlet before putting on or taking off parts.
2. Close supervision is necessary when this product is used by, or near children.
3. Close supervision is necessary when this product is used by, or near pets or animals.
4. Use the Upeasy Power Lifting Seat only for its intended use as described in these instructions.
5. Do not use attachments not recommended by the manufacturer.
6. Never operate the Upeasy Power Lifting Seat if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water, return the Upeasy Power Lifting Seat to a service center for examination and repair.
7. Keep the cord away from heated surfaces.
8. Never drop or insert any foreign object into any opening.
9. Do not use outdoors.
10. Do not use in or near water.
11. To disconnect, turn all controls to the off position, then remove plug from outlet.
12. Do not unplug by pulling on cord. To unplug, grasp the plug, not the cord.
13. Unplug from outlet when not in use and before servicing or cleaning.
14. Use the Upeasy Power Lifting Seat only in armchairs or sofas with at least one armrest for optimum stability when sitting or rising. Upeasy Power Lifting Seat is not intended for use in rocking chairs, wheelchairs, office chairs, or any chair with wheels. Not for use in motor vehicles.
15. Keep fingers clear of the lifting mechanism when in use.
16. Do not allow direct contact with any type of heat source.
17. Not for use in bathrooms.
18. Consult your physician to ensure there are physical limitations or are other medical reasons that may contraindicate safe usage of the Upeasy Power Lifting Seat.

#### SAVE THESE INSTRUCTIONS

The Upeasy Power Lifting Seat is intended for household use only.

#### **Disclaimer**

*Instructions and recommendations for safe and proper use are outlined in this Users Guide. Carex Health Brands is not liable for property damage or personal injury arising out of unsafe or improper use of the Upeasy Power Lifting Seat. Carex Health Brands is not responsible for injuries and/or damage resulting from any person's failure to exercise good judgment and/or common sense. If you are having any difficulty operating your Upeasy Power Lifting Seat, or do not feel capable of operating it safely and comfortably, please contact the place of purchase for assistance.*



## Welcome Power Seat Users

Thank you for purchasing the Upeasy Power Lifting Seat. It has been engineered to provide you with the utmost in comfort and convenience. With proper care, you should enjoy years of trouble-free mobility assistance.

For instructions on correct assembly and proper use of your Upeasy Power Lifting Seat, please read this manual in its entirety. For more information on our lifting seat products, please visit our web site at [www.carex.com](http://www.carex.com).

## Safety Instructions

The following safety instructions and tips will help you operate your Upeasy Power Seat safely. Follow them carefully to avoid personal injury.

- As with any electrical device, do not use the Upeasy Power Seat in or near water.
- Do not operate if the power cord has been damaged.
- Do not grab the side of the lifting cushion when in use. Keep fingers clear of the lifting mechanism.
- Unplug your Upeasy Power Seat before performing any electrical inspection.
- Do not attempt to fix electrical problems yourself. This will void the warranty.
- Do not allow children or pets near the Upeasy Power Seat while it is in motion.
- Do not allow the fabric to come into direct contact with any type of heat source, as personal injury or fire damage may occur.
- Use only in chairs and sofas that have at least one available arm for optimum stability when you are sitting or rising.
- Not for use in rocking chairs, wheelchairs, office chairs, or any chair with wheels.
- Not for use in motor vehicles.
- Not for use in the bathroom.
- Consult your physician to ensure you do not have particular physical limitations or are taking medications which may contraindicate the safe usage of the Upeasy Power Seat.

## WARNINGS

To reduce the risk of fire or electric shock, do not attempt to remove the outer cover. There are no user-serviceable parts inside.

Repair should be done only by authorized service personnel. Note: If the power supply cord is damaged, it should be replaced by qualified personnel. Any other servicing, except as described in this user guide, is to be performed by an authorized service representative.

When not in use, the Upeasy Power Seat should be unplugged and/or the Power Lever(s) be removed. The Upeasy Power Seat may be stored in its original packaging.

## Indications for Use

Those who have lost the ability to get up unassisted from a seated position that may be due to the following:

- Muscular dystrophy
- Lou Gehrig's disease/Amyotrophic lateral sclerosis
- Post-operative conditions
- Degenerative joint disease
- Arthritis
- Parkinson's Disease

## User Benefits

- *Helps keep people active and independent.*
- *Provides safe lift assistance for accident prevention.*
- *Relieves stress and strain on joints.*
- *Provides 100% lift, allowing users to conserve muscular strength.*

## Specifications

CSA C/US CE Approved  
US Patent #6,702,383 and #5,316,370. Other US,  
European and International patents pending.

### Seat Dimensions

Length:	19.5"
Width:	17.1"
Weight:	12.3 lbs
Cover:	100% polyester washable


### Seat Height:

Lifted:	10"
Flat (when lowered):	3"

**Capacity:** For users up to 300 lbs

**Power Supply:** 120Vac, 60 Hz, 0.6 A

**Cord Length:** 10 ft



Please check to see that all parts are included in your package. Contact your dealer if any product parts are missing or damaged. Send in your completed Warranty Card today to ensure that you are registered.

## Assembly & Testing

For correct operation, and to ensure prevention of electrical shock, follow these directions for assembly.

1. Place your lifting cushion on an appropriate chair.
2. Insert the power lever(s) over the steel prongs on either side of the lifting cushion at the front. Angle the levers towards the front as shown in the photo below is operating properly.


## Before Using your Power Seat

1. Check package contents to ensure that all components are enclosed and in good condition (see picture).
2. Test the lifting seat to ensure that it is in working order.

### Package Contents

Your Upeasy Power Lifting Seat comes with the following components:

- Upeasy Power Lifting Seat (in the flat position) with power cord.
- 2 Detachable Power Levers.
- Washable cover.
- User Guide.
- Warranty Registration Card.



3. Plug the power cord into an electrical socket.
4. Test your Upeasy Power Seat before using it to ensure it is operating properly. There are three operating positions: up, down and stop. We recommend that you become familiar with these before using your lifting cushion.

Gently lift the lever up to raise the seat. Select the middle setting to stop. Press down to lower.

**Note:** Your Upeasy Power Seat will also automatically stop when it is at its full height, and when it is completely lowered, unless you move the power lever to the “stop” position.

## Using Your Power Seat

This section explains how to safely use your Upeasy Power Seat.

### Choosing Appropriate Chairs

The Upeasy Power Seat may be used on most chairs, recliners and sofas – but always next to an armrest. This is to ensure the utmost stability when using the Upeasy Power Seat. We also recommend that the chair or sofa you are using it on have a suitably stable seating surface.

### Setting Up For Use

When placing your Upeasy Power Seat on your chair or sofa, the lever(s) should be to the front and free from obstructions. Also, make sure the cord is not under the base of the seat to prevent it from damage or from becoming unplugged. Finally, ensure that the cord is out of the way so that no one can trip over it.

Your Upeasy Power Seat will add about two inches (2"/5cm) to the height of your chair, recliner or sofa when you are seated. In order to use your Upeasy Power Seat safely, you will need to keep both feet on the floor at all times.



There are three ways to keep your feet on the floor at all times during use (*Please try these in the order shown*):

- (1) Stop the Power Seat to control the height. Shift the Power Lever to the middle position to stop at a point where you are able to get out of the chair comfortably with both feet on the floor. When you sit, your Power Seat will be positioned at the perfect level for you.
- (2) If the seat cushion of the chair/sofa is removable, try using your Upeasy Power Seat without the cushion, or substitute a lower cushion.
- (3) Reduce the height of the chair or sofa by adjusting the legs so that your feet remain on the floor.

Now that you have your Upeasy Power Seat placed in your chair or sofa, gently shift the Power Lever and let the seat rise to a height that allows you to keep both of your feet on the floor.

You may stop it at the appropriate height by gently shifting the Power Lever to the center position or let the seat rise to its full height.

Now your Upeasy Power Seat is ready for use!

**Note:** If you have to hop up to get onto the Upeasy Power Lifting Seat, then you have raised it too high. Simply lower your power seat so you can sit down with both feet on the floor.

# Operating Instructions

When operating your Upeasy Power Seat, keep at least one hand on the arm of your chair or sofa to ensure stability.

You only need to use one of the power levers to raise and lower the seat.

## Lowering into a Seated Position

**Note:** Your Upeasy Power Seat should be in the raised position.

- Facing away from the seat, place your hand(s) on the arm(s) of the chair, or sofa.
- Sit on the raised lifting cushion – keeping both of your feet flat on the floor.
- Click one of the levers downward – past the “stop” position – and release.
- The motor will lower the lifting cushion and automatically shut off when it is fully lowered – unless you choose to stop it first.

**Note:** Please remember that it is neither necessary, nor advisable to keep pushing/pulling on the lever once the motor has started.



## Rising from a Seated Position

- Place your hand(s) on the arm(s) of the chair, or sofa and ensure both feet are on the floor.
- Click one of the levers upwards – past the “stop” position – and release.
- The motor will raise the lifting cushion and automatically shut off when it is fully raised – unless you choose to stop it first.



## Using the “Stop” position

You can use the “stop” (center) position on the power lever to stop your Upeasy Power Seat at any time.

You may want to use this:

- When rising – if you have reached the desired height.
- Anytime you feel it is necessary to maintain your stability.

## Care and Maintenance


In order to ensure the proper care and maintenance of your Upeasy Power Seat, please read these instructions carefully.

### Easy Care

Your Upeasy Power Seat comes with a machine washable, fabric cover. The surface of the foam is also coated with a water resistant covering for easy care.

For proper care, please follow these instructions:

- Wipe up any spills on cover immediately.
- Remove cover for washing or to clean the foam cushion surface.
- Machine-wash your fabric cover in **cold** water and **hang to dry**.
- Wipe down the foam cushion surface with a damp sponge. (We recommend using mild cleaners and detergents.)



### Inspection of Electrical Components

We recommend performing routine maintenance checks on your Upeasy Power Seat to ensure your seat gives you years of trouble-free operation. Visually inspect the power cord for any cracks or wear.

If a problem is discovered with any component of the electrical system, contact your authorized Upeasy Power Seat vendor immediately.

**NOTE:** Do not attempt to repair an electrical problem yourself. Attempting to fix electrical components may result in personal injury, and will void the warranty.

## Trouble Shooting

### What if my Upeasy Power Seat does not operate at all when I push/pull on the lever?

- First, ensure the lever is firmly installed. You should hear a click when it is inserted over the metal prong in the base. It should not come out when you tug gently on it.
- Check that the power cord is plugged firmly into a properly wired wall outlet.
- If the above measures fail, check your circuit breaker box as a final step before concluding there is a problem with the product.

### Where do I find the serial number on my lifting cushion?

The serial number is a 10 digit number on a narrow white sticker, located on the base of the seat.


### What if I am having difficulty setting up or using my Upeasy Power Seat?

First read all instructions carefully to ensure you have set up your lifting cushion correctly. If you are still having difficulty, or do not feel capable of operating this lifting cushion safely, please contact your authorized Upeasy Power Seat vendor for assistance.



**Correct Usage: Both Feet on the Floor**

Before sitting, ensure the Upeasy Power Lifting Seat is at the ideal height for you to sit while keeping both feet flat on the floor. Keep at least one hand on an armrest for support while lowering and raising the seat.



**Incorrect Usage: Cushion too High**

If your Upeasy Power Lifting Seat is raised too high, you will have difficulty sitting down without hopping onto the cushion. To maintain stability and reduce your risk of falling, lower the cushion to a height that allows you to keep both feet on the floor.



**Incorrect Usage: Sitting Lopsided**

If your Upeasy Power Lifting Seat is correctly positioned, you should sit down squarely on the seat and keep both feet on the floor. If you find yourself sliding sideways, the seat is lifted too high for you, and you could fall or damage the product.

Lower the cushion to the correct height.



**Incorrect Usage: Too Much Force on Power Lever**  
The Power Lever requires minimum force to operate. Placing excessive force on the switch is unnecessary and could result in the Power Lever breaking.



# Upeasy Power Lifting Seat

## User Guide

### Warranty Returns

If your Upeasy Power Lifting Seat fails to function as warranted, contact the original place of purchase for assistance. You will be required to provide date and proof of purchase.

We encourage all dealers to provide warranty service for Upeasy Power Lifting Seat customers so that you do not experience delays in receiving a replacement unit.

**Note:** No returns will be accepted without a Return Authorization Number (RA). Please have your dealer call 1-800-526-8051 for your RA Number.

We strongly recommend that you use the original packaging material and box for all Warranty returns. These items were designed to withstand the shipping process.

If the Upeasy Power Lifting Seat incurs any damage due to improper packaging, the cost for repair will be charged back to you.

### Warranty

Carex Health Brands warrants your Upeasy Power Lifting Seat against product failure for one year from the date of purchase. Damage due to neglect or misuse is not warranted, nor is this warranty transferable.

If your Upeasy Power Lifting Seat is defective in material or workmanship Carex Health Brands will, at its discretion, either repair or replace your Upeasy Power Lifting Seat throughout the duration of the year at no charge except for the current delivery charge.

This warranty is void if the product has been modified with unauthorized parts or attachments.

**To ensure that we have your Warranty information on file, complete and return your Warranty Registration card today.**