

Posey Alarm Sensors Compatibility Chart (Sorted by Application)

Alarm systems that can be powered by both AC and batteries include the KeepSafe Essential and KeepSafe Deluxe. The other alarms are all battery powered only.

Pull-cord alarms include the following: Personal Alarm, KeepSafe Essential or Sitter Select.

Motion sensor alarms include the Wireless Infrared Alarm or the PIR Sensor matched with KeepSafe, KeepSafe Deluxe, Sitter Select or Sitter II.

Pressure sensor alarms include the KeepSafe, KeepSafe Deluxe, Sitter Select and Sitter II.

Belt sensor alarms include the KeepSafe, KeepSafe Deluxe, Sitter Select and Sitter II.

	Item #	Applications	Description	Personal Alarm 8202L	Door Alarm 8205	KeepSafe 8350	KeepSafe Essential 8373 8373N	KeepSafe Deluxe 8374 7374N	Sitter Select 8361	Sitter II 8281	Wireless Infrared Alarm 8376
Bedside Alarms	POS8250L	Bedside , cord measures 120 inches	Exit Alarm Sensor Mat, 60 x 24 inches			✓		✓	✓	✓	
	POS8351	Bedside , door or room sensor	Passive Infrared (PIR) Sensor			✓		✓	✓	✓	
Bed Mattress Alarms	POS8285L	Bed mattress - Hospital or nursing home beds with mattresses over 38 - 42 pounds, cord measures 96 inches	Under-Mattress Sensor-Light			✓		✓	✓	✓	
	POS8285H	Bed mattress - Hospital or nursing home beds with mattresses over 42 pounds, cord measures 96 inches	Under-Mattress Sensor-Heavy			✓		✓	✓	✓	
	POS8292	Bed mattress - Hospital or nursing home beds, cord measures 96 inches	Under-Mattress Bed Pad Sensor			✓		✓	✓	✓	
	POS8307	Bed mattress sensor for longer term use	Over Mattress Sensor, 6-month			✓		✓	✓	✓	
	POS83073PK	Bed mattress sensor for longer term use, 3 pack	Over Mattress Sensor 3 Pack, 6-month			✓		✓	✓	✓	
	POS8283	Bed mattress - short-term admissions or fall assessment, single patient use, cord measures 96 inches	Bed Pad Alarm Sensor, 30-day, 13 x 31 inches			✓		✓	✓	✓	
	POS8291	Bed mattress - under-mattress sensor support	Bed Alarm Sensor Support, under-mattress								
	POS8292	Bed mattress - Under-mattress sensor	Bed Alarm Sensor, under-mattress			✓		✓	✓	✓	
Door Alarms	POS8351	Door , bedside or room sensor	Passive Infrared (PIR) Sensor			✓		✓	✓	✓	
	POS8250S	Door or room sensor, cord measures 120 inches	Exit Alarm Sensor Mat, 36 x 24 inches			✓		✓	✓	✓	
Room Alarms	POS8351	Room , door or bedside sensor	Passive Infrared (PIR) Sensor			✓		✓	✓	✓	
	POS8250S	Room or door sensor, cord measures 120 inches	Exit Alarm Sensor Mat, 36 x 24 inches			✓		✓	✓	✓	
Wheelchair Alarms	POS8360	Wheelchair for additional posture support, cord measures 24 inches	Chair Belt Sensor			✓		✓	✓	✓	
	POS8308	Wheelchair or other seating	Chair Pad Alarm Sensor, Square 13 x 13 inches			✓		✓	✓	✓	
	POS7223F	Wheelchair or other seating to reduce pressure	Gel Foam Alarm Sensor Cushion			✓		✓	✓	✓	
	POS8358	Wheelchair or other seating, easily transfers from chair to chair, D-ring attachment	EZ Clean Alarm Sensor Belt			✓		✓	✓	✓	
	POS8371	Wheelchair or other seating, easily transfers from chair to chair, D-ring attachment	Mobile Chair Belt Sensor			✓		✓	✓	✓	
	POS8372	Wheelchair or other seating, easily transfers from chair to chair, D-ring attachment	Hook & Loop Alarm Belt Sensor			✓		✓	✓	✓	
	POS8372L	Wheelchair or other seating, easily transfers from chair to chair, D-ring attachment	Hook & Loop Alarm Belt Sensor, Long			✓		✓	✓	✓	
	POS8309	Wheelchair or other seating, single patient use	Chair Pad Alarm Sensor, 30-day, Square 13 x 13 inches			✓		✓	✓	✓	
Bathroom Alarms	POS8332	Bathroom Alarm Sensor - Toilet Seat Sensor	Toilet Seat Alarm Sensor			✓		✓	✓	✓	

To order your customized system call [800-397-5899](tel:800-397-5899).