



# Posey® Bedfellow Positioning Roll

The Posey Bedfellow provides full body support in both lateral and supine positions. The Bedfellow is constructed of three inner pillows covered in a wipe clean vinyl. The outer cover is available in machine washable soft fleece or wipe clean vinyl.

- Three pillow construction allows the Bedfellow to bend and fit body contours
- Center pillow rolls out to provide a headrest or lumbar support
- Nylon ties easily secure the center pillow into position and prevent sliding from under the head
- Name tag allows Bedfellow to be personalized
- Front button secures the Bedfellow around the patient/resident and prevents the Bedfellow from falling out of the chair
- Side support protects ribs and positions the patient in the chair
- Helps prevent leaning or falling to the side
- Reduces use of siderails

#6306 Bedfellow, fleece cover

#6306SC Bedfellow, vinyl cover

#6308 Bedfellow replacement cover, fleece

#6308SC Bedfellow replacement cover, vinyl



Center pillow can be rolled out or fully detached for additional positioning options.



Center pillow provides a headrest while the two end pillows provide lateral positioning support.



70" length provides full body support.



Center pillow rolls out to provide lumbar support and prevent rolling.

