

# Tubigrip® and TSSB®

## Sizing guide

### Tubigrip

Size	Low pressure circumference	Medium pressure circumference	High pressure circumference
A	–	10–12.5 cm/3.9–4.9 in	12.5–15 cm/4.9–5.9 in
B	10–12.5 cm/3.9–4.9 in	12.5–15 cm/4.9–5.9 in	15–24.5 cm/5.9–9.6 in
C	13.5–15 cm/5.3–5.9 in	15–24.5 cm/5.9–9.6 in	24.5–35.5 cm/9.6–14 in
D	15–24.5 cm/5.9–9.6 in	24.5–35.5 cm/9.6–14 in	35.5–45 cm/14–17.7 in
E	24.5–35.5 cm/9.6–14 in	35.5–45 cm/14–17.7 in	45–50.4 cm/17.7–19.8 in
F	35.5–45 cm/14–17.7 in	45–50.4 cm/17.7–19.8 in	50.4–60.7 cm/19.8–23.9 in
G	45–50.4 cm/17.7–19.8 in	50.4–60.7 cm/19.8–23.9 in	60.7–70.3 cm/23.9–27.7 in
J	60.7–70.3 cm/23.9–27.7 in	70.3–75.5 cm/27.7–29.7 in	–
K	70.3–75.5 cm/27.7–29.7 in	–	–

Low pressure = 5–10 mmHg

Medium pressure = 10–20 mmHg

High pressure = 20–30 mmHg

### Tubigrip TSSB

Size	SMALL	MEDIUM	LARGE
<b>Circumference</b>	32–35.5 cm/12.6–14 in	35.5–38.5 cm/14–15.2 in	38.5–41.5 cm/15.2–16.3 in

For more information call toll free: 1-800-882-4582.

Tubigrip has been designed to provide tissue support, not graduated therapeutic compression. If graduated therapeutic compression is required, please consult your physician. The information provided herein is not to be construed as the practice of medicine or substituted for the independent medical judgment of a patient's treating physician. This information, including but not limited to suggestions for product sizing, product selection and suggested use is based on generalizations and does not consider the unique characteristics of an individual's condition. Each patient's physician shall remain solely responsible for assessing the appropriate treatment and managing treatment of the condition. For additional information, please refer to the applicable product insert or contact Mölnlycke Health Care at 1-800-882-4582.

