


# Extremity Pump® Systems

7000 / 7500

## PNEUMATIC SLEEVE SIZING CHART

This sizing chart reflects anatomical measurements that **must** be taken to ensure proper fitting of arm and leg sleeves (both with and without inserts). Simply take measurements as shown by the figure to the right, then determine which pump system (System 7000 Single-chamber Intermittent Pump or System 7500 Three-chamber Sequential Pump) will be used.



### UPPER EXTREMITY PNEUMATIC SLEEVES FOR INTERMITTENT PUMP SYSTEM 7000

ARM CIRCUMFERENCE AT AXILLA							
ITEM #	DESCRIPTION	ARM LENGTH	NO INSERT	W/NARROW INSERT	ITEM #	W/WIDE INSERT	ITEM #
147010	Hand & Wrist	N/A	N/A	N/A	N/A	N/A	N/A
147011	Large Hand & Wrist	N/A	N/A	N/A	N/A	N/A	N/A
147012	Short Arm	17" or shorter	Up to 14"	N/A	N/A	N/A	N/A
147013	Medium Arm	26" or shorter	Up to 16½"	Up to 20½"	147020	N/A	N/A
147014	Long Arm	30" or shorter	Up to 17½"	Up to 21½"	147021	Up to 24½"	147022

### UPPER EXTREMITY PNEUMATIC SLEEVE FOR SEQUENTIAL PUMP SYSTEM 7500

ARM CIRCUMFERENCE AT AXILLA							
ITEM #	DESCRIPTION	ARM LENGTH	NO INSERT	W/NARROW INSERT	ITEM #	W/WIDE INSERT	ITEM #
147514	Medium Arm	30" or shorter	Up to 21½"	N/A	N/A	N/A	N/A

### LOWER EXTREMITY PNEUMATIC SLEEVES FOR INTERMITTENT PUMP SYSTEM 7000

ITEM #	DESCRIPTION	LEG LENGTH	THIGH CIRCUMFERENCE			ANKLE CIRCUMFERENCE				
			NO INSERT	W/NARROW INSERT	ITEM #	W/WIDE INSERT	ITEM #	NO INSERT	W/NARROW INSERT	W/WIDE INSERT
147030	Foot & Ankle	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
147032	Youth Leg	17" or longer	Up to 18½"	N/A	N/A	N/A	N/A	Up to 17"	N/A	N/A
147034	X-Short Leg	24" or longer	Up to 20½"	Up to 24½"	147020	N/A	N/A	Up to 17"	Up to 19½"	N/A
147035	Short Leg	27" or longer	Up to 21¼"	Up to 25¼"	147021	Up to 28¼"	147022	Up to 17"	Up to 19½"	Up to 21"
147036	Medium Leg	30" or longer	Up to 22¼"	Up to 26¼"	147021	Up to 29¼"	147022	Up to 17"	Up to 19½"	Up to 21"
147037	Long Leg	33" or longer	Up to 23"	Up to 27"	147023	Up to 30"	147024	Up to 17"	Up to 19½"	Up to 21"
147038	X-Long Leg	35" or longer	Up to 23½"	Up to 27½"	147023	Up to 30½"	147024	Up to 17"	Up to 19½"	Up to 21"

### LOWER EXTREMITY PNEUMATIC SLEEVES FOR SEQUENTIAL PUMP SYSTEM 7500

ITEM #	DESCRIPTION	LEG LENGTH	THIGH CIRCUMFERENCE			ANKLE CIRCUMFERENCE				
			NO INSERT	W/NARROW INSERT	ITEM #	W/WIDE INSERT	ITEM #	NO INSERT	W/NARROW INSERT	W/WIDE INSERT
147535	Medium Leg	28½" or longer	Up to 25½"	N/A	N/A	N/A	N/A	Up to 19½"	N/A	N/A
147537	Long Leg	33½" or longer	Up to 27"	N/A	N/A	N/A	N/A	Up to 19½"	N/A	N/A



KCI Licensing, Inc., P.O. Box 659508, San Antonio, Texas 78265-9508. 1-800-275-4524.  
All trademarks, registered trademarks and service marks designated herein are property of KCI and its affiliates.  
©2004 KCI USA, Inc. All rights reserved.

Lit. No. 55-A-006 Rev. 4/2004