

Features and Benefits

Patient Comfort and Safety

Sheer Ankle Crest Prevents excess fabric abrading against the skin.

Heel Pocket More room and improved comfort.

Comfort Panel Accommodates changes in limb size and shape.

Top-View Inspection Port Better patient comfort when ambulating.

Rapid Moisture Transfer Moves fluids and body oils from skin to surface of the fabric quickly to speed drying; launders better.

Two-Way Stretch Fabric Allows stocking to conform to shape of limb.

Premium Elastic Top Keeps stocking in place without causing tourniquet effect.

Outward Facing Seams Reduces skin irritation and/or depressions as limb size changes.

Latex-Free Safe for latex-sensitive skin.

Clinician Focused

Sheer Ankle Crest Prevents excess fabric abrading against the skin.

Heel Pocket Eases donning and doffing while assisting proper placement.

Comfort Panel Allows more stretch to position monitors or check incision sights.

Color-Coded Top-View Inspection Port Conveniently access toes and top of foot and identify size.

Alignment Indicator Visible guide to prevent twisting during application.

Available Sizes and Styles

19 Sizes To Choose From Knee and thigh length.

For more information, visit www.DVTSTAR.com

C.A.R.E. Anti-Embolism Stockings

SHORT	CATALOG NUMBER		STYLE	SIZE	COLOR	CASE QTY
	REGULAR	LONG				
	453-01	458-01	Knee	Small	Yellow	12 PR
	453-02	458-02	Knee	Medium	White	12 PR
	453-03	458-03	Knee	Large	Blue	12 PR
	453-04	458-04	Knee	X-Large	Green	12 PR
463-01	463-02	463-03	Thigh	Small	Yellow	12 PR
468-01	468-02	468-03	Thigh	Medium	White	12 PR
473-01	473-02	473-03	Thigh	Large	Blue	12 PR
	473-04	473-05	Thigh	X-Large	Green	12 PR



C.A.R.E.®
Control and Resist Embolisms



Better care starts
with safety and comfort.



All items featured in this brochure are free of natural rubber latex

Albahealth is a subsidiary of Encompass Group LLC

800-284-4540 | www.encompassgroup.net

© 2014 Encompass Group, LLC
All Rights Reserved
Lit. No. AC10180 REV 06/14