











Tool-Free Adjustable-Length Home Style Bed Rail, Brown-Vein Finish

Standard Telescoping Full-Length Side Rail, Brown-Vein Finish

| Item # | UOM | Item # | UOM |
|---------|--------|----------|--------|
| 16500BV | 1pr/cs | 15001ABV | 1pr/cs |

- Rail length can be adjusted without tools to 37", 43" or 57" long. (Figure A - C)
- Crossbars are easy to install and remove. Adjusts between 36" and 72".
- Crossbars come with extensions that will accommodate a KING SIZE BED.
- Adjustment button allows for 3 different length adjustments. (Figure D)
- Fits all home-style beds.
- Easily adjusts up or down with spring-loaded release.
- Attractive brown-vein finish.
- Constructed of 1" steel.
- Fits all home-style beds.

- Attractive brown-vein finish. (Figure A B)
- Constructed of 1" steel.
- Telescopes with movement of bed.
- Easily adjusts up or down with spring-loaded release.
- Spring-loaded crossbars are easy to install and remove.
 (Figure C)

| specifications | pecifications | | | | | |
|----------------|---|----------|--|--|--|--|
| Warranty | Limited | Lifetime | | | | |
| accessories | | | | | | |
| 16500CBV | Cross Bars with King Size Bed Extensions for 16500BV (Comes with 16500BV) | 1pr/cs | | | | |

| specifica | specifications | | | | |
|-------------|---------------------|------------------|--|--|--|
| Warranty | | Limited Lifetime | | | |
| accessories | | | | | |
| 15050BV | Universal Crossbars | 1 pr/cs | | | |