

Lift Scale

- Digital scale easily attaches to Lumex LF1050 and LF1090 Patient Lifts
- Accurately weighs patients during daily lifts
- Combines the latest microprocessor and load cell technology to provide accurate and repeatable weight data
- Easy 4-button operation
- Weight is displayed for 5 minutes
- Auto power-off prolongs battery life
- Low battery indicator
- Measures pounds or kilograms
- Maximum weight capacity 600 lb
- Comes complete with attachment hardware
- 2 year limited warranty

DSC250 1 ea

