





# EZ POOL LIFT

\*MEASUREMENT TAKEN WITH THE LIFT AT FULL LOAD\*



**PRODUCT WEIGHT: 65lbs [30kg]**  
**MAX WEIGHT CAPACITY: 400lbs [181kg]**  
**MIN ANCHOR SETBACK: 12"**  
**MAX ANCHOR SETBACK: 38"**  
**MAX WATER DRAFT: \*18"**



**CONTINUOUS  
360° ROTATION**



**Aqua Creek Products**  
 The Total Solution for Pool and Spa Access

9889 GARRYMORE LANE  
 MISSOULA, MT 59808  
 (406) 549-0769  
 FAX (406) 549-2602

NOTE:  
 \*SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE

PRODUCT:	<b>EZ LIFT</b>	
PART #:	<b>F-03EZPL</b>	
SCALE:	<b>1:24</b>	DATE: <b>2/16/16</b>
DRAWN BY:	<b>KR</b>	REVISION: <b>A</b>