




*Both resident and caregiver benefit from the new Deluxe Assist Handle*

## Enhance Resident Safety with NEW Deluxe Assist Handle

- New shrouded design minimizes limb entrapment and entanglement while featuring two locking positions for additional safety
- Meets the current IEC 60601-2-38 requirements in all positions. Compliant with dimensions for zones 1, 2, 3 and 4 of the March 2006 FDA Entrapment Guidelines. Also addresses areas of concern identified as zones 5, 6, and 7 in the March 2006 FDA Entrapment Guidelines
- Mounted at an ergonomic height to best assist in safe ingress/egress while promoting independent repositioning
- Features ergonomic and sturdy yet soft touch handholds for improved resident grip and comfort
- Easy installation and quick removal without any tools needed, minimizes time demands on staff
- Optional integrated pendant holder on both sides making pendant stowage and retrieval simple for both residents and caregivers
- For use with EasyCare® 5, 7, 9, UltraCare® and UltraCare XT bed models



*Industry's best assist handle for resident egress/ingress promoting safety and comfort. Shown in two locking positions.*



*Deluxe assist handle stays with the bed frame and resident.*



*Easily removed for storage.*



*Deluxe assist handle shown with both pendant holder styles.*



**Enhance safety and comfort for your residents and caregivers while promoting independence.**

**Order the Deluxe Assist Handle today! Call 800.826.0270.**



**Joerns  
Healthcare**

A Tradition of Caring • [www.joerns.com](http://www.joerns.com)

6180078 RevB • 09-0337