

## Deluxe Tilt Overbed Table

- Chrome plated and welded tubular steel column and base make the table durable and stable.
- 1 ½" hooded, ball-bearing casters provide easy maneuverability across all floor surfaces.
- Table top adjusts 13" in infinite positions and is constructed of wood-grained oak Melamine laminate on particle-board core.
- Raised T-molded plastic edge.
- Chrome-plated "H" base and column for aesthetics and durability.
- 2 locking casters
- Designed to meet the requirements of HCPCS code: E0274.
- Supports up to 30 lb of dispersed weight.
- Warranty is Three years

GF8905-1A      1 ea



Item	Table Height	Table Width	Table Depth	Base Width	Base Depth	Weight	Dispersed Weight Capacity
GF8905-1A	30" - 43.8"	30"	15.2"	26.2	16.5	21 lb	30 lb