

# Patient Mobility



MDS806300D

## Excel® 2000 Wheelchairs

These wheelchairs have padded and upholstered armrests and calf pads.

- » Carbon steel frame with rust- and chip-resistant chrome plating
- » Threaded seat-screw inserts reduce stripping from upholstery changes
- » Hemi height-adjustable axles raise and lower seat 2-inches (5.08 cm)
- » Chart pocket on back
- » Easy-to-clean vinyl upholstery
- » Weighs: 38-lbs. (17 kg)
- » Weight capacity: 300-lbs. (135 kg)
- » Warranty: Lifetime on the frame; one year on all parts

Medline has a complete line of wheelchairs to meet every need. Ask your Medline Education Sales Representative for additional information.

Item No.	Description	Pkg
MDS806400	Desk-Length Removable Arm, Detachable Footrest, Seat 20 x 16" (50.8 x 40.6 cm)	1/ea
MDS806250D	Desk-Length Removable Arm, Detachable Footrest, Seat 18 x 16" (45.7 x 40.6 cm)	1/ea
MDS806300D	Desk-Length Removable Arm, Detachable Elev. Legrest, Seat 18 x 16" (45.7 x 40.6 cm)	1/ea

Full-length arms available by adding "FLA" to the end of the item number.  
Navy and Ruby models available; simply replace the "D" with "NAVY" or "RBY".

Item No.	Description	Pkg
<b>Accessories</b>		
WCA6100	Oxygen Holder	1/ea
MDS85196	Anti-Fold Device	1/ea
MDS85183	I.V. Pole	1/ea
MDS85197A	Overhead Anti-Theft Device	1/ea
MDS85190	Oxygen/I.V. Combo	1/ea
MDS85189	Anti-Tip Device	1/ea
WCA806991	Tinnerman Legrest Locks	1/ea



MDS806700

### Extra-Wide Wheelchairs

#### Features:

- » Threaded seat-screw inserts minimize stripping from upholstery changes
- » Padded and upholstered armrests and calf pads are standard
- » Carbon steel frame with rust- and chip-resistant chrome plating
- » Easy-to-clean navy vinyl upholstery
- » Chart pocket on back
- » 500-lb. (226.8 kg) weight capacity

Item No.	Description	Pkg
MDS806700	Desk-Length Removable Arm, Detachable Footrest, 20 x 18" (50.8 x 45.7 cm)	1/ea

\*Full-length arms available, just add "FLA" to the item number.

Item No.	Description	Pkg
<b>Accessories</b>		
MDS85196	Anti-Fold Device	1/ea
MDS85197A	Overhead Anti-Theft Device for 18–22" (45.7–55.9 cm) Models	1/ea
MDS85183	I.V. Pole	1/ea
WCA6100	Oxygen Holder	1/ea
MDS85190	Oxygen/I.V. Combo	1/ea
MDS85189	Anti-Tip Device	1/ea
MDS85192	22" (55.9 cm) Pushbar	1/ea
WCA806991	Tinnerman Legrest Locks	1/ea
MDS851965SH	Anti-Fold, 24" (61 cm) Models	1/ea