



MDS809750

Shuttle Bariatric Wheelchair

- » Our strongest and widest wheelchair
- » Double cross-brace for maximum support and durability
- » Arms and legrests are removable to make exiting easy
- » Comfortable nylon upholstery and padded legrests
- » Desk length arms
- » Hemi-height adjustable with a triple-axle design
- » 850-lb. (385.6 kg) weight capacity
- » Warranty: 12 months on parts, lifetime on the frame and welding

Item No.	Leg Rests/Footrests	Seat Width
MDS80650	Elevating Leg Rests	24" (61 cm)
MDS809750	Elevating Leg Rests	26" (66 cm)
MDS809850	Elevating Leg Rests	28" (71.1 cm)
MDS809950	Elevating Leg Rests	30" (76.2 cm)



MDS808200BAR

Freedom Plus Transport Wheelchair

- » Weighs just 33-lb. (15 kg)
- » Back folds down for easy storage and transport
- » Handbrakes and a seat belt for safety and control
- » Easy-to-clean vinyl upholstery
- » Anti-tippers
- » Push bar with foam grips, cup holder and a chart pocket on the back
- » Full-length permanent arms
- » 400-lb. (181 kg) weight capacity
- » Warranty: 12 months on parts, lifetime on the frame and welding

Item No.	Leg Rests/Footrests	Seat Width
MDS808200BAR	Elevating Leg Rests	22" (55.9 cm)

Our Warranty

We stand behind all our wheelchairs with the industry's best warranty: 12 months on parts and lifetime on the frame and welding.