



***STANDING STRONG
AND READY.***

IV Poles
Heavy-Duty, Smooth-Rolling, and Durable

Designed with a genuine understanding of real-world use.

All Medline IV poles resist tipping and spinning, are easy to maintain and durably constructed to stay in service longer.

Infusion Stand with CuVerro® Copper Alloy

The center pole and optional pump hoop and handle are made with CuVerro copper alloy, which is registered with the EPA to continuously kill the bacteria that cause infection within two hours of exposure.¹

- » Infusion pump holder for up to six infusion pumps can be adjusted rotationally and vertically
- » Bag drop prevention device for controlled descent of the pole
- » Pole height adjustment: 69"-124" (175.3-315 cm)
- » 27-lb. (12.2 kg) low-profile powder-coated steel base
- » Heavy-duty 3" (7.6 cm) casters


Item No.	Description	Pkg.
PEDP10836	Infusion Stand, 5-Leg	ea
PEDP1083F	Infusion Pump Hoop, 29"H x 19"W (73.7 x 48.3 cm)	ea
PEDP3303	Patient Handle	ea



Infusion Pump Hoop (PEDP1083F)



Patient Handle (PEDP3303)



PEDP10836

Heavy-Duty IV Pole

The top-quality 3" (7.6 cm) quick-release ball-bearing casters are built to last, keeping the pole in service longer.

- » Rubber bumpers for scuff prevention
- » Heavy-duty six-leg 24" (61 cm) dia. base with a low center of gravity—recommended for triple-channel pumps
- » D-tube insert for stability and to help prevent the base from spinning
- » Stainless steel telescoping pole, adjusts 73"–99.5" (185.4–253 cm)
- » Removable top with four ram-horn style hooks—easily replaced with an eight-hook rake (sold separately)

Item No.	Description	Base Color	Pkg.
MDS80600	Heavy-Duty IV Poles	Blue	2/cs
MDS80600RED	Heavy-Duty IV Poles	Red	2/cs
MDS80600GRY	Heavy-Duty IV Poles	Gray	2/cs
MDS08600BLK	Heavy-Duty IV Poles	Black	2/cs



Heavy-Duty IV Stand Accessories


Item No.	Description	Pkg.
MDS80600UBH	Urine Bag Hook	10/cs
MDS80600SAT	Push Handle with Tray	6/cs
MDS80600B	Basket	2/cs
MDS806000TH	Oxygen Tank Holder	2/cs
MDS80600C	Replacement Casters	2/bx
MDS80600R	8-Hook Rake Top	1/ea
MDSIVCLIPZ	Barrier Drape Clip	10/bg



Oxygen Tank Holder
(MDS806000TH)



8-Hook Rake Top
(MDR80600R)





Deluxe IV Pole

Corrosion-resistant aluminum provides years of dependable service.

- » Smooth, quiet operation
- » Microphone-style knob for easy adjustments
- » Telescoping pole, adjusts 48.5–82" (123.2–208.28 cm)
- » Five-leg, 24" dia. (61 cm), cast aluminum base
- » 2" (5.1 cm) heavy-duty casters
- » Four ram-horn style hooks

Item No.	Description	Pkg.
MDS80494	Deluxe IV Pole	2/cs



5-Leg IV Pole, Low Center of Gravity

- » Telescoping pole, adjusts 48–80" (122–203 cm)
- » Five-leg, 21.5" (61 cm) dia., cast iron base
- » 2" (5.1 cm) heavy-duty casters
- » Four ram-horn style hooks

Item No.	Description	Pkg.
MDS80500	5-Leg IV Pole	2/cs

For more information, contact your Medline representative, call **1-800-MEDLINE** or visit our Capital Quote Tool at www.medlinecapitalquote.com today!



Medline Industries, Inc. One Medline Place, Mundelein, IL 60060

Medline United States
 1-800-MEDLINE (633-5463)
medline.com | info@medline.com

Medline Canada
 1-800-396-6996
medline.ca | canada@medline.com

Medline México
 01-800-831-0898
medlinemexico.com | mexico@medline.com

FOLLOW US 