



HCPCS  
E0143/E0156



## Rollator with 6" Wheels

Item #	Color	HCPCS	UOM
R800BL	Blue	E0143 / E0156	1/ea
R800RD	Red	E0143 / E0156	1/ea

- Comfortable padded seat.
- Pouch under seat for added privacy and security of personal items.
- Removable, hinged, padded backrest can be folded up or down.
- 6" Black non-marring casters with soft grip tires are ideal for indoor and outdoor use.
- Easy-to-use deluxe loop locks. (Figure A)
- Brakes with serrated edges provide firm hold.

specifications	
Handle Height	32" - 37"
Unit Dimensions	24.5" (W) x 27.5" (L)
Seat Dimensions	12" (W) x 20" (H) x 14" (D)
Weight	16.7 lbs.
Carton Shipping Weight	17 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime