

## PVC Shower Chairs with Commode Opening and Gurney's

Comes standard with locking casters making transportation easier and can eliminate an additional transfer. PVC bath safety equipment is strong, lightweight and easy to clean.

- Anti-slip grips provide secure grip in a slippery environment
- Ideal for use in bathing since PVC won't rust, stain, or fade color
- Products are long lasting with 5 year warranty on the frame
- PVC chairs are warm to the touch for added patient comfort
- Reclining chair has adjustable headrest and footrest for multiple bathing positions



### PVC Shower Chairs

Item No.	Description	Overall Dimensions	Internal	Wt Cap.	Pkg
PVCM1183	Standard Shower Chair	22"W x 18"D x 40"H (55.9 x 45.7 x 101 cm)	18" (45.7 cm)	300lbs (136 kg)	Each
PVCM1223	Wide Seat Shower Chair	26"W x 18"D x 40" H (66 x 45.7 x 101 cm)	22" (55.9 cm)	375lbs (170 kg)	Each
PVCM193	Reclining Shower Chair	24"W x 62"D x 48"H (61 x 157 x 122 cm)	20" (50.8 cm)	325lbs (147 kg)	Each

### PVC Shower Gurney

Item No.	Description	Weight Capacity Lbs (kg)	Overall Length Inch (cm)	Overall Width Inch (cm)	Overall Height Inch (cm)	Height with Rails Up Inch (cm)	Pkg
PVCM910	Standard	450 (204)	79 (201)	30 (76)	33.5 (85)	40 (102)	Each
PVCM910B	Bariatric	600 (272)	79 (201)	32 (81)	31.5 (80)	40 (102)	Each

## Safe-T-Pole

This floor-to-ceiling pole safely assists with standing, sitting, climbing stairs or transferring.

- Pole diameter of 1.5" (4cm)
- Height adjustment 7'8- 8.25 (234-260cm)
- Weight: 14lbs (6.3 kg)
- Order both items to get complete unit

Item No.	Description	Diameter	Pkg
G98120B	Base for the Safe-T-Pole	1.5" (4cm)	Each
G98120P	Pole portion for the Safe-T-Pole	N/A	Each



## FALL STATISTICS

*In 2010, 2.3 million nonfatal injuries among older adults were treated in emergency departments and more than 662,000 of these patients were hospitalized.*

Centers for Disease Control and Prevention.  
Falls Among Older Adults: An Overview.  
Available at: <http://www.cdc.gov/home-andrecreational-safety/falls/adultfalls.html>.  
Accessed January 8, 2014.