



The Monaghan Advantage

- ❖ Proprietary anti-static polymer ensures consistent and optimum medication delivery
- ❖ Comfort sensitive design for special needs patients
- ❖ Easy actuation with minimal inspiratory effort
- ❖ Minimizes static charge and provides consistent aerosol availability
- ❖ Compatible with most commonly prescribed MDIs
- ❖ Not made or manufactured with BPA, natural rubber latex, DEHP or lead
- ❖ Designed and manufactured in North America

AEROTRACH PLUS® aVHC

The AEROTRACH PLUS® anti-static valved holding chamber (aVHC) sets the standard for special-need tracheostomy patients. This cost-effective system quickly delivers medicine where it's needed most – straight to your patient's lungs – with the least discomfort to the patient.

The AeroTrach Plus® aVHC anti-static polymer chamber minimizes static charge and improves aerosol availability. The see-through chamber allows caregivers to visualize the medication dose in the chamber.

THE FACTS

Monaghan is known for extensive research and development of aerosol delivery systems.

The AeroTrach Plus® aVHC aerosol delivery is quick and easy.

Consistent concentrations are delivered on inspiration only.

Exceptional particle suspension time delivers a cost-effective treatment.

Typical 15 mm adapters are designed to connect firmly in place. The connection fits tightly enough for an exceptional aerosol delivery, without excessive torque when connecting or disconnecting the device.

Monaghan's anti-static polymer chamber minimizes static charge which improves aerosol availability.

EFFICIENT,
COMFORTABLE
DELIVERY FOR
TRACHEOSTOMY
PATIENTS



518.561.7330 phone
800.833.9653 toll free
518.561.5088 fax


5 Latour Avenue, Suite 1600
Plattsburgh, New York 12901
monaghanmed.com

THE FACTS (continued)

With a protective design for low resistance and superior responsiveness, even the smallest inspiratory effort will be rewarded. Additionally, the valve features increased durability.

Exhalation valve is designed for patients to experience low resistance to exhaled flow. Additional protection from valve tampering or damage is built in.

The MDI adapter is made of a durable polymer material which increases product durability.



ORDERING INFORMATION

PART #	NAME/DESCRIPTION	UNITS PER CASE
52510	AEROTRACH PLUS® anti-static Valved Holding Chamber (aVHC)	cs/10



Contact your Monaghan Sales Territory Manager for more information or a demonstration.

