

URINARY LEG BAGS

AMSure® Urinary Leg Bags provide discrete and comfortable drainage. The bag is available in two sizes with pre-attached leg straps, drainage port and an anti-reflux device. Push-pull and twist-turn drainage port design options available for patient dexterity.

- Soft and breathable pre-attached leg straps provide patient comfort
- Anti-reflux device helps improve patient safety by reducing urine backflow
- Leak-proof seams ensure fluid containment
- Sterile fluid pathway, latex-free and BPA-free



ITEM CODE	600 mL LEG BAG	900 mL LEG BAG	PUSH-PULL DRAIN PORT	TWIST-TURN DRAIN PORT	18 IN. (45.72 CM) EXTENSION TUBE WITH ADAPTER AND PROTECTIVE CAP	QUANTITY/UNIT
AS306N	•		•			48/Case
AS307N		•	•			48/Case
AS308N	•			•		48/Case
AS309N		•		•		48/Case
LBE306N	•		•		•	48/Case
LBE307N		•	•		•	48/Case
LBE308N	•			•	•	48/Case
LBE309N		•		•	•	48/Case

EXTENSION TUBE FOR LEG BAGS

The 18 in. leg bag extension tubing is essential for comfort and proper placement of the leg bag. The leg bag can be positioned either on the thigh or calf. The extension tubing can be cut to the length the patient requires. It includes an adapter and protective cap.

- Kink-resistant tubing safeguards against fluid buildup
- Flexible tubing permits leg bag to be positioned on the thigh or calf
- Sterile, latex-free and BPA-free

ITEM CODE	DESCRIPTION	QUANTITY/UNIT
AS311	18 in. (45.72 cm) Extension Tube	50/Case



PEDIATRIC URINE COLLECTOR

AMSure® Pediatric Urine Collectors are hypoallergenic and provide graduation marks designed for convenient and simple urine collection.

- Adhesive surface folds to form a spill-proof specimen transport container
- Hypoallergenic feature makes it safe to use with most pediatric patients
- Visible graduation marks designed for ease of use
- Single-use disposable
- Sterile, latex-free, DEHP-free and BPA-free

ITEM CODE	DESCRIPTION	QUANTITY/UNIT
AS409	200 mL, Pediatric Urine Bag	50/Box

