



Aqua Creek Products

The Leaders of Recreation, Fitness and Ability

AquaJogg

Strengthen Your Stride!

The AquaJogg combines the benefits of water buoyancy with the ability to run or walk. Users can rehabilitate lower extremities while maintaining a sense of stability and security. The AquaJogg allows the individual to build muscle strength and endurance at their own pace.

The AquaJogg is designed for intensive use in group class or individual home exercise. It can be used barefoot. Its innovative folding system makes it compact for storage. Anti-scratch roller wheels allow easy transfers into the water without marking the floor.

FEATURES:

- Adjustable resistance
- Barefoot use
- Easily folds for storage
- Excellent stability
- Lightweight, easy to move
- Scratch resistant wheels
- Two year warranty on metal frame, six months on all other components



Item# F-WXAQJG

SPECIFICATIONS:

- 100% stainless steel frame
- Anti-corrosion treatment
- Recommended pool depth:
3 ft - 5 ft (0.9 m - 1.5 m)
- Weight: 58 lb (26 kg)
- 10.2° incline
- Open dimensions: L 53" x W 27" x H 51"
(135 cm x 67 cm x 128 cm)
- Folded dimensions: L 25" x W 22" x H 51"
(62 cm x 56 cm x 128 cm)

PROUDLY SOLD BY:

AquaJogg