



Summer 2021
INCONTINENCE CARE GUIDE





Tranquility:

Having freedom from unwanted distractions, allowing you to focus on the things that matter most – a calm and confident space to live well.

Our products help you manage even the most challenging incontinence care needs with higher-performing absorbent protection – so you can enjoy the comfort, confidence and freedom to live well and focus on the things that matter most. Inspired by personal experiences with loved ones who deserved better incontinence solutions, our family-owned company relentlessly perfects the most absorbent products – **delivering Real Life Protection for your real, everyday demands.**

Why Tranquility?

A Personal Focus

We offer a breadth of products that uniquely support specific medical conditions, lifestyle needs and budgets, so you can find a product that feels like it was made just for you.



Compassion-Driven Innovation

Through rigorous research and testing, we never stop refining better ways to support healthier, more fulfilled lives and higher quality care.



Commitment to Care

Family-minded to the core, we are committed to your needs, ensuring you may live with the confidence, comfort and continual support you deserve.



Catalog Guide

OVERVIEW

Tranquility® Overview	2
Real Life Protection	4
Real Life Stories	5
Product Brands	6
Product Guide.....	7
Features and Benefits.....	8

PRODUCTS

Briefs.....	9
Bariatric Briefs.....	12
Protective Underwear	15
Pads & Liners	19
Boosters.....	22
Underpads	25
Extra Care	28
Disposable Swimwear	30

REFERENCE

Sizing Guide.....	33
Education & Programs	34
Tranquility Essential Cross-Reference	36
Quick Order Guide.....	38



Tranquility® products absorb and contain their stated fluid capacity, no matter how much pressure is applied – whether you're sleeping, sitting or on-the-move. **We guarantee it!** You can rely on our clinically proven design to protect you during your real, everyday activities, so you can live freely and confidently, with fewer worries.

We Hold Cups, Not Drops

We get it, using drops to show absorbency can be confusing.

Tranquility values transparency. This is why we state our product's full retention capacity on every package, in cups.

No more misleading raindrops. No more second-guessing how much each product can hold.

We make it easy for you to understand which product you need and feel confident in the product you choose.

 **TRANQUILITY.**



VS

Leading Brands



We measure our products using **Capacity Under Pressure**, a proprietary test method that measures the product's retention capacity (the amount of fluid absorbed by the product that cannot be forced out with pressure), providing a more accurate representation of the product's absorbency under daily use. International Organization for Standardization testing (ISO 11948-1:1996) measures capacity as the amount of fluid a product can absorb without any pressure applied.

Real Life Stories

Tranquility® products are designed to make a real difference in the lives of wearers, caregivers and healthcare professionals. Read a few real stories below describing how people like you have benefited from Tranquility over the last 40 years.

“I have tried every product out there. None work as good as Tranquility. ... I no longer feel I need to take a change of clothes with me. I can go on long trips with peace of mind, and that alone is worth it!”

Karen, WEARER, KANSAS



“My wife is an early onset Alzheimer's sufferer of 11 years and has no control of bladder or bowels. **I had tried every adult diaper on the market, and they all failed to work, causing major cleanup, until I found Tranquility Premium OverNight™ Underwear.** ... Before Tranquility I was using three to four diapers per day and ... now I only use two with no additional cleaning.”

Chuck, CAREGIVER, MONTANA



“I started ordering Tranquility Premium Overnight Underwear for all my patients with significant incontinence problems. **EVERY patient or caregiver would call me and let me know how thrilled they were with the Tranquility products!**”

Staff member, PA-C HOME-BASED PRIMARY CARE, VA MEDICAL CENTER



To read more stories, visit TranquilityProducts.com/Real-Stories

Our simplified product lineup and new packaging

are designed to make it easier for you to select the right product, in the correct size, with the best features to manage your individual needs.



premium

Our highest-performing range of incontinence products. Find the ideal balance of protection, comfort and design to fit your unique needs, with a diverse selection of sizes, types and absorbency levels.



essential NEW!

High-quality products with **reliable performance** to support comfort and confidence for those with moderate to heavy incontinence.

Previously Known As



specialty

Products specifically innovated for unique use cases and **specialized needs**, empowering wearers and caregivers to achieve a new level of freedom.




connect
Smart Diaper Technology
COMING SOON!

Incontinence Product Guide


We believe you deserve an incontinence solution that addresses your unique journey, so we offer a wide selection of products, styles and sizes to meet your individual needs.

Incontinence Products




 **Pads, Liners and Guards** are absorbent pads designed for light to moderate incontinence. An adhesive strip holds the pad securely in your regular underwear to provide discreet protection.




 **Booster Pads** are flow-through pads designed to add absorbent capacity to your primary incontinence product. The flexible design allows for personalized positioning to add protection when and where you need it – ideal for pairing with protective underwear or brief-style products.




 **Protective Underwear** are pull-on products that fit like regular underwear. The soft, elastic waistband provides a discreet and comfortable fit that moves with you.




 **Underpads** provide an additional layer of protection to safeguard surfaces, such as bedding, furniture and wheelchairs, from heavy leakage and odor.



 **Briefs** are tab-style products that provide you with the highest level of protection and coverage. Available in either a poly (plastic) backsheet or a clothlike backsheet for added comfort and breathability. Briefs are designed to contain both urinary and bowel incontinence.



 **Disposable Swimwear** are worn under a bathing suit to help contain bowel incontinence. Designed for minimal swelling and will not break apart when wet. Not suited for containing urinary incontinence.

Incontinence Types

Stress:

Leakage of small amounts of urine due to sudden pressure on the bladder, such as coughing, sneezing, laughing or exercising.

Urge:

A sudden urge to urinate that is so strong it cannot be held long enough to reach the bathroom.

Mixed:

A combination of the symptoms of stress and urge incontinence.

Overflow:

Involuntary urination. Bladder muscles are unable to contain urine under pressure. Feeling of inability to empty bladder completely.

Reflex:

Involuntary urination without sensation of a full bladder. Unaware of the need to urinate.

Functional:

Inability to get to the bathroom fast enough due to mobility or cognitive impairments.

Features and Benefits

Compassion-driven innovation starts with you. Our research team is committed to learning what matters most to you and pairing that with the latest technology available to develop solutions for your unique needs.



LEAKAGE CONTROL

Absorbency target zones combined with other features, such as Kufguards®, gentle leg elastics and moisture-proof materials, are designed to secure leaks and maintain your confidence.



ODOR REDUCTION

Superabsorbent core wicks away fluid, neutralizing pH and reducing odor for your discretion.



SKIN HEALTH PROTECTION

Breathable sides combined with our superabsorbent core support an ideal microclimate that neutralizes pH and keeps your skin dry.



ADDITIONAL CAPACITY

Additional capacity extends the wear time of your primary product.




SOFT COMFORT

Soft, clothlike layer is gentle against skin.



SECURE PLACEMENT

Adhesive strip secures product in place to prevent sliding.



KUFGUARDS®

Inner leg cuffs channel fluids into the core, helping to contain high-volume episodes.



CLOSURE TABS

Advanced tape system securely fastens the brief while allowing for multiple openings.



BREATHABLE SIDES

Allows air to circulate around the body, promoting skin integrity.



OVERNIGHT

Maximum absorbent capacity for extended or nighttime protection.




WETNESS INDICATOR

Yellow lines turn blue when product reaches capacity to indicate change is needed.



LATEX-FREE

All Tranquility products are latex-free for your skin's safety and comfort.



Briefs

TOP
SELLER



ABSORBENCY:
34.0 oz/1005 mL (Size Large)



SmartCore® Briefs

MAXIMUM

Breathable sides help eliminate heat and moisture, providing the ideal micro-climate for skin integrity.

Breathable Sides		Kufguards*	Clothlike Backsheet	Micro-Hook Closure Tabs	Color-Coded by Size	Wetness Indicator	Latex-Free
Item	Size	Product Color	Waist/Hip	Weight	Capacity	Case Count	HCPC
2311	Small	Green	24" - 32"	85 - 100 lbs	22.0 oz/650 mL	100 (10 bags of 10)	T4521
2312	Medium	White	32" - 44"	100 - 150 lbs	28.0 oz/828 mL	96 (8 bags of 12)	T4522
2313	Large	Blue	45" - 58"	150 - 200 lbs	34.0 oz/1005 mL	96 (8 bags of 12)	T4523
2314	X-Large	Beige	56" - 64"	200 - 250 lbs	34.0 oz/1005 mL	72 (6 bags of 12)	T4524
2315	2XL/Bariatric	White	60" - 80"	250+ lbs	34.0 oz/1005 mL	32 (4 bags of 8)	T4543

TOP
SELLER



ABSORBENCY:
33.0 oz/975 mL (Size Large)



ATN™ Briefs (All-Through-the-Night)

MAXIMUM

A full-fit, poly-backed brief with expanded surface area providing longer wear time by absorbing multiple voids.

Kufguards®		Refastenable Tapes	Wetness Indicator	Latex-Free		
Item	Size	Waist/Hip	Weight	Capacity	Case Count	HCPC
2183	XS/Youth 6-7*	18" - 26"	42 - 90 lbs	18.5 oz/547 mL	100 (10 bags of 10)	T4533
2184	Small	24" - 32"	85 - 100 lbs	22.2 oz/657 mL	100 (10 bags of 10)	T4521
2185	Medium	32" - 44"	100 - 150 lbs	27.5 oz/812 mL	96 (8 bags of 12)	T4522
2186	Large	45" - 58"	150 - 200 lbs	33.0 oz/975 mL	96 (8 bags of 12)	T4523
2187	X-Large	56" - 64"	200 - 250 lbs	34.0 oz/1005 mL	72 (6 bags of 12)	T4524

*Clothlike backsheet, micro-hook closure tabs



ABSORBENCY:
21.4 oz/632 mL (Size Large)



SlimLine® Original Briefs

SUPER*

Classic full-fit, poly-backed brief with a slimmer profile makes it versatile for everyday use.

Kufguards*		Refastenable Tapes	Wetness Indicator	Latex-Free		
Item	Size	Waist/Hip	Weight	Capacity	Case Count	HCPC
2112	Youth 5-6**	16" - 22"	28 - 42 lbs	10.2 oz/302 mL	120 (10 bags of 12)	T4530
2166	XS/Youth 6-7**	18" - 26"	42 - 90 lbs	10.9 oz/322 mL	100 (10 bags of 10)	T4533
2120	Small	24" - 32"	85 - 100 lbs	12.7 oz/374 mL	100 (10 bags of 10)	T4521
2122	Medium	32" - 44"	100 - 150 lbs	20.4 oz/602 mL	96 (8 bags of 12)	T4522
2132	Large	45" - 58"	150 - 200 lbs	21.4 oz/632 mL	96 (8 bags of 12)	T4523
2134	X-Large	56" - 64"	200 - 250 lbs	23.7 oz/702 mL	72 (6 bags of 12)	T4524

*This absorbency category will change from Heavy to Super when Tranquility Premium packaging is updated in 2022.

**Clothlike backsheet, micro-hook closure tabs

Essential Breathable Briefs

HEAVY

Clothlike, breathable, tab-style brief offers reliable protection for heavy incontinence.

Breathable Sides		Kufguards®	Clothlike Backsheet	Adjustable Tabs	Color-Coded by Size	Wetness Indicator	Latex-Free
Item	Size	Product Color	Waist/Hip	Weight	Capacity	Case Count	HCPC
2743	XS/Youth 6-7	White	18" - 26"	42 - 90 lbs	12.0 oz/355 mL	100 (10 bags of 10)	T4533
2744	Small	Green	24" - 32"	85 - 100 lbs	12.0 oz/355 mL	100 (10 bags of 10)	T4521
2745	Medium	White	32" - 44"	100 - 150 lbs	15.8 oz/467 mL	96 (8 bags of 12)	T4522
2746	Large	Blue	44" - 58"	150 - 200 lbs	19.3 oz/571 mL	72 (6 bags of 12)	T4523
2747	X-Large	Beige	56" - 64"	200 - 250 lbs	23.0 oz/680 mL	64 (8 bags of 8)	T4524
2748	2XL/Bariatric	White	60" - 80"	250+ lbs	23.0 oz/680 mL	32 (4 bags of 8)	T4543



ABSORBENCY:
19.3 oz/571 mL (Size Large)




Essential Breathable Briefs

MODERATE


Clothlike, breathable, tab-style brief offers reliable protection for moderate incontinence.

Breathable Sides		Kufguards®	Clothlike Backsheet	Adjustable Tabs	Color-Coded by Size	Wetness Indicator	Latex-Free
Item	Size	Product Color	Waist/Hip	Weight	Capacity	Case Count	HCPC
2965-100	Medium	White	32" - 44"	100 - 150 lbs	13.5 oz/400 mL	100 (4 bags of 25)	T4522
2966-100	Large	Blue	44" - 58"	150 - 200 lbs	13.5 oz/400 mL	100 (4 bags of 25)	T4523
2967-100	X-Large	Beige	56" - 64"	200 - 250 lbs	13.5 oz/400 mL	100 (4 bags of 25)	T4524



ABSORBENCY:
13.5 oz/400 mL (Size Large)





Bariatric Briefs

AIR-Plus™ Bariatric Briefs **MAXIMUM**

100% breathable allowing for heat and moisture vapor to escape, enhancing skin integrity. Fits up to a 108" waist.

Breathable	Kufguards®	Clothlike Backsheet	Micro-Hook Closure Tabs	Wetness Indicator	Latex-Free	
Item	Size	Waist/Hip	Weight	Capacity	Case Count	HCPC
2195	4-5XL	70" - 108"	250+ lbs	34.0 oz/1005 mL	32 (4 bags of 8)	T4543



ABSORBENCY:
34.0 oz/1005 mL



HI-Rise™ Bariatric Briefs **MAXIMUM**

Full-fit brief with higher coverage in the front and back for a secure fit. Fits up to a 96" waist. Breathable, stretchy side panels.

Breathable Sides	Kufguards®	Clothlike Backsheet	Micro-Hook Closure Tabs	Wetness Indicator	Latex-Free	
Item	Size	Waist/Hip	Weight	Capacity	Case Count	HCPC
2192	3XL	64" - 96"	250+ lbs	34.0 oz/1005 mL	32 (4 bags of 8)	T4543



ABSORBENCY:
34.0 oz/1005 mL



Bariatric Briefs **MAXIMUM**

Full-fit brief with breathable, stretchy side panels comfortably fits up to a 96" waist.







Breathable Sides	Kufguards®	Clothlike Backsheet	Micro-Hook Closure Tabs	Wetness Indicator	Latex-Free	
Item	Size	Waist/Hip	Weight	Capacity	Case Count	HCPC
2190	3XL	64" - 96"	250+ lbs	34.0 oz/1005 mL	32 (4 bags of 8)	T4543



ABSORBENCY:
34.0 oz/1005 mL



Briefs Comparison	 SmartCore® Briefs Sizes: S - 2XL	 ATN™ Briefs (All-Through-the-Night) Sizes: XS - XL	 SlimLine® Original Briefs Sizes: Youth - XL	 Essential Breathable Briefs - Heavy Sizes: XS - 2XL	 Essential Breathable Briefs - Moderate Sizes: M - XL
Capacity Under Pressure*	 34.0 oz/1005 mL	 33.0 oz/975 mL	 21.4 oz/632 mL	 19.3 oz/571 mL	 13.5 oz/400 mL
Overnight/Extended Wear	●	●			
Clothlike Backsheet	●	● XS Only	● Youth & XS Only	●	●
Poly Backsheet		● S - XL Only	● S - XL Only		
Breathable Sides	●			●	●
Color-Coded by Size	●			●	●

Bariatric Briefs Comparison	 AIR-Plus™ Bariatric Briefs Size: 4-5XL	 HI-Rise™ Bariatric Briefs Size: 3XL	 Bariatric Briefs Size: 3XL
Capacity Under Pressure*	 34.0 oz/1005 mL	 34.0 oz/1005 mL	 34.0 oz/1005 mL
Overnight/Extended Wear	●	●	●
Clothlike Backsheet	●	●	●
Breathable Sides	●	●	●
100% Breathable	●		
Stretchable Side Panels	●	●	●
High-Rise Waist	●	●	

All Tranquility Briefs feature:

Kufguards®
 Refastenable Tabs
 Wetness Indicator



All Tranquility Products offer Real Life Protection:



LEAKAGE CONTROL



ODOR CONTROL




SKIN HEALTH PROTECTION



LATEX-FREE

*Capacity Under Pressure is a proprietary Tranquility test method that measures the product's retention capacity (the amount of fluid absorbed by the product that cannot be forced out with pressure), providing a more accurate representation of the product's absorbency under daily use. Capacity Under Pressure varies by product size. For comparison purposes, all capacities listed in this chart are for large-sized products.



Protective Underwear

#1
Preferred
Product



ABSORBENCY:
34.0 oz/1005 mL (Size Large)



Premium OverNight™ Underwear

MAXIMUM

Maximum absorbency provides extended wear time for a full night of uninterrupted sleep.

Breathable Sides Kufguards* Clothlike Backsheet Latex-Free

Item	Size	Waist/Hip	Weight	Capacity	Case Count	HCPC
2113	XS/Youth Large	17" - 28"	65 - 85 lbs	20.3 oz/600 mL	88 (4 bags of 22)	T4534
2114	Small/Youth X-Large	22" - 36"	80 - 125 lbs	20.3 oz/600 mL	80 (4 bags of 20)	T4525/T4534
2115	Medium*	34" - 48"	120 - 175 lbs	34.0 oz/1005 mL	72 (4 bags of 18)	T4526
2116	Large*	44" - 54"	170 - 210 lbs	34.0 oz/1005 mL	64 (4 bags of 16)	T4527
2117	X-Large*	48" - 66"	210 - 250 lbs	34.0 oz/1005 mL	56 (4 bags of 14)	T4528
2118	2XL/Bariatric*	62" - 80"	250+ lbs	34.0 oz/1005 mL	48 (4 bags of 12)	T4544

*Blue leg elastics denote OverNight from DayTime Underwear.



ABSORBENCY:
25.0 oz/740 mL (Size Large)



Premium DayTime™ Underwear

SUPER*

Discreet daytime protection for active individuals with heavier urine voids.

Breathable Sides Kufguards* Clothlike Backsheet Latex-Free

Item	Size	Waist/Hip	Weight	Capacity	Case Count	HCPC
2105	Medium	34" - 48"	120 - 175 lbs	25.0 oz/740 mL	72 (4 bags of 18)	T4526
2106	Large	44" - 54"	170 - 210 lbs	25.0 oz/740 mL	64 (4 bags of 16)	T4527
2107	X-Large	48" - 66"	210 - 250 lbs	25.0 oz/740 mL	56 (4 bags of 14)	T4528
2108	2XL/ Bariatric	62" - 80"	250+ lbs	25.0 oz/740 mL	48 (4 bags of 12)	T4544

*This absorbency category will change from Heavy to Super when Tranquility Premium packaging is updated in 2022. For sizes X-Small and Small, try OverNight Underwear.

Essential Underwear

HEAVY

Soft, pull-on product offers discreet and reliable protection for heavy incontinence.

Breathable Sides Kufguards® Clothlike Backsheet Stretchable Waistband Latex-Free

Item	Size	Waist/Hip	Weight	Capacity	Case Count	HCPC
2602	Youth Medium	15" - 25"	38 - 65 lbs	14.5 oz/429 mL	96 (8 bags of 12)	T4534
2603	XS/Youth Large	17" - 28"	65 - 85 lbs	14.5 oz/429 mL	96 (4 bags of 24)	T4534
2604	Small/Youth X-Large	22" - 36"	80 - 125 lbs	14.5 oz/429 mL	88 (4 bags of 22)	T4525/T4534
2605	Medium	34" - 48"	120 - 175 lbs	19.0 oz/562 mL	80 (4 bags of 20)	T4526
2606	Large	44" - 54"	170 - 210 lbs	19.0 oz/562 mL	72 (4 bags of 18)	T4527
2607	X-Large	48" - 66"	210 - 250 lbs	19.0 oz/562 mL	56 (4 bags of 14)	T4528
2608	2XL/Bariatric	62" - 80"	250+ lbs	19.0 oz/562 mL	48 (4 bags of 12)	T4544



TOP
SELLER

ABSORBENCY:
19.0 oz/562 mL (Size Large)



Essential Underwear

MODERATE

Soft, pull-on product offers discreet and reliable protection for moderate incontinence.


Breathable Sides Kufguards® Clothlike Backsheet Stretchable Waistband Latex-Free

Item	Size	Waist/Hip	Weight	Capacity	Case Count	HCPC
2974-100	Small/Youth XL	22" - 36"	80 - 125 lbs	10.1 oz/299 mL	100 (4 bags of 25)	T4525/T4534
2975-100	Medium	34" - 48"	120 - 175 lbs	13.5 oz/400 mL	100 (4 bags of 25)	T4526
2976-100	Large	44" - 54"	170 - 210 lbs	13.5 oz/400 mL	100 (4 bags of 25)	T4527
2977-100	X-Large	48" - 66"	210 - 250 lbs	13.5 oz/400 mL	100 (4 bags of 25)	T4528



ABSORBENCY:
13.5 oz/400 mL (Size Large)



Protective Underwear Comparison	 Premium OverNight™ Underwear Sizes: XS - 2XL	 Premium DayTime™ Underwear Sizes: M - 2XL	 Essential Underwear - Heavy Sizes: Youth - 2XL	 Essential Underwear - Moderate Sizes: S - XL
Capacity Under Pressure*	 34.0 oz/1005 mL	 25.0 oz/740 mL	 19.0 oz/562 mL	 13.5 oz/400 mL
Overnight	●			
Extended Wear	●	●		
Clothlike Backsheet	●	●	●	●
Breathable Sides	●	●	●	●
Kufguards*	●	●	●	●
Stretchable Waistband	●	●	●	●



All Tranquility Products offer Real Life Protection:



LEAKAGE CONTROL



ODOR CONTROL




SKIN HEALTH PROTECTION



LATEX-FREE

*Capacity Under Pressure is a proprietary Tranquility test method that measures the product's retention capacity (the amount of fluid absorbed by the product that cannot be forced out with pressure), providing a more accurate representation of the product's absorbency under daily use. Capacity Under Pressure varies by product size. For comparison purposes, all capacities listed in this chart are for large-sized products.



Pads & Liners



ABSORBENCY:
9.5 oz/280 mL



Male Guard

MAXIMUM

Absorbent core locks in fluid and controls odor. Contoured shape with side elastics forms a cup that conforms to a man's body.

Clothlike
Backsheet

Adhesive
Strip

Latex-Free

Item	Size	Length/Width	Capacity	Case Count	HCPC
2385	One Size	12.25" x 5.25"	9.5 oz/280 mL	104 (2 bags of 52)	T4535



ABSORBENCY:
15.2 oz/450 mL (Size Ultimate)



Personal Care Pads

MAXIMUM

Features two embossed channels to direct fluid into the pad and elastics to form a cup shape. Worn in regular underwear.

Clothlike
Backsheet

Adhesive
Strip

Latex-Free

Item	Size	Length/Width	Capacity	Case Count	HCPC
2380	Super*	10.5" x 5.5"	12.6 oz/375 mL**	96 (4 bags of 24)	T4535
2381	Ultimate*	13.5" x 6.5"	15.2 oz/450 mL***	96 (4 bags of 24)	T4535
2382	OverNight	16.5" x 7.25"	16.9 oz/500 mL***	96 (4 bags of 24)	T4535

*Individually wrapped for hygienic protection.

**This absorbency category will change from Moderate Protection to Super when Tranquility Premium packaging is updated in 2022.

*** This absorbency category will change from Moderate Protection to Maximum when Tranquility Premium packaging is updated in 2022.



ABSORBENCY:
15.5 oz/458 mL



Adult Liner

MAXIMUM

Thin, moisture-proof, backed liner is worn in regular underwear for a discreet appearance.

Adhesive
Strip

Latex-Free

Item	Size	Length/Width	Capacity	Case Count	HCPC
2078	One Size	24" x 9"	15.5 oz/458 mL	120 (4 bags of 30)	T4535

Essential Personal Care Pads

HEAVY

Clothlike, contoured pad designed for discreet placement in regular underwear.

Clothlike Backsheet Adhesive Strip Contoured Shape Latex-Free

Item	Size	Length/Width	Capacity	Case Count	HCPC
2881	Regular*	13.5" x 6.5"	11.8 oz/349 mL	96 (4 bags of 24)	T4535
2882	Long	16.5" x 7.25"	12.6 oz/375 mL	96 (4 bags of 24)	T4535

*Individually wrapped for hygienic protection.



ABSORBENCY:
11.8 oz/349 mL (Size Regular)



Essential Personal Care Pads

MODERATE

Clothlike, contoured pad designed for discreet placement in regular underwear.

Clothlike Backsheet Adhesive Strip Contoured Shape Latex-Free

Item	Size	Length/Width	Capacity	Case Count	HCPC
2880	Regular*	10.5" x 5.5"	10.1 oz/299 mL	96 (4 bags of 24)	T4535

*Individually wrapped for hygienic protection.





ABSORBENCY:
10.1 oz/299 mL (Size Regular)



Did You Know?

Pads and liners have a moisture-proof backing and are meant to be worn in regular underwear. They should not be used as a liner inside another incontinence product, as this may actually cause moisture to run out of the product and leak. **Do not use with other incontinence products.**





Boosters

What is a Booster Pad?

Booster pads have a flow-through design allowing them to fill to capacity before passing additional moisture into a primary incontinence product. Booster pads should always be placed inside another disposable incontinence product that has a moisture-proof backing. **Do not use with regular underwear.**

For added capacity



BOOSTERS



PULL-ON
UNDERWEAR



TAB-STYLE
BRIEFS



2XL-5XL
BARIATRIC
BRIEFS

TopLiner® Booster Pads **MAXIMUM**

Thin design allows for placement in target areas where additional absorption is needed.

Flow-Through Design | Clothlike | Adhesive Strip | Latex-Free

Item	Size	Length/Width	Capacity	Case Count	HCPC	Best Suited for Host Garment Sizes
2072	Mini*	10.5" x 2.75"	5.7 oz/169 mL	200 (8 bags of 25)	T4535	Youth, XS
2070	Regular**	14" x 4"	10.9 oz/322 mL	200 (8 bags of 25)	T4535	S, M, L
2060	Super**	15" x 4.25"	14.2 oz/420 mL	200 (8 bags of 25)	T4535	XL, 2XL

*This absorbency category will change from Added Protection to Super when Tranquility Premium packaging is updated in 2022.

**This absorbency category will change from Added Protection to Maximum when Tranquility Premium packaging is updated in 2022.



ABSORBENCY:

14.2 oz/420 mL (Size Regular)



TopLiner® Booster Contour Pads **MAXIMUM**

Wide hourglass shape helps contain bowel. Soft texture and no adhesive strip makes cleanup easier.

Flow-Through Design | Clothlike | Latex-Free

Item	Size	Length/Width	Capacity	Case Count	HCPC	Best Suited for Host Garment Sizes
3096	Contour	21.5" x 13.5"	13.6 oz/403 mL	120 (10 bags of 12)	T4535	M, L, XL
3097	Super-Plus	32" x 14"	27.5 oz/813 mL	96 (8 bags of 12)	T4535	2XL, 3XL, 4-5XL



ABSORBENCY:

13.6 oz/403 mL (Size Contour)





ABSORBENCY:
6.7 oz/198 mL (Size Regular)



Essential Booster Pads

HEAVY

Flow-through pad designed to add capacity to any disposable incontinence product.

Flow-Through Design Clothlike Adhesive Strip Latex-Free

Item	Size	Length/Width	Capacity	Case Count	HCPC	Best Suited for Host Garment Sizes
2770	Youth	12" x 3.25"	5.5 oz/163 mL	180 (6 bags of 30)	T4535	Youth, XS
2760	Regular	12" x 4.25"	6.7 oz/198 mL	200 (8 bags of 25)	T4535	S, M, L
2762	Long	15" x 4.25"	8.0 oz/236 mL	200 (8 bags of 25)	T4535	XL, 2XL



ABSORBENCY:
5.3 oz/157 mL (Size Regular)



Essential Booster Pads

MODERATE

Flow-through pad designed to add capacity to any disposable incontinence product.

Flow-Through Design Clothlike Adhesive Strip Latex-Free

Item	Size	Length/Width	Capacity	Case Count	HCPC	Best Suited for Host Garment Sizes
19244	Regular	12" x 4"	5.3 oz/157 mL	200 (8 bags of 25)	T4535	S, M, L, XL

Did You Know?

Booster pad placement is important. **The pads can be placed in targeted areas where additional absorption is needed.**





Booster pads can be folded and even doubled for added protection, but should never be cut.



Host Garment – The pad's adhesive side should face out, away from the body. Make sure the pad is fully contained within the host incontinence product; no part should extend beyond the edges.

For Women – Place the pad more toward the bottom.

For Men – Place the pad more toward the front.



Underpads

TOP
SELLER



ABSORBENCY:
34.0 oz/1005 mL



AIR-Plus™ Extra-Strength Breathable Underpads

MAXIMUM

Strong backsheet for boosting, turning and repositioning.*
100% breathable; allows heat and moisture vapor to escape.

Breathable Clothlike Backsheet Latex-Free

Item	Length/Width	Capacity	Case Count	HCPC
2711	30" x 36"	34.0 oz/1005 mL	40 (4 bags of 10)	T4541/A4554

*Product should be used in accordance with Safe Patient Handling Guidelines. Lay printed side down against bed, chair or other surface.



ABSORBENCY:
34.0 oz/1005 mL (Size 30" x 36")



AIR-Plus™ Breathable Underpads

MAXIMUM

100% breathable backsheet allows heat and moisture vapor to escape. Ideal for use with low-air-loss bed systems.

Breathable Clothlike Backsheet Latex-Free

Item	Length/Width	Capacity	Case Count	HCPC
2709	23" x 36"	29.0 oz/850 mL	40 (4 bags of 10)	T4541/A4554
2710	30" x 36"	34.0 oz/1005 mL	40 (4 bags of 10)	T4541/A4554

Lay printed side down against bed, chair or other surface.



ABSORBENCY:
34.0 oz/1005 mL



Heavy-Duty Underpads

MAXIMUM

Strong, clothlike backsheet resists tearing as an individual moves on the underpad.

Clothlike Backsheet Latex-Free

Item	Length/Width	Capacity	Case Count	HCPC
2088	30" x 36"	34.0 oz/1005 mL	60 (6 bags of 10)	T4541/A4554



ABSORBENCY:
34.0 oz/1005 mL



Peach Sheet Underpads

MAXIMUM

Poly-backed underpad has four adhesive tapes to secure it to wheelchairs, bedding and furniture.

Adhesive Tabs Latex-Free

Item	Length/Width	Capacity	Case Count	HCPC
2074	21.5" x 32.5"	34.0 oz/1005 mL	96 (8 bags of 12)	T4541/A4554

Essential Underpads MODERATE

Soft and smooth underpad is comfortable against skin and protects surfaces from leakage and odor.

Latex-Free



ABSORBENCY:
18.5 oz/550 mL (Size 36" x 36")



Item	Length/Width	Product Color	Capacity	Case Count	HCPC
2676	22" x 36"	Blue	2.0 oz/59 mL	150 (5 bags of 30)	T4541/A4554
2677	28" x 30"	Green	10.6 oz/314 mL	100 (10 bags of 10)	T4541/A4554
2679	36" x 36"	Beige	18.5 oz/550 mL	50 (5 bags of 10)	T4541/A4554

Did You Know?

Underpads are typically used to protect mattresses and furniture from urine damage and reduce linen laundry. They also can protect other surfaces, including wheelchairs, vehicle seats and the floor beneath a portable or bedside toilet.

In fact, they have many uses – to temporarily absorb light leakage of kitchen pipes and refrigerator drips, soak up car oil or other spills in the garage, or as a protective cover for baby diaper changes when using a public changing station.



“As a nurse, I know the challenge of avoiding skin breakdown in elderly patients. I am very pleased with the Tranquility® products for my own Dad. **They really help to keep urine away from his skin,** especially during the longer periods between diaper changes.”

Lenneal, NURSE



“

“I am a retired registered nurse and a full-time caregiver for my husband of 53 years. ... I have tried all the incontinence products available at the drugstore, supermarket and healthcare stores ... and **I am more than pleased with the absorption, odor control and versatility of your products ...**

These Tranquility® products have been a Godsend to me.”

Rosemary,

RETIRED REGISTERED NURSE
AND FULL-TIME CAREGIVER



EXTRA
CARE



Extra Care

Cleansing Wipes

Paraben- and alcohol-free, hypoallergenic with a mild scent. Generous 9" x 13" size allows for quicker cleanup.



Latex-Free

Item	Size	Length/Width	Case Count	HCPC
3101	Full Case	9" x 13"	672 (12 packs of 56)	A4335
3121	Mini Case	9" x 13"	224 (4 packs of 56)	A4335

ThinLiner Moisture Management® Sheets

100% breathable, soft, clothlike sheets reduce skin-on-skin friction. Collects, absorbs and retains moisture.



Breathable

Clothlike

Latex-Free

Item	Length/Width	Absorbency*	Case Count	HCPC
3190	6" x 10"	6 oz/177 mL	200 (20 bags of 10)	A9999
3191	6" x 14"	8.5 oz/251 mL	200 (20 bags of 10)	A9999


*Moisture absorbency; not capacity under pressure.

Did You Know?


ThinLiner Moisture Management Sheets offer multiple applications and uses. They are ideal for tucking, folding and wrapping in or around skin folds, the groin area, under breasts – wherever light-to-moderate moisture is present.




ThinLiner Sheets can be folded to fit various sizes and applications, but should never be cut.




Relieves tightened digits



Collects wound exudate




Wicks away moisture



Absorbs loose stool*

*ThinLiner Sheets should always be placed inside another disposable incontinence product that has a moisture-proof backing. **Do not use with regular underwear.**



Disposable Swimwear

Swimmates®

These discreet, unisex, disposable underwear are worn under a bathing suit. Features Kufguards® (inner leg cuffs) and moisture barrier to help contain bowel incontinence. Designed for minimal swelling and does not break apart when wet.




Kufguards® Latex-Free

Item	Size	Waist/Hip	Weight	Case Count	HCPC
2844	Small/Youth XL	22" - 36"	80 - 125 lbs	88 (4 bags of 22)	T4525/T4534/A9999
2845	Medium	34" - 48"	120 - 175 lbs	80 (4 bags of 20)	T4526/A9999
2846	Large	44" - 54"	170 - 210 lbs	72 (4 bags of 18)	T4527/A9999
2847	X-Large	48" - 66"	210 - 250 lbs	56 (4 bags of 14)	T4528/A9999
2848	XXL-Plus/Bariatric	62" - 80"	250+ lbs	48 (4 bags of 12)	T4544/A9999

“

“I am a physical education teacher at an elementary school with a therapeutic pool ... **Swimmates® is the reliable swim diaper that I count on to provide protection to my older students.** Having the proper attire keeps the pool safe for all students and allows the students the opportunity to excel in the aquatic environment.”

Paul, PHYSICAL EDUCATION TEACHER, MINNESOTA



“

“We recommend Tranquility® products because it works! We tell customers 'you get what you pay for' and it really applies to incontinence challenges.”

P.W.,

HOME-CARING DEALER



Reference Materials



How to Measure

Use a tape measure to measure in inches the **widest part of your body between waist and hip**. Use this number to determine product size.



Briefs Size Guide


Size	Waist/Hip	Weight	Women's Pant Size	Men's Pant Size
Youth	16" - 22"	28 - 42 lbs	N/A	N/A
X-Small	18" - 26"	42 - 90 lbs	0 - 3	18 - 26
Small	24" - 32"	85 - 100 lbs	4 - 6	24 - 32
Medium	32" - 44"	100 - 150 lbs	8 - 12	32 - 44
Large	44" - 58"	150 - 200 lbs	14 - 18	44 - 58
X-Large	56" - 64"	200 - 250 lbs	18 - 24	56 - 64
2XL	60" - 80"	250+ lbs	26 - 28	60+
3XL	64" - 96"	250+ lbs	28+	64+
4-5XL	70" - 108"	250+ lbs	28+	70+

Good Fit



Brief should fit snugly with the placement of the tape tabs **on the hips** and no gaps around the legs.

Too Big




Too Small



Protective Underwear Size Guide


Size	Waist/Hip	Weight	Women's Pant Size	Men's Pant Size
Youth	15" - 25"	38 - 65 lbs	N/A	N/A
X-Small	17" - 28"	65 - 85 lbs	0 - 3	17 - 28
Small	22" - 36"	80 - 125 lbs	4 - 6	22 - 36
Medium	34" - 48"	120 - 175 lbs	6 - 12	34 - 48
Large	44" - 54"	170 - 210 lbs	12 - 18	44 - 54
X-Large	48" - 66"	210 - 250 lbs	18 - 24	48 - 66
2XL	62" - 80"	250+ lbs	24+	62+

Good Fit



Underwear should fit snugly **touching the crotch** with no gaps around the legs.

Too Big



Too Small



Sizing is not always perfect; the sizes shown are recommendations, not a guarantee. Please call **Tranquility Customer Care at 1-866-865-6101** if you have questions.



Tranquility® Youth

Our youth products bridge the gap when your child has outgrown toddler diapers and youth pull-ons.

Toddler Diapers

Commonly known as Huggies® or Pampers®

Size 5/6

Waist/Hip 16" - 22"
Weight 28 - 42 lbs



**SlimLine®
Original Briefs**
2112
10.2 oz/302 mL

Size 6/7

Waist/Hip 18" - 26"
Weight 42 - 90 lbs



ATN™ Briefs
2183
18.5 oz/547 mL



**SlimLine®
Original Briefs**
2112
10.2 oz/302 mL



**Essential Breathable
Briefs - Heavy**
2743
12.0 oz/355 mL

Youth Pull-ons

Commonly known as GoodNites®

Youth Size M

Waist/Hip 15" - 25"
Weight 38 - 65 lbs



**Essential Underwear
- Heavy**
2602
14.5 oz/429 mL

Youth Size L

Waist/Hip 17" - 28"
Weight 65 - 85 lbs



**Premium OverNight™
Underwear**
2113
20.3 oz/600 mL



**Essential Underwear
- Heavy**
2603
14.5 oz/429 mL

Youth Size XL

Waist/Hip 22" - 36"
Weight 80 - 125 lbs



**Premium OverNight™
Underwear**
2114
20.3 oz/600 mL



**Essential Underwear
- Heavy**
2604
14.5 oz/429 mL



**Essential Underwear
- Moderate**
2974-100
10.1 oz/299 mL

Other Youth Products



**Essential Booster
Pads* - Heavy**
2770
Youth
12" x 3.25"



Swimmates®
2884
Small/Youth XL
Waist/Hip 22" - 36"
Weight 80 - 125 lbs



For complete product details,
see individual items listed in this catalog.

Newly Designed with Options for Medicaid



Tranquility® Essential Breathable Briefs in Moderate and Heavy protection may be available – **at little or no cost** – for those with special healthcare needs, starting at age 3, under state Medicaid benefits.*



*Incontinence products are not covered by Medicaid in all states. We may be able to help you find a Medicaid provider in your state that offers covered Tranquility products.

INTRODUCING Tranquility® Essential

Coming
Summer
2021

The Select® and ComfortCare™ brands will now be known as Tranquility Essential. Use this chart to find the equivalent new Tranquility Essential product:

Previous Items

essential Tranquility Essential Equivalent



Select Briefs and Soft n' Breathable Briefs

3665	Youth
3666	X-Small
2620	Small
2624	Medium
2634	Large
2635	X-Large
2627	Medium
2628	Large
2629	X-Large



Tranquility Essential Breathable Briefs – Heavy

2743	XS/Youth 6-7
2743	XS/Youth 6-7
2744	Small
2745	Medium
2746	Large
2747	X-Large
2745	Medium
2746	Large
2747	X-Large
2748	2XL/Bariatric NEW SIZE



Select Absorbent Disposable Underwear

2602	Youth
2603	X-Small
2604	Small
2605	Medium
2606	Large
2607	X-Large
2608	2XL/Bariatric



Tranquility Essential Underwear – Heavy

2602	Youth Medium
2603	XS/Youth Large
2604	Small/Youth X-Large
2605	Medium
2606	Large
2607	X-Large
2608	2XL/Bariatric



Select Personal Care Pads

2881	Regular
2882	Extra



Tranquility Essential Personal Care Pads – Heavy

2881	Regular
2882	Long



Select Booster Pads

2760	Medium
2762	Extended
2770	Small



Tranquility Essential Booster Pads – Heavy

2760	Regular
2762	Long
2770	Youth



Previous Items

essential Tranquility Essential Equivalent



Select Personal Care Pads

2880 Light



Tranquility Essential Personal Care Pads – Moderate

2880 Regular



Select Underpads

2675 Large
2677 Extra Large
2679 Ultra Large



Tranquility Essential Underpads – Moderate

2676 22" x 36"
2677 28" x 30"
2679 36" x 36"



ComfortCare Briefs

2965-100 Medium
2966-100 Large
2967-100 X-Large



Tranquility Essential Breathable Briefs – Moderate

2965-100 Medium
2966-100 Large
2967-100 X-Large



ComfortCare Disposable Absorbent Underwear

2974-100 Small/Youth XL
2975-100 Medium
2976-100 Large
2977-100 X-Large



Tranquility Essential Underwear – Moderate

2974-100 Small/Youth XL
2975-100 Medium
2976-100 Large
2977-100 X-Large



ComfortCare Booster Pads

19244 Booster



Tranquility Essential Booster Pads – Moderate

19244 Regular



Have questions or need assistance?

Call 1-800-467-3224 or visit TranquilityProducts.com/TQEssential

Quick Order Reference Guide

premium

Item	Product Description	Size	Unit/Case	Pkg/Case	HCPC	Page #
2060	Tranquility® TopLiner® Booster Pads	Super	200	8/25	T4535	23
2070	Tranquility® TopLiner® Booster Pads	Regular	200	8/25	T4535	23
2072	Tranquility® TopLiner® Booster Pads	Mini	200	8/25	T4535	23
2074	Tranquility® Peach Sheet Underpads	21.5" x 32.5"	96	8/12	T4541/A4554	26
2078	Tranquility® Adult Liner	One Size	120	4/30	T4535	20
2088	Tranquility® Heavy-Duty Underpads	30" x 36"	60	6/10	T4541/A4554	26
2105	Tranquility® Premium DayTime™ Underwear	Medium	72	4/18	T4526	16
2106	Tranquility® Premium DayTime™ Underwear	Large	64	4/16	T4527	16
2107	Tranquility® Premium DayTime™ Underwear	X-Large	56	4/14	T4528	16
2108	Tranquility® Premium DayTime™ Underwear	2XL/Bariatric	48	4/12	T4544	16
2112	Tranquility® SlimLine® Original Briefs	Youth 5-6	120	10/12	T4530	10
2113	Tranquility® Premium OverNight™ Underwear	X-Small/Youth Large	88	4/22	T4534	16
2114	Tranquility® Premium OverNight™ Underwear	Small/Youth X-Large	80	4/20	T4525/T4534	16
2115	Tranquility® Premium OverNight™ Underwear	Medium	72	4/18	T4526	16
2116	Tranquility® Premium OverNight™ Underwear	Large	64	4/16	T4527	16
2117	Tranquility® Premium OverNight™ Underwear	X-Large	56	4/14	T4528	16
2118	Tranquility® Premium OverNight™ Underwear	2XL/Bariatric	48	4/12	T4544	16
2120	Tranquility® SlimLine® Original Briefs	Small	100	10/10	T4521	10
2122	Tranquility® SlimLine® Original Briefs	Medium	96	8/12	T4522	10
2132	Tranquility® SlimLine® Original Briefs	Large	96	8/12	T4523	10
2134	Tranquility® SlimLine® Original Briefs	X-Large	72	6/12	T4524	10
2166	Tranquility® SlimLine® Original Briefs	X-Small/Youth 6-7	100	10/10	T4533	10
2183	Tranquility® ATN™ Briefs (All-Through-the-Night)	X-Small/Youth 6-7	100	10/10	T4533	10
2184	Tranquility® ATN™ Briefs (All-Through-the-Night)	Small	100	10/10	T4521	10
2185	Tranquility® ATN™ Briefs (All-Through-the-Night)	Medium	96	8/12	T4522	10
2186	Tranquility® ATN™ Briefs (All-Through-the-Night)	Large	96	8/12	T4523	10
2187	Tranquility® ATN™ Briefs (All-Through-the-Night)	X-Large	72	6/12	T4524	10
2190	Tranquility® Bariatric Briefs	3XL	32	4/8	T4543	13
2192	Tranquility® HI-Rise™ Bariatric Briefs	3XL	32	4/8	T4543	13
2195	Tranquility® AIR-Plus™ Bariatric Briefs	4-5XL	32	4/8	T4543	13
2311	Tranquility® SmartCore® Briefs	Small	100	10/10	T4521	10
2312	Tranquility® SmartCore® Briefs	Medium	96	8/12	T4522	10
2313	Tranquility® SmartCore® Briefs	Large	96	8/12	T4523	10
2314	Tranquility® SmartCore® Briefs	X-Large	72	6/12	T4524	10
2315	Tranquility® SmartCore® Briefs	2XL/Bariatric	32	4/8	T4543	10
2380	Tranquility® Personal Care Pads	Super	96	4/24	T4535	20
2381	Tranquility® Personal Care Pads	Ultimate	96	4/24	T4535	20
2382	Tranquility® Personal Care Pads	OverNight	96	4/24	T4535	20
2385	Tranquility® Male Guard	One Size	104	2/52	T4535	20
2709	Tranquility® AIR-Plus™ Breathable Underpads	23" x 36"	40	4/10	T4541/A4554	26
2710	Tranquility® AIR-Plus™ Breathable Underpads	30" x 36"	40	4/10	T4541/A4554	26
2711	Tranquility® AIR-Plus™ Extra-Strength Breathable Underpads	30" x 36"	40	4/10	T4541/A4554	26
3096	Tranquility® TopLiner® Booster Contour Pads	Contour	120	10/12	T4535	23
3097	Tranquility® TopLiner® Booster Contour Pads	Super-Plus	96	8/12	T4535	23
3101	Tranquility® Cleansing Wipes	Full Case	672	12/56	A4335	29
3121	Tranquility® Cleansing Wipes	Mini Case	224	4/56	A4335	29
3190	Tranquility® ThinLiner Moisture Management® Sheets	6" x 10"	200	20/10	A9999	29
3191	Tranquility® ThinLiner Moisture Management® Sheets	6" x 14"	200	20/10	A9999	29

Quick Order Reference Guide

essential

Item	Product Description	Size	Unit/Case	Pkg/Case	HCPC	Page #
2602	Tranquility® Essential Underwear - Heavy	Youth Medium	96	8/12	T4534	17
2603	Tranquility® Essential Underwear - Heavy	X-Small/Youth Large	96	4/24	T4534	17
2604	Tranquility® Essential Underwear - Heavy	Small/Youth X-Large	88	4/22	T4525/T4534	17
2605	Tranquility® Essential Underwear - Heavy	Medium	80	4/20	T4526	17
2606	Tranquility® Essential Underwear - Heavy	Large	72	4/18	T4527	17
2607	Tranquility® Essential Underwear - Heavy	X-Large	56	4/14	T4528	17
2608	Tranquility® Essential Underwear - Heavy	2XL/Bariatric	48	4/12	T4544	17
2676	Tranquility® Essential Underpads - Moderate	22" x 36"	150	5/30	T4541/A4554	27
2677	Tranquility® Essential Underpads - Moderate	28" x 30"	100	10/10	T4541/A4554	27
2679	Tranquility® Essential Underpads - Moderate	36" x 36"	50	5/10	T4541/A4554	27
2743	Tranquility® Essential Breathable Briefs - Heavy	X-Small/Youth 6-7	100	10/10	T4533	11
2744	Tranquility® Essential Breathable Briefs - Heavy	Small	100	10/10	T4521	11
2745	Tranquility® Essential Breathable Briefs - Heavy	Medium	96	8/12	T4522	11
2746	Tranquility® Essential Breathable Briefs - Heavy	Large	72	6/12	T4523	11
2747	Tranquility® Essential Breathable Briefs - Heavy	X-Large	64	8/8	T4524	11
2748	Tranquility® Essential Breathable Briefs - Heavy	2XL/Bariatric	32	4/8	T4543	11
2760	Tranquility® Essential Booster Pads - Heavy	Regular	200	8/25	T4535	24
2762	Tranquility® Essential Booster Pads - Heavy	Long	200	8/25	T4535	24
2770	Tranquility® Essential Booster Pads - Heavy	Youth	180	6/30	T4535	24
2880	Tranquility® Essential Personal Care Pads - Moderate	Regular	96	4/24	T4535	21
2881	Tranquility® Essential Personal Care Pads - Heavy	Regular	96	4/24	T4535	21
2882	Tranquility® Essential Personal Care Pads - Heavy	Long	96	4/24	T4535	21
2965-100	Tranquility® Essential Breathable Briefs - Moderate	Medium	100	4/25	T4522	11
2966-100	Tranquility® Essential Breathable Briefs - Moderate	Large	100	4/25	T4523	11
2967-100	Tranquility® Essential Breathable Briefs - Moderate	X-Large	100	4/25	T4524	11
2974-100	Tranquility® Essential Underwear - Moderate	Small/Youth X-Large	100	4/25	T4525/T4534	17
2975-100	Tranquility® Essential Underwear - Moderate	Medium	100	4/25	T4526	17
2976-100	Tranquility® Essential Underwear - Moderate	Large	100	4/25	T4527	17
2977-100	Tranquility® Essential Underwear - Moderate	X-Large	100	4/25	T4528	17
19244	Tranquility® Essential Booster Pads - Moderate	Regular	200	8/25	T4535	24

specialty

Item	Product Description	Size	Unit/Case	Pkg/Case	HCPC	Page #
2844	Swimmates®	Small/Youth X-Large	88	4/22	T4525/T4534/A9999	31
2845	Swimmates®	Medium	80	4/20	T4526/A9999	31
2846	Swimmates®	Large	72	4/18	T4527/A9999	31
2847	Swimmates®	X-Large	56	4/14	T4528/A9999	31
2848	Swimmates®	XXL-Plus/Bariatric	48	4/12	T4544/A9999	31

connect
Smart Diaper Technology

COMING SOON!



GUARANTEE

You deserve real-life protection that works, whether you're sleeping, sitting or on-the-move. Tranquility products are clinically proven to absorb and contain their stated cup capacity, no matter how much pressure is applied. We stand behind our products. If you are not satisfied for any reason or need additional support, let us know. Our care team is here to help you.



CUSTOMER CARE

Consumers may contact the Tranquility Care Center for FREE product samples, private consultation and product information.

1-866-865-6101

CustomerService@pbenet.com

TranquilityProducts.com/contact



DEALERS & DISTRIBUTORS

Call your Tranquility Representative or contact us for information on Tranquility products, marketing services and support, and consumer sampling programs.

1-800-467-3224 ext. 6

CustomerService@pbenet.com

TranquilityProducts.com/dealer



P.O. Box 129
Dunbridge, OH 43414
419-352-1551
PrincipleBusinessEnterprises.com

Follow us:  

Family Owned and Operated for 3 Generations

We are inspired by experiences with loved ones who deserve better incontinence solutions. Our team combines engineering expertise and clinical best practices to deliver products that provide dignity, freedom and confidence for wearers and caregivers.

– Andy, CEO

