

# BD Saf-T-Intima™ IV Catheters

## BD Angiocath™ IV Catheters



### BD Saf-T-Intima™

Closed IV Catheter System with BD Vialon™ Biomaterial

Catalog #	Gauge	Length	Flow Rate	O.D.	I.D.
383312	24	0.75" (0.7 mm x 19 mm) PRN	22 mL/min	.0260-.0290	.0190-.0220
383313	24	0.75" (0.7 mm x 19 mm) Y	22 mL/min	.0260-.0290	.0190-.0220
383322	22	0.75" (0.9 mm x 19 mm) PRN	36 mL/min	.0330-.0360	.0240-.0270
383323	22	0.75" (0.9 mm x 19 mm) Y	36 mL/min	.0330-.0360	.0240-.0270
383335	20	1.00" (1.1 mm x 25 mm) PRN	55 mL/min	.0410-.0440	.0300-.0330
383336	20	1.00" (1.1 mm x 25 mm) Y	55 mL/min	.0410-.0440	.0300-.0330
383346	18	1.00" (1.3 mm x 25 mm) Y	82 mL/min	.0500-.0530	.0370-.0400

### BD Angiocath™

IV Catheter with FEP Polymer

Catalog #	Gauge	Length	Flow Rate	O.D.	I.D.
381112	24	0.75" (0.7 mm x 19 mm)	17 mL/min	.0275-.0285	.0175-.0205
381123	22	1.00" (0.9 mm x 25 mm)	28 mL/min	.0325-.0335	.0215-.0255
381137	20	1.88" (1.1 mm x 48 mm)	42 mL/min	.0425-.0435	.0275-.0315
381147	18	1.88" (1.3 mm x 48 mm)	79 mL/min	.0514-.0526	.0345-.0385
381157	16	1.88" (1.7 mm x 48 mm)	147 mL/min	.0633-.0677	.0443-.0487
382258	16	3.25" (1.7 mm x 83 mm)	127 mL/min	.0633-.0677	.0443-.0487
382259	16	5.25" (1.7 mm x 133 mm)	108 mL/min	.0633-.0677	.0443-.0487
381167	14	1.88" (2.1 mm x 48 mm)	277 mL/min	.0822-.0838	.0602-.0648
382268	14	3.25" (2.1 mm x 83 mm)	249 mL/min	.0822-.0838	.0602-.0648
382269	14	5.25" (2.1 mm x 133 mm)	219 mL/min	.0822-.0838	.0602-.0648
382277	12	3.00" (2.8 mm x 76 mm)	449 mL/min	.1060-.1140	.0880-.0920
382287	10	3.00" (3.4 mm x 76 mm)	609 mL/min	.1300-.1380	.1040-.1080



Helping all people live healthy lives

Product Attributes & Features	BD Angiocath™ Autoguard™	BD Insyte™ Autoguard™	BD Saf-T-Intima™	BD Angiocath™
Product Description	Shielded IV Catheter	Shielded IV Catheter	Closed IV Catheter System	IV Catheter
Catheter Material	FEP Polymer	BD Vialon™ Biomaterial	BD Vialon™ Biomaterial	FEP Polymer
BD Instaflex™ Needle Technology	20, 22, 24 G	20, 22, 24 G	NA	NA
Needle Protection	Needle Encapsulation	Needle Encapsulation	Needle Encapsulation	NA
Safety Activation	Clinician Controlled	Clinician Controlled	Passive	NA

## BD Medical Systems

BD, the world's leading manufacturer of IV access devices, knows the importance to both clinician and patient of thoughtfully designed, effective, and safe medical devices. The BD family of IV catheters provides you with products designed to make IV access simple, safe, and reliable.

BD Customer Service: 1.888.237.2762

Visit us at: [www.bd.com/infusion](http://www.bd.com/infusion)



**BD Medical**  
9450 South State Street  
Sandy, UT 84070  
[www.bd.com](http://www.bd.com)

# BD Autoguard™ Shielded IV Catheters



## BD Insyte™ Autoguard™

Shielded IV Catheter with BD Vialon™ Biomaterial

Catalog #	Gauge	Length	Flow Rate	O.D.	I.D.
381412	24	0.75" (0.7 mm x 19 mm)	20 mL/min	.0260-.0290	.0190-.0220
381423	22	1.00" (0.9 mm x 25 mm)	35 mL/min	.0330-.0360	.0240-.0270
381433	20	1.00" (1.1 mm x 25 mm)	65 mL/min	.0410-.0440	.0300-.0330
381434	20	1.16" (1.1 mm x 30 mm)	60 mL/min	.0410-.0440	.0300-.0330
381437	20	1.88" (1.1 mm x 48 mm)	55 mL/min	.0410-.0440	.0300-.0330
381444	18	1.16" (1.3 mm x 30 mm)	105 mL/min	.0500-.0530	.0370-.0400
381447	18	1.88" (1.3 mm x 48 mm)	95 mL/min	.0500-.0530	.0370-.0400
381454	16	1.16" (1.7 mm x 30 mm)	220 mL/min	.0670-.0700	.0520-.0550
381457	16	1.77" (1.7 mm x 45 mm)	205 mL/min	.0670-.0700	.0520-.0550
381467	14	1.75" (2.1 mm x 45 mm)	330 mL/min	.0830-.0860	.0665-.0705



## BD Insyte-N™ Autoguard™

Shielded IV Catheter with BD Vialon™ Biomaterial

Catalog #	Gauge	Length	Flow Rate	O.D.	I.D.
381411	24	0.56" (0.7 mm x 14 mm)	20 mL/min	.0278-.0288	.0207-.0214
381511	24	0.56" (0.7 mm x 14 mm) Wings	20 mL/min	.0278-.0288	.0207-.0214



## BD Angiocath-N™ Autoguard™

Shielded IV Catheter with FEP Polymer

Catalog #	Gauge	Length	Flow Rate	O.D.	I.D.
381720	24	0.56" (0.7 mm x 14 mm)	20 mL/min	.0278-.0288	.0207-.0214



## BD Insyte™ Autoguard™ Winged

Shielded IV Catheter with BD Vialon™ Biomaterial

Catalog #	Gauge	Length	Flow Rate	O.D.	I.D.
381512	24	0.75" (0.7 mm x 19 mm)	20 mL/min	.0260-.0290	.0190-.0220
381523	22	1.00" (0.9 mm x 25 mm)	35 mL/min	.0330-.0360	.0240-.0270
381533	20	1.00" (1.1 mm x 25 mm)	65 mL/min	.0410-.0440	.0300-.0330
381534	20	1.16" (1.1 mm x 30 mm)	60 mL/min	.0410-.0440	.0300-.0330
381537	20	1.88" (1.1 mm x 48 mm)	55 mL/min	.0410-.0440	.0300-.0330
381544	18	1.16" (1.3 mm x 30 mm)	105 mL/min	.0500-.0530	.0370-.0400
381547	18	1.88" (1.3 mm x 48 mm)	95 mL/min	.0500-.0530	.0370-.0400
381554	16	1.16" (1.7 mm x 30 mm)	220 mL/min	.0670-.0700	.0520-.0550
381557	16	1.77" (1.7 mm x 45 mm)	205 mL/min	.0670-.0700	.0520-.0550



## BD Angiocath™ Autoguard™

Shielded IV Catheter with FEP Polymer

Catalog #	Gauge	Length	Flow Rate	O.D.	I.D.
381700	24	0.75" (0.7 mm x 19 mm)	20 mL/min	.0260-.0290	.0190-.0220
381701	22	1.00" (0.9 mm x 25 mm)	35 mL/min	.0330-.0360	.0240-.0270
381702	20	1.00" (1.1 mm x 25 mm)	65 mL/min	.0410-.0440	.0300-.0330
381703	20	1.16" (1.1 mm x 30 mm)	60 mL/min	.0410-.0440	.0300-.0330
381704	20	1.88" (1.1 mm x 48 mm)	55 mL/min	.0410-.0440	.0300-.0330
381705	18	1.16" (1.3 mm x 30 mm)	105 mL/min	.0500-.0530	.0370-.0400
381706	18	1.88" (1.3 mm x 48 mm)	95 mL/min	.0500-.0530	.0370-.0400
381707	16	1.16" (1.7 mm x 30 mm)	220 mL/min	.0670-.0700	.0520-.0550
381708	16	1.77" (1.7 mm x 45 mm)	205 mL/min	.0670-.0700	.0520-.0550
381709	14	1.75" (2.1 mm x 45 mm)	330 mL/min	.0830-.0860	.0665-.0705

# BD IV Catheter Family



## BD Insyte™ Autoguard™

Shielded IV Catheter with BD Vialon™ Biomaterial



## BD Insyte™ Autoguard™ Winged

Shielded IV Catheter with BD Vialon™ Biomaterial



## BD Insyte-N™ Autoguard™

Shielded IV Catheter with BD Vialon™ Biomaterial



## BD Angiocath™ Autoguard™

Shielded IV Catheter with FEP Polymer



## BD Angiocath-N™ Autoguard™

Shielded IV Catheter with FEP Polymer



## BD Saf-T-Intima™

Closed IV Catheter System with BD Vialon™ Biomaterial



## BD Angiocath™

IV Catheter with FEP Polymer