

Cunningham Incontinence Clamp

Easy adjustable
to give proper fit
and comfort.

- For the effective control of enuresis (bed wetting) in males of all ages
- Constructed of soft sponge rubber and can be adjusted to fit comfortably
- Available in 3 sizes

JUVENILE | REGULAR | LARGE



BARD

Cunningham Incontinence Clamp



Easy adjustable to give proper fit and comfort.



A very light pressure on the urethral canal on the underside of the penis will prevent any leakage. The necessary pressure is achieved by the hump in the under part of the Cunningham Clamp. For proper fit and comfort the upper part of the clamp can be easily shaped by the fingers.



The clamp should be replaced every 3 months, or sooner if the foam deteriorates.

- For the effective control of enuresis (bed wetting) in males of all ages
 - Constructed of soft sponge rubber and can be adjusted to fit comfortably
 - Available in 3 sizes

JUVENILE | **REGULAR** | **LARGE**
for Boys for Adults

Description	Size	Order-No.
Cunningham Incontinence Clamp, juvenile	1.5 Inch / 38 mm	004052
Cunningham Incontinence Clamp, regular	2 Inch / 51 mm	004053
Cunningham Incontinence Clamp, large	3 Inch / 76 mm	004054

Ordering amount: 1 piece

Please consult product labels and inserts for any indications, contraindications, hazards, warnings, cautions and instructions for use.

BARD

C. R. Bard GmbH
Wachhausstraße 6
D-76227 Karlsruhe
TEL: +49 721 9445-124
FAX: +49 721 9445-100
www.bard.de

Bard Medica S.A.
Seestrasse 64
CH-8942 Oberrieden/Zürich
TEL: +41 44 722 53 60
FAX: +41 44 722 53 70
www.crbard.ch

Bard Medica S.A.
Thaliastraße 125a/1/5
A-1160 Wien
TEL: +43 1 49 49 130
FAX: +43 1 49 49 130-30
www.bard.at

**Bard Czech Republic
s.r.o.**
Táboršká 619/46
CZ-140 00 Praha 4
TEL: +420 242 408 630
FAX: +420 242 410 185
www.bard.com

Bard Poland Sp.z o.o.
ul. Cybernetyki 7B
PL-02-677 Warszawa
TEL: +48 22 321 09 30
FAX: +48 22 321 09 38
www.crbard.pl