

# Personal Protective Apparel

Head-to-Toe Protection for Healthcare



# Protection at every level

You want to meet your team's PPE needs consistently. That's why we empower you with every level of protection. From our comprehensive line of protective apparel to next-level supply chain reliability and PPE expertise, you can feel confident that your team has the high-quality protection they need.

When you partner with us, we support you with a higher level of service, making it easier to keep up with industry standards and follow best practices. That way, you're covered no matter what lies ahead.

**Recommended Gowns 4-5**

**AAMI Isolation Gowns 6-7**

**Cover Gowns 8-9**

**Protective Gowns 10-11**

**Footwear 12-13**

**Headwear and Beard Cover 14-15**

**Coveralls 16-17**

**Lab Coats 18-20**

**Scrubs 21**

**Patient Wear 22**

# Understanding gown materials

## Spunbond polypropylene

For basic coverage, this non-woven fabric bonds fibers together to form a single layer.

## SMS material

Spunbond/meltblown/spunbond (SMS) is a multi-layer fabric composed of inner layers of meltblown polypropylene between outer layers of spunbond polypropylene.

## Coated polypropylene (poly-coated)

Soft, spunbond polypropylene is coated with a layer of polyethylene (plastic) film.

## Polyethylene film

Polyethylene (plastic) film is composed of a single layer of film and depending on thickness can be appropriate for low to moderate fluid exposure.

## Microporous

Microporous fabric is composed of a soft nonwoven inner layer, laminated with breathable microporous film.



## Choose the right gown

Our recommendations for your workplace environment.

### A Housewide, Med/Surg and ICU

Our QuickComply™ AAMI Level 2 gowns feature an over-the-head style which improves compliance while reducing donning and removing time. Side ties eliminate the difficulties of tying behind the back, and thumb loop wrists allow for a good gown and glove fit.

Item No.	Description
NONLV240	Regular/Large
NONLV240XL	X-Large

*Additional items on page 7*

### B Endoscopy and Emergency

AAMI Level 3 gowns are made from comfortable premium multi-ply material. Features a classic tie waist, hook and loop neck and knit cuffs.

Item No.	Description
NONLV325	Regular/Large
NONLV325XL	X-Large

*Additional items on page 7*

### C Decontamination

Made from microporous material, these AAMI Level 3 gowns are suitable for use in decontamination. Features tape tab neck, knit cuff wrists and strong-sealed seams designed to withstand tough conditions.

Item No.	Description
NONLV350	Regular/Large
NONLV350XL	X-Large

*Additional items on page 11*

### D Chemotherapy

USP 800 Compliant Gown with protective polyethylene coated materials and seams tested against top 11 most hazardous chemotherapy drugs per ASTM F739 and ASTM D6978. Features hook and loop neck, waist ties, knit cuffs and taped seams.

Item No.	Description
NON25457M	Medium
NON25457	Regular/Large
NON25457XL	X-Large
NON25457XXL	XX-Large

*Additional items on page 11*



**A** Housewide, Med/Surg and ICU



**B** Endoscopy and Emergency



**C** Decontamination



**D** Chemotherapy



# Why AAMI?

## AAMI Guidelines: A Reliable Benchmark for Protection

Under AAMI Guidelines, protective apparel must meet minimum performance standards for strength, barrier protection and fluid resistance to ensure the safety of healthcare workers and patients.



## Levels of Barrier Protection – AAMI PB70:2022

AAMI Level	Test Method	Requirements at 4% AQL	Anticipated Risk of Fluid Exposure
Level 1	AATCC 42 - Water Resistance: Impact Penetration	≤ 4 g	Low
Level 2	AATCC 42 - Water Resistance: Impact Penetration	≤ 1.0 g	Low to Moderate
	AATCC127 - Water Resistance: Hydrostatic Pressure	≥ 20 cm	
Level 3	AATCC 42 - Water Resistance: Impact Penetration	≤ 1.0 g	Moderate to High
	AATCC127 - Water Resistance: Hydrostatic Pressure	≥ 50 cm	

### AAMI Level 1 Isolation Gown

Recommended for: Med/Surg, Visitors, Contact Precautions

- Made from multi-ply material
- Side ties eliminate the difficulties of tying behind the back
- Choose between over-the-head or neck tie styles

Item No.	Description	Pkg.
<b>Medium Weight, Neck Ties, Elastic Wrists</b>		
NONLV100	Regular/Large, Yellow	100/cs
NONLV100XL	X-Large, Yellow	100/cs
<b>Medium Weight, Over-the-Head, Thumb Loop Wrists</b>		
NONLV140	Regular/Large, Yellow	100/cs
NONLV140XL	X-Large, Yellow	100/cs

### AAMI Level 2 Isolation Gown

Recommended for: Med/Surg, ICU, Emergency, Airborne and Contact Precautions

- Made from multi-ply material with side ties to eliminate the difficulties of tying behind the back
- Choose between over-the-head or tape tab neck styles
- Available with elastic wrists or thumb loops

Item No.	Description	Pkg.
<b>Medium Weight, Tape Tab Neck, Elastic Wrists</b>		
NONLV200M	Medium, Yellow	100/cs
NONLV200	Regular/Large, Yellow	100/cs
NONLV200XL	X-Large, Yellow	100/cs
NONLV200XXL	XX-Large, Yellow	100/cs
NONLV200XXXL	XXX-Large, Yellow	100/cs
<b>Medium Weight, Over-the-head, Overlapping Back Panels, Elastic Wrists, Thumb Loops</b>		
NONLV235	Regular/Large, Yellow	100/cs
NONLV235XL	X-Large, Yellow	100/cs
<b>Medium Weight, Over-the-Head, Thumb Loops</b>		
NONLV240	Regular/Large, Yellow	100/cs
NONLV240XL	X-Large, Yellow	100/cs
NONLV240XXL	XX-Large, Yellow	100/cs
NONLV240XXXL	XXX-Large, Yellow	100/cs

### AAMI Level 2 Polyethylene-Coated Isolation Gown

Recommended for: Endoscopy, Trauma and ER, Airborne and Contact Precautions

- Available over-the-head
- Available with thumb loops or knit cuffs

Item No.	Description	Pkg.
<b>Polyethylene Coated, Over-the-Head, Knit Cuff</b>		
NON27340	Regular/Large, Yellow	100/cs
NON27340XL	X-Large, Yellow	100/cs
<b>Polyethylene Coated, Over-the-Head, Thumb Loops</b>		
NON27360	Regular/Large, Yellow	100/cs
NON27360XL	X-Large, Yellow	100/cs

### AAMI Level 3 Isolation Gown

Recommended for: ICU, Decontam/SPD, Emergency/Trauma, Airborne, Contact and Droplet Precautions

- Comfortable materials with classic tie waist
- Choose between tape tab, hook and loop and over the head styles
- Available with elastic wrists, knit cuffs, or thumb loops

Item No.	Description	Pkg.
<b>Heavyweight Multi-Ply, Tape Tab Neck, Elastic Wrists</b>		
NONLV315	Regular/Large, Yellow	100/cs
NONLV315XL	X-Large, Yellow	100/cs
<b>Premium Multi-Ply, Hook and Loop Neck, Knit Cuffs</b>		
NONLV325	Regular/Large, Blue	50/cs
NONLV325XL	X-Large, Blue	50/cs
<b>Medium Weight Multi-Ply, Over-The-Head, Thumb Loop Wrists</b>		
NONLV340	Regular/Large, Yellow	50/cs
<b>Breathable Microporous, Tape Tab Neck, Knit Cuffs</b>		
NONLV350	Regular/Large, White	50/cs
NONLV350XL	X-Large, White	50/cs

## Gown Features



#### Side ties

Available on AAMI Level 1, 2 and 3 gowns



#### Neck tie style

Available on AAMI Level 1 gowns



#### Thumb loop wrists

Available on AAMI Level 1, 2 and 3 gowns



#### Tape tab neck style

Available on AAMI Level 2 and 3 gowns



CRI4000



Neck ties help ensure the gown remains securely in place.

### Classic Cover Gown

Recommended for: Visitors, Food Services

- Available in lightweight polypropylene or multi-layer materials
- Generously sized for full coverage and flexibility
- Side ties eliminate difficulties of tying behind the back
- Neck ties and elastic wrists

Item No.	Description	Pkg.
<b>Lightweight Polypropylene</b>		
CRI4000	Regular/Large, Yellow	50/cs
CRI4001	X-Large, Yellow	50/cs
CRI4003	XXX-Large, Yellow	50/cs
<b>Medium Weight Multi-Layer</b>		
NON27SMS2	Regular/Large, Yellow	100/cs
NON27SMS2XL	X-Large, Yellow	100/cs
NON27SMS2XXXL	XXX-Large, Yellow	100/cs
<b>Lightweight Multi-Layer</b>		
NON27236B	Regular/Large, Blue	50/cs
NON27236LTXL	X-Large, Yellow, Extra Long Center Tie	50/cs
NON27236LT2XL	XX-Large, Yellow, Extra Long Center Tie	50/cs





NON27SMS5



Convenient over-the-head, open back style reduces donning and removing time.

### Over-the-Head, Open Back Cover Gown

Recommended for: Visitors, Food Services

- Medium weight multi-layer
- Generously sized for coverage and flexibility
- Over-the-head style reduces donning and removing time
- Thumb loop wrists

Item No.	Description	Pkg.
<b>Medium Weight Multi-Layer</b>		
NON27SMS5	Regular/Large	100/cs
NON27SMS5XL	X-Large	100/cs

### Polyethylene Gown

Recommended for: Emergency, ICU, Food Services, Contact and Droplet Precautions

- Over-the-head style reduces donning and doffing time
- Thumb loop wrists help to prevent sleeves from sliding
- Available in bulk “B” or folded in a dispenser box

Item No.	Description	Pkg.
<b>Polyethylene, Blue</b>		
CRI5000	Regular/Large	75/cs
CRI5000B	Regular/Large	100/cs
CRI5010	Universal	75/cs
CRI5001	X-Large	75/cs
CRI5001B	X-Large	100/cs
CRI5003	XXX-Large	75/cs
<b>Premium Polyethylene, Blue</b>		
NONTH150	Regular/Large	75/cs
NONTH180	Universal	75/cs
NONTH200	X-Large	75/cs
NONTH100B	Regular/Large	100/cs
NONTH200B	X-Large	100/cs

### Over-the-Head Microporous Gown

Recommended for: Endoscopy, GI Lab, ICU, Emergency, Maternity

- Made of microporous breathable material
- Over-the-head style reduces donning and doffing time
- Thumb loop wrists prevent sleeve slide
- Individually folded in a dispenser box

Item No.	Description	Pkg.
NONTH350	Regular/Large	40/cs
NONTH400	X-Large	40/cs



CRI5010

NONTH350

Learn more about bouffant caps and shoe covers.

See pages 12–14

### Polyethylene-Coated Gown

Recommended for: ICU, Endoscopy, GI Lab

- Neck and side ties
- Available with elastic wrists or knit cuffs

Item No.	Description	Pkg.
<b>Polyethylene Coated, Yellow</b>		
CRI4010	Elastic Wrists, Regular/Large	50/cs
CRI4011	Elastic Wrists, X-Large	50/cs
CRI4011CF	Knit Cuffs, X-Large	50/cs
CRI4012CF	Knit Cuffs, XX-Large	50/cs
CRI4040PCC	Knit Cuffs, Yellow, Large	100/cs

### Chemotherapy Gown

Recommended for: Oncology, Radiology, Lab, Pharmacy

- CRI23457 is TAA Compliant and features a hook and loop neck, elastic cuffs and side waist ties
- NON25457 series features taped seams, hook and loop neck, knit cuffs and waist ties and offers AAMI Level 2 fluid resistance
- NON27457 series is made from breathable, bi-laminate material and features hook and loop neck, knit cuffs and waist ties

Item No.	Description	Pkg.
<b>Essential Protective Polyethylene Coated</b>		
CRI23457	Regular/Large	100/cs
CRI23457XL	X-Large	100/cs
<b>Protective Polyethylene Coated AAMI Level 2</b>		
NON25457M	Medium	30/cs
NON25457	Regular/Large	30/cs
NON25457XL	X-Large	30/cs
NON25457XXL	XX-Large	30/cs
<b>Protective Premium Breathable Material</b>		
NON27457M	Medium	30/cs
NON27457	Regular/Large	30/cs
NON27457XL	X-Large	30/cs
NON27457XXL	XX-Large	30/cs

### Over-the-Head Polyethylene-Coated Gown

Recommended for: Emergency, ICU, Maternity

- Over-the-head style reduces donning and doffing time
- Available with thumb loop wrists, elastic wrists and thumb loops, or knit cuffs

Item No.	Description	Pkg.
CRI4500	Thumb Loops, Regular/Large, Yellow	100/cs
<b>Polyethylene Coated</b>		
NON26278	Knit Cuffs, Regular/Large, Blue	50/cs
NON26278XL	Knit Cuffs, X-Large, Blue	50/cs
NON26278XXL	Knit Cuffs, XX-Large, Blue	50/cs
NON26276	Thumb Loops, Regular/Large, Blue	50/cs
NON26276XL	Thumb Loops, X-Large, Blue	50/cs
NON26276XXL	Thumb Loops, XX-Large, Blue	50/cs

USP800 Compliant Gowns tested against the top 11 most hazardous chemotherapy drugs per ASTM F739 and ASTM D6978:

- Dacarbazine
- Carmustine (BiCNU)
- Cyclophosphamide (Cytosan)
- Doxorubicin Hydrochloride
- 5-Fluorouracil
- Cisplatin
- Etoposide (Toposar)
- Paclitaxel (Taxol)
- Thio-Tepa
- Mitomycin C
- Vincristine Sulfate





NON29859

### Multi-Layer Shoe Covers

- Multi-layer material features stripe non-skid pattern for increased traction
- Range of sizes to fit all shoe types

Item No.	Description	Pkg.
NON29858	Heavyweight, Regular/Large	100/bx, 300/cs
NON29859	Heavyweight, X-Large	100/bx, 200/cs
NON28858	Medium Weight, Regular/Large	100/bx, 300/cs
NON28859	Medium Weight, X-Large	100/bx, 200/cs
NON28852	Medium Weight, XX-Large	100/bx, 200/cs



CRI2002



CRI2000

### Spunbond Polypropylene Shoe Covers

- Made of spunbond polypropylene
- Choice of non-skid pattern or smooth bottom

Item No.	Description	Pkg.
NON28758	Non-Skid, Regular/Large	100/bx, 300/cs
NON28759	Non-Skid, X-Large	100/bx, 200/cs
NON28752	Non-Skid, XX-Large	100/bx, 200/cs
CRI2002	Non-Skid, Regular/Large	100/bx, 300/cs
CRI2003	Non-Skid, X-Large	100/bx, 200/cs
CRI2000	Smooth Bottom, Regular/Large	100/bx, 300/cs
CRI2001	Smooth Bottom, X-Large	100/bx, 200/cs

Size	Shoe Size Conversion	
	Women's	Men's
<b>Regular/Large</b>	Fits up to size 13.5	Fits up to size 12
<b>X-Large</b>	Fits up to size 16.5	Fits up to size 15
<b>XX-Large</b>	Fits up to size 18.5	Fits up to size 17



CRI2010

### Polyethylene Shoe Covers

- Smooth bottom for situations where traction is not essential
- Made from polyethylene material

Item No.	Description	Pkg.
CRI2010	One Size Fits Most	1,000/cs

### Polyethylene-Coated Shoe Covers

- Medium-weight coated spunbond polypropylene
- Available with seam on bottom or seamless
- Non-skid treated for better traction
- Bulk packaged for Cleanroom use

Item No.	Description	Pkg.
<b>Non-Skid/Bottom Seam</b>		
NON30858	Regular/Large	300/cs
NON30859	X-Large	300/cs
NON30852	XX-Large	300/cs
<b>Non-Skid/Seamless</b>		
NON31858	Regular/Large	300/cs
NON31859	X-Large	300/cs
NON31852	XX-Large	300/cs



NON31852



## Knee-High Boot Covers

- Choose between non-skid bottom or foam strips for increased traction
- Available with hook and loop closure or elastic ankle

Item No.	Description	Pkg.
NON27348P	Premium Breathable, Bi-laminate, Foam Strip Bottom, Hook and Loop Closure, Regular/Large	50/bx, 150/cs
NON27348PXL	Premium Breathable, Bi-laminate, Foam Strip Bottom, Hook and Loop Closure, X-Large	50/bx, 150/cs
NON27349PXL	Thigh High Premium Breathable, Bi-laminate, Foam Strip Bottom, Hook and Loop Closure, X-Large	50/bx, 150/cs
NON27200	Laminated Polypropylene, Foam Strip Bottom, Hook and Loop Closure, Regular/Large	50/bx, 150/cs
NON27200XL	Laminated Polypropylene, Foam Strip Bottom, Hook and Loop Closure, X-Large	50/bx, 150/cs
NON27165XL	Medium-Weight Poly-Coated, Non-skid, Hook and Loop Closure, X-Large	200/cs
NON27150	Poly-coated, Foam Strip Bottom, Hook and Loop Closure, Regular/Large	50/bx, 150/cs
NON27150XL	Poly-coated, Foam Strip Bottom, Hook and Loop Closure, X-Large	50/bx, 150/cs
NON27144	Poly-coated, Non-skid, Hook and Loop Closure, Regular/Large	50/bx, 150/cs
NON27144XL	Poly-coated, Non-skid, Hook and Loop Closure, X-Large	50/bx, 150/cs
NON27143	Poly-coated to the Ankle, SMS Material to the Knee, Non-skid, Elastic, Regular/Large	50/bx, 150/cs
NON27143XL	Poly-coated to the Ankle, SMS Material to the Knee, Non-skid, Elastic, X-Large	50/bx, 150/cs



NON27200



NON27144



NON27348P




NON27143



NON28229




CRI1009




NON30233B




NON28230R




CRI1004



CRI1005



CRI1009



CRI1207

### Multi-Layer Bouffant Caps

- Variety of prints and sizes, all come 100/bx, 500/cs
- Choice of heavy or medium weight multi-ply or heavyweight single-ply

Item No.	Description	Pkg.
NON30232B	21", Blue	500/cs
NON30233B	24", Blue	500/cs
NON30223B	28", Blue	500/cs
NON28229	24", Brush Stroke	500/cs
NON28230R	24", Breast Cancer Ribbon	500/cs
NON28533B	24", Blue, Heavyweight	500/cs

### Polypropylene Bouffant Caps

- Variety of colors and sizes, all come 100/bx, 500/cs
- Choice of medium weight polypropylene or lightweight polypropylene

Item No.	Description	Pkg.
CRI1000	21", Green	500/cs
CRI1001	21", Blue	500/cs
CRI1002	21", White	500/cs
CRI1003	24", White	500/cs
CRI1004	24", Blue	500/cs
CRI1005	24", Green	500/cs
CRI1009	24", Red	500/cs
CRI1202	24", Blue, Medium Weight	500/cs
CRI1203	24", White, Heavyweight	500/cs
CRI1207	24", Purple, Medium Weight	500/cs
CRI1404AZ	24", Blue, Clip Cap	100/cs

**Surgical headwear that meets AORN Guidelines for Surgical Headwear**

### Surgical Headwear

- Better coverage from the back of the neck to sideburn area
- One size fits most—provides all-day comfort and a secure fit
- No elastic band in front means cap won't leave a mark
- Quality, longer-lasting design creates less waste than bouffant caps

Item No.	Description	Pkg.
NON61950	Premium Nonwoven Material	600/cs
NON61990	Scrim Material	500/cs
NON61980	Multi-layer Material	500/cs



NON61980



NON61950



NON61990

### Surgeon Cap

- Available in scrim, multi-layer or spunlace materials
- One size fits most
- Ties in the back
- Polypropylene crown for breathability

Item No.	Description	Pkg.
NON28625	Scrim Material	500/cs
NON28626	Multi-layer Material	500/cs
NON28628	Spunlace Material	500/cs



NON28625

### Surgeon Hood and Beard Cover

- Made of multi-layer or polypropylene material
- One size fits most

Item No.	Description	Pkg.
NONSH100C	Surgeons Hood with Multi-Layer Material, Polypropylene Crown	300/cs
NONSH200	Multi-layer Hood	300/cs
NONSH600	Polypropylene Hood	300/cs
NONSH400	Polypropylene Beard Cover	1,000/cs



NONSH100C



NONSH200



NONSH600



NONSH400

## Multi-Layer Coveralls

Recommended for: Staff, Visitors

- Made from multi-layer polypropylene

Item No.	Description	Pkg.
<b>Elastic Wrists, Straight Ankles</b>		
NONCV250XL	X-Large	25/cs
NONCV250XXL	XX-Large	25/cs
<b>Elastic Wrists and Ankles</b>		
NONCV280XL	X-Large	25/cs
NONCV280XXL	XX-Large	25/cs
NONCV280XXXL	XXX-Large	25/cs
NONCV2804XL	XXXX-Large	25/cs
NONCV2805XL	XXXXX-Large	25/cs

## Spunbond Coveralls

Recommended for: Staff, Visitors

- Available in heavy weight polypropylene material
- Choose between elastic wrist and straight ankle or straight wrist and straight ankle

Item No.	Description	Pkg.
<b>Heavyweight, Elastic Wrists, Straight Ankles</b>		
NONCV150L	Large	25/cs
NONCV150XL	X-Large	25/cs
NONCV150XXL	XX-Large	25/cs
NONCV150XXXL	XXX-Large	25/cs
NONCV150XXXXL	XXXX-Large	25/cs
<b>Elastic Wrists and Ankles</b>		
NONCV120M	Medium	25/cs
NONCV120L	Large	25/cs
NONCV120XL	X-Large	25/cs
NONCV120XXL	XX-Large	25/cs
NONCV120XXXL	XXX-Large	25/cs
NONCV1204XL	X-Large	25/cs
NONCV1205XL	X-Large	25/cs



NONCV250XL



NONCV150L



### Breathable Coveralls

Recommended for: Vivariums, Pharmaceutical Research, Manufacturing, Cleanrooms

- Microporous coveralls provide protection for extended periods of wear
- Full length zipper in front
- Elastic wrists with thumb loops and elastic ankles
- Meets ASTM F1670

Item No.	Description	Pkg.
NONCV940S	Small	25/cs
NONCV940M	Medium	25/cs
NONCV940L	Large	25/cs
NONCV940XL	X-Large	25/cs
NONCV940XXL	XX-Large	25/cs
NONCV940XXXL	XXX-Large	25/cs
NONCV940XXXXL	XXXX-Large	25/cs
NONCV9405XL	XXXXX-Large	25/cs

### Breathable Coveralls

- The front and sleeve are made from microporous breathable material, which meets ASTM F1670
- The back is made from spunbonded material that enhances ventilation
- Elastic wrist and ankles
- Thumb loop wrists prevent sleeve slide

Item No.	Description	Pkg.
NONCV960S	Small	25/cs
NONCV960M	Medium	25/cs
NONCV960L	Large	25/cs
NONCV960XL	X-Large	25/cs
NONCV960XXL	XX-Large	25/cs
NONCV960XXXL	XXX-Large	25/cs
NONCV9604XL	XXXX-Large	25/cs
NONCV9605XL	XXXXX-Large	25/cs

### Static-Dissipative Microporous Breathable Coveralls

Recommended for: Vivariums, Pharmaceutical Research, Manufacturing, Cleanrooms

- Made from static-dissipative microporous breathable material
- Full zipper front
- Elastic wrists and thumb loops
- Attached hood and boots
- Meets ASTM F1670

Item No.	Description	Pkg.
NONCV980M	Medium	25/cs
NONCV980L	Large	25/cs
NONCV980XL	X-Large	25/cs
NONCV980XXL	XX-Large	25/cs
NONCV980XXXL	XXX-Large	25/cs
NONCV9804XL	XXXX-Large	25/cs
NONCV9805XL	XXXXX-Large	25/cs



NONCV960M

NONCV940M

### ComfortSure™ Lab Jacket with Knit Collar and Cuffs

- Premium Comfortsure™ material provides soft, comfortable protection
- Two reinforced pockets
- Convenient loop for hanging

Item No.	Description	Pkg.
NONRP602S	Small	30/cs
NONRP602M	Medium	30/cs
NONRP602L	Large	30/cs
NONRP602XL	X-Large	30/cs
NONRP602XXL	XX-Large	30/cs



NONRP602M

### Breathable Fluid-Proof Lab Coat

Recommended for: Vivariums, Pharmaceutical Research, Manufacturing, Cleanrooms

- Breathable, static dissipative microporous material
- Knit collar, front zipper closure and knit cuffs
- Three reinforced pockets
- Meets ASTM F1670

Item No.	Description	Pkg.
NONSW1750S	Small	30/cs
NONSW1750M	Medium	30/cs
NONSW1750L	Large	30/cs
NONSW1750XL	X-Large	30/cs
NONSW1750XXL	XX-Large	30/cs
NONSW1750XXXL	XXX-Large	30/cs
NONSW17504XL	XXXX-Large	30/cs



NONSW1750L

### Breathable Fluid-Proof Snap Lab Coat

- Breathable, static-dissipative microporous material
- Snap front with knit cuffs and collar
- Suitable for use in nonsterile cleanrooms
- Meets ASTM F1670

Item No.	Description	Pkg.
NONSW1760S	Small	30/cs
NONSW1760M	Medium	30/cs
NONSW1760L	Large	30/cs
NONSW1760XL	X-Large	30/cs
NONSW1760XXL	XX-Large	30/cs
NONSW1760XXXL	XXX-Large	30/cs



NONSW500M

### Classic Lab Coat with Traditional Collar and Knit Cuffs

- Multi-layer material provides comfort and protection
- Three reinforced pockets
- Convenient loop for hanging

Item No.	Description	Pkg.
<b>White, Snap Front</b>		
NONSW100S	Small	30/cs
NONSW100M	Medium	30/cs
NONSW100L	Large	30/cs
NONSW100XL	X-Large	30/cs
NONSW100XXL	XX-Large	30/cs
NONSW100XXXL	XXX-Large	30/cs
<b>Blue, Snap Front</b>		
NONSW400S	Small	30/cs
NONSW400M	Medium	30/cs
NONSW400L	Large	30/cs
NONSW400XL	X-Large	30/cs
NONSW400XXL	XX-Large	30/cs

### Disposable Spunbond Lab Coat

- Made with spunbond polypropylene material for breathability and comfort
- Coats fall just below the knee and include a convenient loop for hanging
- Chest and waist pockets for storage

Item No.	Description	Pkg.
<b>White, Elastic Wrists, Traditional Collar, Snap Front</b>		
NONSW150M	Medium	30/cs
NONSW150L	Large	30/cs
NONSW150XL	X-Large	30/cs
NONSW150XXL	XX-Large	30/cs
NONSW150XXXL	XXX-Large	30/cs
<b>Blue, Elastic Wrists, Traditional Collar, Snap Front</b>		
NONSW450XL	X-Large	30/cs
NONSW450XXL	XX-Large	30/cs
NONSW450XXXL	XXX-Large	30/cs

### Classic Lab Coat with Knit Collar and Cuffs

- Multi-layer material provides comfort and protection
- Three reinforced pockets
- Convenient loop for hanging

Item No.	Description	Pkg.
<b>White</b>		
NONSW500S	Small	30/cs
NONSW500M	Medium	30/cs
NONSW500L	Large	30/cs
NONSW500XL	X-Large	30/cs
NONSW500XXL	XX-Large	30/cs
NONSW500XXXL	XXX-Large	30/cs
<b>Blue</b>		
NONSW600S	Small	30/cs
NONSW600M	Medium	30/cs
NONSW600L	Large	30/cs
NONSW600XL	X-Large	30/cs
NONSW600XXL	XX-Large	30/cs



### Classic Lab Jacket with Knit Collar and Cuffs

- Multi-layer material provides comfort and protection
- Three reinforced pockets
- Convenient loop for hanging

Item No.	Description	Pkg.
<b>White</b>		
NONRP500S	Small	30/cs
NONRP500M	Medium	30/cs
NONRP500L	Large	30/cs
NONRP500XL	X-Large	30/cs
NONRP500XXL	XX-Large	30/cs
<b>Blue</b>		
NONRP600S	Small	30/cs
NONRP600M	Medium	30/cs
NONRP600L	Large	30/cs
NONRP600XL	X-Large	30/cs
NONRP600XXL	XX-Large	30/cs

### Breathable Frock

Recommended for: Vivariums, Pharmaceutical Research, Manufacturing, Cleanrooms

- Breathable, static-dissipative microporous material
- Snap collar and snap front closure
- Elastic wrists
- Made from ASTM F1670 passing material

Item No.	Description	Pkg.
NONSW1765S	Small	30/cs
NONSW1765M	Medium	30/cs
NONSW1765L	Large	30/cs
NONSW1765XL	X-Large	30/cs
NONSW1765XXL	XX-Large	30/cs
NONSW1765XXXL	XXX-Large	30/cs
NONSW17654XL	XXXX-Large	30/cs

### Spunbond Polypropylene Frock

- Comfort and breathable multi-layer, heavyweight material
- Ideal for less critical areas, such as visitor areas and animal research
- Traditional collar, hook-and-loop closure and elastic wrists

Item No.	Description	Pkg.
NONSW025S	Small	30/cs
NONSW025M	Medium	30/cs
NONSW025L	Large	30/cs
NONSW025XL	X-Large	30/cs
NONSW025XXL	XX-Large	30/cs
NONSW025XXXL	XXX-Large	30/cs
NONSW0254XL	XXXX-Large	30/cs

### Protective Sleeve Covers

- Tested against ASTM F1670 to protect against blood penetration
- Available in premium laminate or microporous material
- Elastic bands at both ends secure the covers at the wrists and upper arms



NON18SLEEVE1

Item No.	Description	Pkg.
NON18SLEEVE1	20" (50.8 cm) Long White Microporous Breathable Material	200/cs



NONSW1765L



NONSW025M



## ComfortSure™ Scrubs

- V-neck top with 3 pockets
- Drawstring pants with pockets
- Made from premium ComfortSure™ material that provides soft, comfortable protection

Item No.	Description	Size	Pkg.
<b>ComfortSure™ Disposable Scrubs, V-Neck, Blue</b>			
NON57203XS	Scrub Pants	XS	50/cs
NON57202S	Scrub Shirt	S	50/cs
NON57203S	Scrub Pants		50/cs
NON57202M	Scrub Shirt	M	50/cs
NON57203M	Scrub Pants		50/cs
NON57202L	Scrub Shirt	L	50/cs
NON57203L	Scrub Pants		50/cs
NON57202XL	Scrub Shirt	XL	50/cs
NON57203XL	Scrub Pants		50/cs
NON57202XXL	Scrub Shirt	XXL	50/cs
NON57203XXL	Scrub Pants		50/cs



NON57202M

NON27212L

NON37202M

## Traditional Scrubs

- Choice of v-neck or round neck top with chest pocket
- Choice of elastic waist or drawstring pants with two pockets
- Made from multi-ply material

### Scrub Top, Round Neck, Blue

NON27212S	Small	30/cs
NON27212M	Medium	30/cs
NON27212L	Large	30/cs
NON27212XL	X-Large	30/cs
NON27212XXL	XX-Large	30/cs
NON27212XXXL	XXX-Large	30/cs

### Scrub Top, V-Neck, Blue

NON27202S	Small	30/cs
NON27202M	Medium	30/cs
NON27202L	Large	30/cs
NON27202XL	X-Large	30/cs
NON27202XXL	XX-Large	30/cs

### Scrub Pant, Elastic Waist, Blue

NON27213S	Small	30/cs
NON27213M	Medium	30/cs
NON27213L	Large	30/cs
NON27213XL	X-Large	30/cs
NON27213XXL	XX-Large	30/cs
NON27213XXXL	XXX-Large	30/cs

### Scrub Pant, Drawstring Waist, Blue

NON27203S	Small	30/cs
NON27203M	Medium	30/cs
NON27203L	Large	30/cs
NON27203XL	X-Large	30/cs
NON27203XXL	XX-Large	30/cs

### Scrub Top, Round Neck, Wine Red

NON37212S	Small	30/cs
NON37212M	Medium	30/cs
NON37212L	Large	30/cs
NON37212XL	X-Large	30/cs
NON37212XXL	XX-Large	30/cs
NON37212XXXL	XXX-Large	30/cs

### Scrub Top, V-Neck, Wine Red

NON37202S	Small	30/cs
NON37202M	Medium	30/cs
NON37202L	Large	30/cs
NON37202XL	X-Large	30/cs
NON37202XXL	XX-Large	30/cs

### Scrub Pant, Elastic Waist, Wine Red

NON37213S	Small	30/cs
NON37213M	Medium	30/cs
NON37213L	Large	30/cs
NON37213XL	X-Large	30/cs
NON37213XXL	XX-Large	30/cs
NON37213XXXL	XXX-Large	30/cs

### Scrub Pant, Drawstring Waist, Wine Red

NON37203S	Small	30/cs
NON37203M	Medium	30/cs
NON37203L	Large	30/cs
NON37203XL	X-Large	30/cs

### Multi-Layer Patient Gown

- Disposable material has a fabric-like feel for comfort
- Opening in front or back for easy access
- Attached ties for a custom fit

Item No.	Description	Pkg.
NON27046	Gown, Sleeveless, Waist Tie, Regular/Large, Blue	50/cs
NON27048	Gown, Sleeveless, Waist Tie, Regular/Large, Pink	50/cs
NON27146	Gown, Sleeveless, Neck and Waist Tie, Regular/Large, Blue	50/cs
NON27146SL	Gown, Short Sleeves, Neck and Waist Tie, Regular/Large, Blue	50/cs
NON27146SLXL	Gown, Short Sleeves, Neck and Waist Tie, XL, Blue	50/cs
NON27346SL	3 Arm Hole Gown, Sleeveless, Regular/Large, Blue	50/cs



### Reinforced Tissue Patient Gown

- 3-ply tissue-poly-tissue (TPT) material aids in fluid absorption and protection
- Scrim material provides added strength and absorption
- Opening in front or back for easy access
- Adjustable ties for a custom fit

Item No.	Description	Material	Pkg.
NON24356	Gown, Waist Tie, 30" x 42", Blue	TPT	50/cs
NON24355	Gown, Waist Tie, 30" x 42", White	TPT	50/cs
NON24252	Gown, Waist Tie, 30" x 42", Mauve	TPT	50/cs
NON24354	X-Ray Gown, Waist Tie, 30" x 42", Blue	TP	50/cs
NON24247	Gown, Waist Tie, 30" x 42", Blue	Scrim	50/cs
NON24247B	Bariatric Gown, Waist Tie, 72" x 45", Blue	Scrim	25/cs



### Disposable Exam Capes

- Modest design gives patients privacy
- 3-ply tissue-poly-tissue material aids in fluid absorption and protection
- Elite multi-layer material provides additional comfort and protection
- Opening in front or back for easy access

Item No.	Description	Pkg.
NON24248	Tissue-Poly-Tissue, 30" x 21", White	100/cs
NON24259	Tissue-Poly-Tissue, 30" x 21", Mauve	100/cs
NON25249	Tissue-Poly-Tissue, 30" x 21", Blue	100/cs
NON27036	Multi-Layer Material, Sleeveless, Tab Tape Closure, 30" x 21", Blue	100/cs
NON27038	Multi-Layer Material, Sleeveless, Tab Tape Closure, 30" x 21", Pink	100/cs



### Multilayer Disposable Exam Shorts

- Disposable exam shorts are available in multilayer SMS material for comfort and fluid resistance
- Elastic waist provides a snug fit

Item No.	Description	Pkg.
NON27209S	Small	30/cs
NON27209M	Medium	30/cs
NON27209L	Large	30/cs
NON27209XL	X-Large	30/cs
NON27209XXL	XX-Large	30/cs



### Reinforced Tissue Pediatric Gown

- Made of tissue-poly-tissue (TPT) material
- Disposable

Item No.	Description	Pkg.
NON25363	Blue, 21" x 36"	50/cs



### Disposable Pediatric Patient Gown

- Made from soft, polyester material
- Short sleeve with hook-and-loop closure in back

Item No.	Description	Pkg.
NON28266	Purple/Green, 6-12 Months	100/cs
NON28267	Purple/Green, 1-4 Years	50/cs
NON28268	Purple/Green, 5-8 Years	50/cs
NON28269	Blue, 9-12 Years	50/cs



# Medline PPE Isolation Cart

- Arrives fully assembled and can be rolled to the point of care to streamline caregiving
- At-a-glance supplies of face shields, ear loop masks, N95 respirators and more
- Heavy-duty steel foldable tabletop for space to write on or additional storage
- Up to 3 sizes of gloves can be faced with additional storage behind
- Room to store 50+ gowns, dispensed one at a time



IPC100-4DU

# Disposable Isolation Station

- Keeps hallways free from congestion and PPE close on hand
- Organizer stores flat-packed isolation gowns, three gloves boxes and one box of face masks
- 3 rows of pockets with eyelets for hanging
- Made of SMS material



NONISOCADDY

**Know you're protected.** Contact your Medline representative for more information.



**We make  
healthcare  
run better™**

**Medline Industries, LP**  
Three Lakes Drive, Northfield, IL 60093  
Medline United States | 1-800-MEDLINE (633-5463)  
medline.com | info@medline.com

**Medline Canada**  
1-800-268-2848 | medline.ca | canada@medline.com  
**Medline México**  
01-800-831-0898 | medlinemexico.com | mexico@medline.com

Follow us    

#### REFERENCE:

1. Considerations for selecting protective clothing used in healthcare for protection against microorganisms in blood and body fluids. Centers for Disease Control and Prevention website. Table 3: ANSI/AAMI PB 70:12 Classification of barrier performance of surgical gowns, other protective apparel, surgical drapes and drape accessories. Available at: <http://www.cdc.gov/niosh/nppt/topics/protectiveclothing/>. Accessed May 3, 2016.

Some products may not be available for sale in Mexico or Canada. We reserve the right to correct any errors that may occur within this brochure.

© 2023 Medline Industries, LP. Companion and QuickComply are trademarks and Medline is a registered trademark of Medline Industries, LP. MKT19W4555184 / LIT171 / 3