

REVOLUTION XL

REVOLUTION SERIES

ITEM# F-REVXL

U.S. Patent No. D507,769S



ABOUT

Now using a direct drive motor, the new and improved Revolution XL is stronger and more sleek than ever before! Other upgrades include an increased anchor setback, as well as an increased water draft. For versatility that doesn't compromise strength, you can't beat the Revolution XL.



FEATURES:

- ADA Compliant & UL Safety Certified
- 500* lb (226 kg) weight capacity
- 360 degree rotation
- Adjustable seat height and foot rest
- Wider, more comfortable 22" seat
- Stainless steel with powder coat finish
- Limited five (5) year warranty
- Ships common carrier, some assembly required
- Made in USA

*350 lb (158 kg) capacity with spine board attachment



STANDARD COLOR OPTIONS:

PLASTICS



GRANITE
WHITE



GRANITE
TAN



GRANITE
GREY



GREY



BLUE



TAN

POWDER
COAT



GREY



WHITE



TAN

PROUDLY SOLD BY:

Access, therapy and fitness equipment for people of all abilities.

1.888.687.3552 | aquacreek.com | sales@aquacreek.com

REVOLUTION XL


APPLICATIONS:

This lift works on in ground and some partially raised pools and spas. To see if this lift works with your pool, visit our website and download the application guide for this product, or contact our friendly sales team. This lift can also be anchored to pool decks made of concrete, pavers and even wood. Choose the appropriate anchor system for your needs. Anchors sold separately.

Standard Anchor
F-808SA




Semi-Recessed Anchor
F-77SSA



Paver Anchor
F-808SA-10



Wood Deck Anchor
F-842SWDA



For access to CADs, Revit drawings and more visit our Partner Login page on our website and become an Aqua Creek partner.

OPTIONS & ACCESSORIES:

- Anchor options for concrete, paver, and wood deck applications
- Anchor adapters for many existing pool decks anchors
- Chest strap for additional upper body support
- Cover to protect your investment from the elements*
- Cycle attachment to convert to an exercise bike
- Headrest for comfortable head support
- Pull out leg rest for additional wall clearances
- Rescue spine board option
- Solar charging station to conserve energy
- Transport cart for easy removal from pool side for storage**

*Color and material options available

**Transport cart not usable with semi recessed anchor


SCAN ME



Use your mobile device to scan the code and view the product webpage.



Match your branding or decor with the custom color of your choice. Additional charges may apply.



Solar chargers can help extend the life of your battery and reduce the need to remove it for charging.



Protect your investment from the elements with a premium cover.