

SideStream reusable and disposable nebulizers



Venturi system

provides a faster delivery rate resulting in shorter treatment times than other conventional nebulizers

Fast treatment times

- Fast treatment times (delivers 2.5 ml of drug output within 7 minutes compared to up to 15 minutes for other conventional nebulizers)

Data on file at Philips Respirationics.

Efficient drug delivery

- Highly efficient (80% of drug output volume below 5 microns in size)*
- The SideStream nebulizer's passive venturi system provides air flow in addition to that supplied by the compressor, resulting in faster drug delivery and shorter treatment times
- Low residual volume
- Fewer parts to assemble or clean than other conventional nebulizers
- Provides consistent delivery of medication time after time

* Ref: Collis G.G. et al. Lancet 1990:336:341



Detachable masks are available separately in two sizes, adult and pediatric, including Tucker the Turtle pediatric mask.

Economical

	HCPSC code
Reusable SideStream	A7005
Disposable SideStream	A7003
SideStream masks	A7015
InnoSpire compressor nebulizer systems	E0570

- Reusable SideStream is recommended for up to 6 months of use
- Disposable SideStream is recommended for up to 2 weeks of use

Ordering information

Product description	Reorder number
SideStream disposable nebulizer	HS800
SideStream reusable nebulizer	HS860
SideStream pediatric mask	1025529
SideStream Tucker the Turtle pediatric mask	1044382
SideStream adult mask	1025531
InnoSpire Elegance compressor nebulizer system with SideStream disposable and reusable nebulizers	1099969
InnoSpire Elegance compressor nebulizer system with SideStream disposable nebulizer	1100134
InnoSpire Essence compressor nebulizer system with SideStream disposable and reusable nebulizers	1100312
InnoSpire Essence compressor nebulizer system with SideStream disposable nebulizer	1099966

Respirationics, InnoSpire, Tucker the Turtle and SideStream are trademarks of Koninklijke Philips N.V. and its affiliates. All rights reserved.



© 2014 Koninklijke Philips N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician or licensed healthcare professional.

Geyer JH 03/31/14 MCI 4106092 PN 1116348

Philips Healthcare is part of Royal Philips

Philips Respirationics Asia Pacific
+65 6882 5282

Philips Respirationics Australia
+61 (2) 9947 0440

Philips Respirationics China
+86 400 828 6665

Philips Respirationics United Kingdom
+44 800 1300 845

Philips Respirationics
1010 Murry Ridge Lane
Murrysville, PA 15668

Customer Service
+1 724 387 4000
800 345 6443 (toll free, US only)
www.philips.com/respirationics