



OR Products  
1-800-962-9888





## Table of Contents

### OR Suction

Argyle™ Yankauer Suction Instruments, Rigid Type .....	2
Argyle™ Yankauer Suction Instruments, Flexible Type .....	3
Argyle™ Poole and Suction Sets .....	4
Curity™ Suction Instruments and Sets .....	4
Argyle™ Suction Tubing .....	5
Argyle™ Universal Bubble Tubing .....	6
Argyle™ Tubing Connectors .....	6
Solidification Products .....	7

### OR Essentials

Devon™ Light Handle Systems .....	9
Devon™ Needle Counters .....	10
Devon™ Surgical Markers .....	11
Devon™ Utility Markers and Labels .....	12
OR Towels .....	13
Laparotomy Sponges .....	13
Devon™ Fog-Out Solution .....	14
Devon™ Instrument & Cautery Tip Cleaners .....	14
Devon™ Drape Clamps .....	14
Devon™ Surgical Vessel Loops .....	15
Devon™ Tube Holder .....	15
Devon™ Instrument Pouch .....	16
Devon™ Scope Adjustment Bag .....	16
Devon™ Suture Bag .....	16
Curity™ Disposable Scalpels .....	16

### OR Safety

Sponge Counting Systems .....	18
Solidification Products .....	18
Devon™ Hands Free Transfer Magnetic Drapes .....	19
Devon™ Magnetic Drapes with Neutral Zones .....	19
Devon™ Trays and Basins .....	19
Devon™ Reusable Magnetic Drapes .....	20
Devon™ Blade Shield Scalpel Holder .....	20
Devon™ Magnetic Needle Finder .....	21
Devon™ Safe-T-Strainer Instrument Handling Products .....	21

### Patient Prep

Merlin™ One Step Prep Applicator .....	23
Curity™ Sterile Wet Skin Scrub Prep Products .....	23
Curity™ Sterile Vaginal Prep, Wet Skin Scrub Products .....	24
Curity™ Sterile Dry Skin Scrub Prep Products .....	25

## OR Suction




### Argyle™ Yankauer Suction Instruments, Rigid Type

- Transparent with smooth tip and eyes
- Sterile; single-use
- Latex Free
- DEHP Free

### Argyle™ Rigid Yankauer Suction Instruments

Item Code	Description	Qty/Case
8888505016	Bulbous tip, without vent	50
8888505024	Bulbous tip, with vent	50
8888505115	Open tip, without vent	50
8888505123	Open tip, with vent	50



Bulbous Tip




Open Tip

### Argyle™ Yankauer Selec-Trol™, Rigid Type

- Sliding switch on handle turns vacuum on and off
- Smooth tip and eyes; sterile; single-use
- Ribbed, easy-to-grip handle
- Latex Free
- DEHP Free

### Argyle™ Yankauer Selec-Trol™, Rigid Type

Item Code	Description	Qty/Case
8888505164	Selec-Trol Yankauer, bulb tip	20
8888505172	Selec-Trol Yankauer, open tip	20



Bulbous Tip



Open Tip

### Argyle™ Rigid Yankauer Suction Instruments with Non-Conductive Connective Tubing (molded connectors)


- Transparent with smooth tip and eyes
- Sterile; single-use
- Latex Free

Item Code	Description	Qty/Case
8888505032	Bulb tip, without vent and ¼" x 6' (6mm x 1.8m) tube	50
8888505040	Bulb tip, without vent and ¼" x 10' (6mm x 3.1m) tube	20



Bulbous Tip

## OR Suction



Regular Capacity, Open Tip



Regular Capacity, Sump Tip



High Capacity, Open Tip



Fine Capacity, Open Tip

### Argyle™ Yankauer Suction Instruments, Flexible Type

- Nonbreakable, malleable vinyl construction with smooth inner lumen
- Optional Tip Trof™ vent for fingertip control of suction
- Sterile; single-use
- Latex Free
- DEHP Free

#### Argyle™ Flexible Yankauer Suction Instruments

Item Code	Description	Qty/Case
8888501007	Regular tip capacity, without vent	50
8888501015	Sump tip, regular tip capacity, without vent	50
8888501023	Regular tip capacity, with vent	50
8888501031	Sump tip, regular tip capacity, with vent	50
8888502005	High capacity tip, without vent	50
8888504001	Fine capacity tip, without vent	50

### Argyle™ Yankauer Suction Instruments, Flexible Type with Non-Conductive Connective Tubing

- Nonbreakable, malleable vinyl construction with smooth inner lumen
- Optional Tip Trof™ vent for fingertip control of suction
- Sterile; single-use
- Latex Free

#### Argyle™ Flexible Yankauer Suction Instruments

##### (regular tip capacity, without vent) with Non-Conductive Tubing

Item Code	Description	Qty/Case
8888501619	Includes 1/4" x 6' (6mm x 1.8m) tubing	50
8888501742	Includes 1/4" x 12' (6mm x 3.7m) tubing	20

#### Argyle™ Sump Tip Flexible Yankauer Suction Instruments

##### (regular tip capacity, without vent) with Non-Conductive Tubing

Item Code	Description	Qty/Case
8888501635	Includes 1/4" x 6' (6mm x 1.8m) tubing	50
8888501775	Includes 1/4" x 12' (6mm x 3.7m) tubing	20

#### Argyle™ Flexible Yankauer Suction Instruments

##### (regular tip capacity, without vent) with Conductive Tubing

Item Code	Description	Qty/Case
8888501411	Includes 3/16" x 6' (5mm x 1.8m) tubing	50
8888501429	Sump tip, includes 3/16" x 6' (5mm x 1.8m) tubing	50

## OR Suction



### Argyle™ Poole and Suction Sets

- Two-part tube locks for use as an assembled unit
- Rigid outer sleeve has 108 smooth eyes with one sump eye at proximal end
- Sterile; single-use
- Latex Free

### Argyle™ Poole and Suction Sets

Item Code	Description	Qty/Case
8888509703	Poole with vent, DEHP Free	50
8888509737	Poole with vent & ¼" x 6' (6mm x 1.8m) tubing	50
8888509745	Poole with vent & ¼" x 12' (6mm x 3.7m) tubing	20



### Curity™ Suction Instruments and Sets

- Sigmoids with semi-rigid construction, transparent for visual monitoring
- Available with suction tubing
- Single-use; sterile
- Latex Free

### Curity™ Sigmoid Instruments and Sets

Item Code	Description	Qty/Case
166030	Sigmoid, 18 French, DEHP Free	40
166035	Sigmoid, 23 French, DEHP Free	40
166031	Sigmoid, 18 French & ¼" x 6' (6mm x 1.8m) tubing	20

### Curity™ Frazier Instruments

Item Code	Description	Qty/Box	Qty/Case
166024	Frazier, 8 French, DEHP Free	10	40
166025	Frazier, 10 French, DEHP Free	10	40
166026	Frazier, 12 French, DEHP Free	10	40

## OR Suction

### Argyle™ Suction Tubing

- Available with female molded connectors (male connector also included) or integral connectors
- Package has a sterile inner pouch
- Non-conductive type only
- Sterile;
- Latex Free



### Argyle™ Suction Tubing with Molded Connectors

Item Code	Description	Qty/Case
8888301507	$\frac{3}{16}$ " x 1.5' (5mm x .5m)	50
8888301515	$\frac{3}{16}$ " x 6' (5mm x 1.8m)	50
8888301509	$\frac{3}{16}$ " x 9' (5mm x 2.7m)	50
8888301523	$\frac{3}{16}$ " x 10' (5mm x 3.1m)	50
8888301531	$\frac{3}{16}$ " x 12' (5mm x 3.7m)	20
8888301549	$\frac{3}{16}$ " x 20' (5mm x 6.1m)	10
8888301606	$\frac{1}{4}$ " x 6' (6mm x 1.8m)	50
8888301609	$\frac{1}{4}$ " x 9' (6mm x 2.7m)	50
8888301614	$\frac{1}{4}$ " x 10' (6mm x 3.1m)	50
8888301622	$\frac{1}{4}$ " x 12' (6mm x 3.7m)	20
8888301630	$\frac{1}{4}$ " x 20' (6mm x 6.1m)	10
8888301645	$\frac{1}{4}$ " x 25' (6mm x 7.6m)	10
8888301705	$\frac{9}{32}$ " x 6' (7mm x 1.8m)	50
8888301713	$\frac{9}{32}$ " x 10' (7mm x 3.1m)	50
8888301721	$\frac{9}{32}$ " x 12' (7mm x 3.7m)	20
8888301739	$\frac{9}{32}$ " x 20' (7mm x 6.1m)	10



### Argyle™ Suction Tubing with Funnel/Funnel Integral Connectors

Item Code	Description	Qty/Case
8888284505	$\frac{3}{16}$ " x 1.5' (5mm x .5m)	50
8888284513	$\frac{3}{16}$ " x 6' (5mm x 1.8m)	50
8888284539	$\frac{3}{16}$ " x 10' (5mm x 3.1m)	50
8888284547	$\frac{3}{16}$ " x 12' (5mm x 3.7m)	20
8888284554	$\frac{3}{16}$ " x 20' (5mm x 6.1m)	10
8888284604	$\frac{1}{4}$ " x 6' (6mm x 1.8m)	50
8888284612	$\frac{1}{4}$ " x 10' (6mm x 3.1m)	50
8888284620	$\frac{1}{4}$ " x 12' (6mm x 3.7m)	20
8888284638	$\frac{1}{4}$ " x 20' (6mm x 6.1m)	10
8888284703	$\frac{9}{32}$ " x 6' (7mm x 1.8m)	50
8888284737	$\frac{9}{32}$ " x 12' (7mm x 3.7m)	20
8888284745	$\frac{9}{32}$ " x 20' (7mm x 6.1m)	10

### Argyle™ Suction Tubing with Funnel/Tapered Integral Connectors

Item Code	Description	Qty/Case
8888284521	$\frac{3}{16}$ " x 6' (5mm x 1.8m)	50
8888284711	$\frac{9}{32}$ " x 6' (7mm x 1.8m)	50

## OR Suction



### Argyle™ Universal Bubble Tubing

- Clean bulk tubing that can be cut to the exact length needed
- Conductive type meets N.F.P.A. specifications for use with flammable agents
- "Bubble" every 3 feet
- Latex Free

### Argyle™ Non-Conductive Bubble Tubing

Item Code	Description	Qty/Case
8888280214	3/16" x 100' (5mm x 30.5m)	1
8888280412	1/4" x 100' (6mm x 30.5m)	1
8888280610	9/32" x 100' (7mm x 30.5m)	1

### Argyle™ Conductive Bubble Tubing

Item Code	Description	Qty/Case
8888290213	3/16" x 100' (5mm x 30.5m)	1
8888290411	1/4" x 100' (6mm x 30.5m)	1
8888290619	9/32" x 100' (7mm x 30.5m)	1

### Argyle™ Tubing Connectors

- Wide range of types and sizes
- Available sterile or clean



### Argyle™ Tubing Connectors

Item Code	Type	Fits Tubing Diameter	Qty/Case
<b>Sterile - Individually Packaged, Polyethylene</b>			
8888270603	5-in-1 (bubble)	3/16" - 7/16" (4.8-11.1mm)	100
8888271411	5-in-1 (barbed)	5/16" - 9/16" (7.9-14.3mm)	100
8888271502	5-in-1 (barbed)	3/16" - 9/16" (4.8-14.3mm)	100
8888272013	Large Y (barbed)	1/4" - 1/2" (6.4-12.7mm)	50
8888275008	Female Luer Adapter	Small catheter	100
<b>Sterile - Individually Packaged, Vinyl</b>			
8888270207	Sims	3/16" - 7/16" (4.8-11.1mm)	100
8888271007	5-in-1 (bubble)	3/16" - 7/16" (4.8-11.1mm)	100
<b>Clean - Bulk Packaged, Polyethylene</b>			
8888270702	5-in-1 (bubble)	3/16" - 7/16" (4.8-11.1mm)	100
8888271403	5-in-1 (barbed)	5/16" - 9/16" (7.9-14.3mm)	100
8888271510	5-in-1 (barbed)	3/16" - 9/16" (4.8-14.3mm)	100
8888272005	Large Y (barbed)	1/4" - 1/2" (6.4-12.7mm)	50
<b>Clean - Bulk Packaged, Vinyl</b>			
8888270306	Sims	3/16" - 7/16" (4.8-11.1mm)	100
8888271106	5-in-1 (bubble)	3/16" - 7/16" (4.8-11.1mm)	100



## OR Suction



### Solidification Products

- Designed to absorb and encapsulate blood, vomit or any fluid with at least 6% water content
- Solidifiers only
- Latex Free

#### Covidien™ Solidification Products Cone Cap

Item Code	Description	Qty/Case
31161006	1,200 ml single dose bottle	75
31146163	1,500 ml single dose bottle	75
31146148	2,000 ml single dose bottle	75
31163473	3,000 ml single dose bottle	75
31146130	16,000 ml multi dose bottle	20

#### Devon™ Solidification Products Flat Cap

Item Code	Description	Qty/Case
54901200	1,200 ml single dose bottle	75
54901500	1,500 ml single dose bottle	75
54902000	2,000 ml single dose bottle	75
54903000	3,000 ml single dose bottle	75
54916000	16,000 ml multi dose bottle	20

1-800-962-9888



## OR Essentials

### Devon™ Light Handle Systems

- Select one system to standardize all of your overhead lighting needs
- Convenient and cost effective
- Disposable, sterile products eliminate time and labor of reprocessing
- Latex Free



#### Devon™ Light Gloves

Item Code	Description	Qty/Case
31140208	1/pack, 25/box	100
31140257	2/pack, 40/box	160
31140216	3/pack, 48/box	192

#### Devon™ Light Handle (Rigid)

Item Code	Description	Qty/Case
31140133	1/pack, 16/box	128
31140125	2/pack, 28/box	112
31140117	3/pack, 21/box	84

#### Devon™ Light Sleeve (Semi-Rigid)

Item Code	Description	Qty/Case
31140273	1/pack	70
31140240	2/pack	120

### Devon™ Universal Light Handles (ULH) and Universal Light Adapters (ULA)

- ULH is designed for use with the Light Glove or the Light Sleeve
- ULA is designed for use with the Light Handle



#### Devon™ Universal Light Handles (ULH) and Universal Light Adapters (ULA)

Available for direct order through your Covidien Representative

OR Essentials

**Devon™ Needle Counters**

- Available in foam block, foam strip and magnetic combinations with needle counts from 10-100
- BoxLocks feature provides secure closure and disposal, but may be reopened as needed for recounts
- Blade Removal System available
- Latex Free

**Devon™ Foam Block Needle Counters**

**(with BoxLocks feature) includes Blade Removal**

Item Code	Description	Qty/Box	Qty/Case
31142352	100 count/100 capacity, double black magnet	8	64
31142378	50 count/50 capacity, double black magnet	12	96
31142360	20 count/40 capacity, adhesive backing	12	96
31142386	20 count/40 capacity, single black magnet	12	96
31142394	10 count/20 capacity, single black magnet	12	96

**Devon™ Foam Block Needle Counters (with BoxLocks feature) no Blade Removal**

Item Code	Description	Qty/Box	Qty/Case
31142469	80 count/80 capacity, double foam	8	64
31142402	40 count/70 capacity, adhesive backing	8	64
31142311	40 count/70 capacity, single black magnet	8	64
31142501	20 count/40 capacity, adhesive backing	12	96
31142493	20 count/40 capacity, single black magnet	12	96
31142485	10 count/20 capacity, adhesive backing	12	96
31142428	10 count/20 capacity, single black magnet	12	96

**Devon™ Foam Strip Needle Counters**

**(with BoxLocks feature) includes Blade Removal**


Item Code	Description	Qty/Box	Qty/Case
31142451	30 count/30 capacity, double foam	12	96
31142295	15 count/30 capacity, black Magnet	12	96

**Devon™ Foam Strip Needle Counters (with BoxLocks feature) no Blade Removal**

Item Code	Description	Qty/Box	Qty/Case
31142204	60 count/60 capacity, double foam	8	64
31142303	30 count/30 capacity, double foam	12	96
31142188	30 count/60 capacity, black magnet	8	64
31142170	15 count/30 capacity, black magnet	12	96
31142212	15 count/15 capacity (clear top)	12	96



## OR Essentials



### Devon™ Magnetic Needle Counters (with BoxLocks feature) includes Blade Removal

Item Code	Description	Qty/Box	Qty/Case
31142238	30 count/30 capacity, double white magnet	12	96
31142113	30 count/30 capacity, double black magnet	12	96

### Devon™ Magnetic Needle Counters (with BoxLocks feature) no Blade Removal

Item Code	Description	Qty/Box	Qty/Case
31142154	60 count/60 capacity, double black magnet	8	64
31142147	30 count/30 capacity, double black magnet	12	96
31142139	30 count/30 capacity, single black magnet (clear top)	8	64

### Devon™ Magnetic Needle Counters (no BoxLocks feature) no Blade Removal

Item Code	Description	Qty/Box	Qty/Case
31142196	20 count/20 capacity, double black magnet	21	168
31142410	10 count/10 capacity, single black magnet (clear top)	21	168



## Devon™ Surgical Markers

- Non-toxic, non-smearing, non-irritating Gentian violet ink
- Choice of regular, fine or dual tip markers
- Sterile labels and/or flexible rulers included with markers as indicated
- Latex Free

### Devon™ Surgical Markers

Item Code	Description	Qty/Box	Qty/Case
31145819	Regular tip	25	100
31145785	Regular tip, flexible ruler	25	100
31145777	Regular tip, flexible ruler, 9 labels	25	100
31177366	Fine tip	25	100
31145827	Fine tip, flexible ruler	25	100
31145793	Fine tip, flexible ruler, 9 labels	25	100



## OR Essentials



### Devon™ Surgical Markers with Ruler Cap

Item Code	Description	Qty/Box	Qty/Case
31145926	Regular tip	25	100
31145918	Regular tip, 9 labels	25	100
31145942	Fine tip	25	100
31145892	Fine tip, flexible ruler	25	100
31145843	Fine tip, flexible ruler, 9 labels	25	100
31145959	Regular tip (*Non-Sterile item)	25	100
31145900	Dual tip	25	100
31145868	Dual tip, 9 labels	25	100
31160560	Dual tip, flexible ruler	25	100
31145934	Dual tip, flexible ruler, 9 labels	25	100

### Devon™ Utility Markers and Labels

#### *Sterile Marker*

- Sterile, single-use utility marker
- Nine labels included
- Non-smearing, fast drying, non-toxic ink
- Latex Free

#### *Non-Sterile Marker*

- Non-sterile, multi-use utility marker
- May be used on most surfaces including non-woven CSR wrap, sterilization tape, plastic, glass or stainless steel
- Non-toxic ink may be sterilized via steam or gas
- Latex Free



### Devon™ Utility Markers

Item Code	Description	Qty/Box	Qty/Case
31145884	Sterile utility marker with 9 labels	25	100
51000177	Non-sterile utility marker	25	50

### Devon™ Utility Labels

Item Code	Description	Qty/Case
31145751	9 labels (25 pouches of 9 per box)	900
31145769	Bulk packaged labels	3,000



**OR Essentials**



**OR Towels**

- Designed for maximum absorbency
- Manufactured from high quality 100% cotton
- Pre-washed, sterile and color fast
- Latex Free

**OR Towels**

Item Code	Description	Qty/Case
77702	Blue OR towel, 2 pack	80
77704	Blue OR towel, 4 pack	80
77706	Blue OR towel, 6 pack	72
77708	Blue OR towel, 8 pack	80



**Laparotomy Sponges, X-ray Detectable**

- Sterile, pre-washed 100% cotton, disposable
- Kerlix™ sponges have dual-colored radiopaque handles that remain visible when saturated with blood
- Latex Free

**Curity™ Lap Sponges — Sterile Tray Pack**

Item Code	Description	Qty/Case
6022	18" x 18"	100
6034	12" x 12"	100
6045	4" x 18"	100
6056	8" x 36"	100



**Curity™ Lap Sponges — Sterile Soft Pouch Pack**

Item Code	Description	Qty/Case
6522	18" x 18"	200
6534	12" x 12"	200
6545	4" x 18"	200
6556	8" x 36"	200

**Curity™ Lap Sponges — Non-Sterile Bulk Pack**

Item Code	Description	Qty/Case
6822	18" x 18"	400
6834	12" x 12"	750

**Kerlix™ Lap Sponges — Sterile Tray Pack, Made in the USA**

Item Code	Description	Qty/Case
4031	17" x 26"	100
4034	16" x 17"	100
4037	8" x 32"	100
4048	3" x 17"	100
4069	8" x 11"	100



## OR Essentials




### Devon™ Fog-Out Solution

- Anti-fog solution is sterile, non-toxic, laboratory tested and safe to use in a body cavity
- Inhibits condensation and clouding
- May be used on endoscopic lenses, surgical cameras, protective eyewear, goggles or glasses
- Latex Free

### Devon™ Fog-Out Solution

Item Code	Description	Qty/Box	Qty/Case
31142527	3910 Anti-fog solution with foam pad	12	48
31142808	3919 Anti-fog solution with #1910 needle counter	6	48



31142790  
Cautery Tip Polisher

### Devon™ Instrument & Cautery Tip Cleaners

- Cautery Tip Polisher — Safely and effectively cleans electro-surgical pencil tips; strong adhesive backing; X-ray detectable
- Micro Tip Wipe — A sterile, lint and particulate-free foam wipe for cleaning the fine tips of micro instruments and lenses
- Latex Free



31142535  
Micro Tip Wipe

### Devon™ Instrument & Cautery Tip Cleaners

Item Code	Description	Qty/Box	Qty/Case
31142790	3200 cautery tip polisher (5cm x 5cm)	25	100
31142717	3204 cautery tip polisher (5cm x 5cm)	80	240
31142535	3300 micro tip wipe	25	100



31140315  
Tube  
Occluding  
Clamp

### Devon™ Drape Clamps

- Select from a metal reusable design or sterile or non-sterile disposable clamp
- Non-sterile, tube occluding, disposable clamp will not cut fine tubes

31140331  
Reusable  
Drape Clamp



31140299  
Disposable  
Drape Clamp

### Devon™ Disposable Drape Clamps

Item Code	Description	Qty/Box	Qty/Case
31140315	Non-sterile tube occluding clamp, bulk pack		50
31140299	Non-sterile drape clamp, bulk pack		100
31140281	Sterile drape clamp, 1 pack	10	80
31140307	Sterile drape clamp, 2 pack	15	120

### Devon™ Reusable Drape Clamps

Item Code	Description	Qty/Case
31140331	Non-sterile reusable drape clamp	10



## OR Essentials



Red=Arteries  
Blue=Veins  
White=Nerves  
Yellow=Ureters

### Devon™ Surgical Vessel Loops

- Silicone, radiopaque
- Available in four colors and two sizes
- Latex Free

### Devon™ Mini Surgical Vessel Loops

- Cross Section 1.3mm x 0.9mm
- Length 17 inches

Item Code	Description	Qty/Box	Qty/Case
31145728	White	20	200
31145710	Red	20	200
31145702	Blue	20	200
31145660	Yellow	20	200

### Devon™ Maxi Surgical Vessel Loops

- Cross Section 2.4mm x 1.15mm
- Length 17 inches

Item Code	Description	Qty/Box	Qty/Case
31145686	White	20	200
31145744	Red	20	200
31145678	Blue	20	200
31145694	Yellow	20	200



### Devon™ Tube Holder

- Organizes and securely holds tubes and cords in place
- Available with removable or fixed touch fasten strap
- Sterile, single-use
- Latex Free

### Devon™ Tube Holder

Item Code	Description	Qty/Box	Qty/Case
31142667	Tube holder, removable fasten strap	25	100
31142659	Tube holder, removable fasten strap	12	96
31142675	Tube holder, fixed fasten strap	25	100

**OR Essentials**



**Devon™ Instrument Pouch**

- Stores a variety of instruments connected to tubes and cords
- Latex Free

**Devon™ Instrument Pouch**

Item Code	Description	Qty/Box	Qty/Case
31142550	3221 instrument pouch (double pocket)	20	80



**Devon™ Scope Adjustment Bag**

- Sterile microscope cover allows surgeon to adjust microscope during procedures
- Latex Free

**Devon™ Scope Adjustment Bag**

Item Code	Description	Qty/Box	Qty/Case
31142543	3406 scope adjustment bag	120	480



**Devon™ Suture Bag**

- Adheres securely to any surface for convenient, quick access
- Latex Free

**Devon™ Suture Bag**

Item Code	Description	Qty/Box	Qty/Case
31142691	3570 suture bag	25	100




**Curity™ Disposable Scalpels**

- Sterile, surgical blades attached to plastic handles
- Individually wrapped; available loose or in boxes of 10 scalpels
- Handles have metric measurements for incision recording
- Latex Free

**Curity™ Disposable Scalpels**

Item Code	Description	Qty/Case
139090	#10, bulk packaged	100
139091	#11, bulk packaged	100
139093	#15, bulk packaged	100
131610	#10, 10 per box	10 boxes
131611	#11, 10 per box	10 boxes
131615	#15, 10 per box	10 boxes



## OR Safety

### Sponge Counting Systems

- Economical way to develop a standard sponge counting protocol
- Assists in the prevention of retained sponges
- Provides a reliable system for counting sponges during surgical procedure
- Latex Free

### Sponge Counting Systems

Item Code	Description	Qty/Box	Qty/Case
31144119	Devon™ Bag-It dispenser box, clear back	50	250
31144127	Devon™ Bag-It dispenser box, white back	50	250
5096	Devon™ Bag-It dispenser box, blue back	50	250

### Sponge Counting Racks and Brackets

Item Code	Description	Qty/Case
31144168	Devon™ rack for I.V. pole, single	1
31144184	Devon™ rack for I.V. pole, triple	1
31144176	Devon™ rack for wall, single	1
31144192	Devon™ rack for wall, triple	1
5097	Curity™ rack with I.V. pole mounting bracket	1



### Solidification Products

- Designed to absorb and encapsulate blood, vomit or any fluid with at least 6% water content
- Solidifiers only
- Latex Free

### Covidien™ Solidification Products Cone Cap

Item Code	Description	Qty/Case
31161006	1,200 ml single dose bottle	75
31146163	1,500 ml single dose bottle	75
31146148	2,000 ml single dose bottle	75
31163473	3,000 ml single dose bottle	75
31146130	16,000 ml multi dose bottle	20

### Devon™ Solidification Products Flat Cap

Item Code	Description	Qty/Case
54901200	1,200 ml single dose bottle	75
54901500	1,500 ml single dose bottle	75
54902000	2,000 ml single dose bottle	75
54903000	3,000 ml single dose bottle	75
54916000	16,000 ml multi dose bottle	20



## OR Safety

### Devon™ Hands Free Transfer Magnetic Drapes

- Disposable magnetic drapes designed to attract and hold surgical instruments on the sterile field
- Helps reduce risk of percutaneous injuries and helps enhance safe, effective instrument transfer in the OR
- Flexible drapes conform to any shape with non-slip foam backing for secure positioning
- Latex Free

### Devon™ Hands Free Transfer Magnetic Drapes

Item Code	Description	Qty/Case
31140570	10" x 16"	30
31140562	16" x 20"	20
31140596	10" x 16", with stays	30
31140547	16" x 20", with stays	20

### Devon™ Magnetic Drapes with Neutral Zones

- Neutral Zone magnetic drapes have a convenient, built-in, non-magnetic area to hold or safely transfer sharp instruments
- Latex Free

### Devon™ Magnetic Drapes with Neutral Zones

Item Code	Description	Qty/Case
31140612	10" x 16", one neutral zone	30
31140539	16" x 20", one neutral zone	20
31140505	10" x 16", two neutral zones	26
31140430	16" x 20", two neutral zones	16

### Devon™ Trays and Basins

- Kidney shape basin and rectangular tray easily attach to magnetic drapes
- Latex Free

### Devon™ Magnetic Drapes with Neutral Zones

Item Code	Description	Qty/Case
31140398	3¾" x 9" x 1" tray with adhesive strip	50
31140380	3¾" x 9" x 1" tray with instrument holder, adhesive strip	50
31140372	4" x 8¼" x 2" basin with metal feet	40



**OR Safety**



**Devon™ Reusable Magnetic Drapes**

- Constructed from high quality medical grade siliconized rubber for strength and durability
- Round ceramic magnets help prevent perforations after repeated handling and use
- May be ETO or steam sterilized
- Latex Free

**Devon™ Reusable Magnetic Drapes**

Item Code	Description	Qty/Case
31140588	14.75x11.75	3

**Devon™ Blade Shield Scalpel Holder**

- Sterile, disposable blade shield holds sharp edges of scalpels in the “down” position, helping to prevent accidental cuts or stabs
- Holds up to three scalpels in a convenient, ready-to-use position
- Adhesive tab on back secures blade shield to work surface
- Latex Free

**Devon™ Blade Shield Scalpel Holder**

Item Code	Description	Sterile	Qty/Case
31142576	3593 blade shield scalpel holder	24	96



## OR Safety



### Devon™ Magnetic Needle Finder

- Rolling, magnetic and reusable needle finder provides a convenient, safe and effective method to find and pick up needles or blades that may have fallen on the floor
- Convenient handle length of 44" for ease of use
- Latex Free

### Devon™ Magnetic Needle Finder

Item Code	Description	Qty/Case
31140422	100 magnetic needle finder	1



### Devon™ Safe-T-Strainer Instrument Handling Products

- Disposable basin liner eliminates need to reach by hand into containers filled with sharps that are contaminated with blood or other potentially infectious material
- An optional flexible lid is available for transporting filled safe-t-strainer product outside of the OR
- Easy grip handles provide a secure grip to lift heavy instrument loads
- Latex Free

### Devon™ Safe-T-Strainer Instrument Handling Products

Item Code	Description	Sterile	Qty/Case
31144234	2701 Safe-T-Strainer with flexible lid	Yes	50
31144291	2705 Safe-T-Strainer (bulk)	No	100
31144283	2706 Safe-T-Strainer w/flexible lid (bulk)	No	100
31144226	2707-B flexible lid (bulk)	No	100





## Patient Prep



### Merlin™ One Step Prep Applicator

#### Merlin™ One Step Prep Applicator

- Viscous solution reduces dripping and pooling, increasing OR safety
- Directions for use that allow for scrubbing, creating friction to remove transient bacteria from the skin
- Unique sponge head design allows for easier access to fingers and toes
- 50% more antiseptic than other applicators on the market (39ml versus 26ml)

Item Code	Description	Color	Qty/Case
605039	39 ml applicator PVP/IPA (povidone iodine 7.5% w/isopropyl alcohol 72% v/v)	Blue	20

#### Merlin™ One Step Prep Applicator (Alcohol Free)

- Alcohol Free (Povidone iodine 10%) in an easy to use applicator

Item Code	Description	Color	Qty/Case
606139	39 ml applicator alcohol-free (povidone iodine 10%)	Purple	20

### Curity™ Sterile Wet Skin Scrub Prep Products

- Your choice of antimicrobials — PVP-I solution, PVP-I gel
- One-step prep PVP-I gel packs combine scrubbing and painting in one simple, convenient step; no alcohol

#### Curity™ Sterile Wet Skin Scrub Trays with PVP-I Solutions

*Wet Skin Scrub Trays include:*

1 CSR wrap	1 Removable tray
1 Tray (3 compartments)	2 Cotton tipped applicators
2 Blue blotting towels	2 White absorbent towels
1 Pair nitrile latex-free, powder-free gloves	1 PVP-I Scrub solution (3.7 oz. bottle)
1 PVP-I Paint solution (3 oz. bottle)	3 Sponge sticks



Item Code	Description	Qty/Case
41589	Unwrapped gloves, 6 small winged grip sponges	20
41590	Wrapped gloves, 2 large winged grip sponges	20
41591	Wrapped gloves, 6 small winged grip sponges	20

## Patient Prep



### Curity™ Sterile Wet Skin Scrub Packs with PVP-I Solutions

*Wet Skin Scrub Packs include:*

- 1 CSR wrap
- 2 Cotton tipped applicators
- 1 Pair nitrile latex-free, powder-free gloves
- 2 Blue blotting towels
- 2 White absorbent towels

Item Code	Description	Qty/Case
41529	Unwrapped gloves, 2 sponge sticks*, 2 large winged grip sponges†	20
41530	Wrapped gloves, 2 sponge sticks*, 2 large winged grip sponges†	20
41531	Wrapped gloves, 2 sponge sticks*, 2 small winged grip sponges†	20
41535	Wrapped gloves, 2 large winged grip sponges†	20

\* *Sponge sticks contain PVP-I Paint Solution*

† *Winged grip sponges contain PVP-I Scrub Solution*



### Curity™ Sterile Wet Skin Scrub Packs with One Step Prep PVP-I Gel

*Wet Skin Scrub Packs include:*

- 1 CSR wrap
- 2 Cotton tipped applicators
- 1 Pair wrapped nitrile latex-free, powder-free gloves
- 2 Blue blotting towels
- 1 One-step PVP-I topical gel (1% available iodine, 3.7 oz. bottle)

Item Code	Description	Qty/Case
41560	2 large winged grip sponges	20
41561	4 small winged grip sponges	20

### Curity™ Sterile Wet Skin Scrub Components

Item Code	Description	Qty/Case
41550	2 sponge sticks saturated with PVP-I paint solution	20
41556	2 large winged grip sponges saturated with PVP-I scrub solution	20
41461	3.7 oz. bottles of One-Step Prep PVP-I Gel	48
41401	Iodine PVP Scrub 3.7 oz.	72
41411	Iodine PVP Paint 3 oz.	72
41561	Pack Gel Scrub Small Winged Sponge	20



## Patient Prep

### Curity™ Sterile Dry Skin Scrub Prep Products

#### Dry Skin Scrub Trays include:

- 1 CSR wrap
- 2 White absorbent towels
- 1 Ring forceps (5-inch long)
- 2 Cotton tipped applicators
- 1 Tray (2 compartments)
- 1 Tray (3 compartments)
- 1 Pair wrapped nitrile latex-free, powder-free gloves



#### Curity™ Sterile Dry Skin Scrub Trays

Item Code	Description	Qty/Case
41514	8 super sponges	20
41516	12 4"x 4" gauze sponges	20
41522	12 4"x 4" filled sponges	20
41523	3 sponge sticks, 6 winged grip sponges, 2 blue blotting towels	20
7402	10 super sponges	20
7405	6 super sponges, 6 radiopaque round sponges	20

### Curity™ Sterile Vaginal Prep, Wet Skin Scrub Products

#### Curity™ Sterile Vaginal Prep, Wet Skin Scrub Tray

Item Code	Description	Qty/Case
41536	1 CSR wrap	20
	2 cotton tipped applicators	
	6 dry contoured sponges	
	1 PVP-I scrub solution (3.7 oz. bottle)	
	1 PVP-I paint solution (3 oz. bottle)	
	2 large winged grip sponges	
	2 white absorbent towels	
	2 blue blotting towels	
	2 pairs wrapped latex-free, powder-free nitrile gloves	
	1 tray (2 compartments)	
	1 tray (3 compartments)	



#### Curity™ Sterile Vaginal Prep, Wet Skin Scrub Pack

Item Code	Description	Qty/Case
41533	1 CSR wrap	20
	2 cotton tipped applicators	
	2 large winged grip sponges saturated with PVP-I Scrub solution	
	2 sponge sticks saturated with PVP-I Scrub solution	
	2 sponge sticks saturated with PVP-I Paint solution	
	1 sponge stick	
	2 white absorbent towels	
	2 blue blotting towels	
	1 pair wrapped latex-free, powder-free nitrile gloves	

## Index

Code	Product	Page	Code	Product	Page
AWC Index			31142147	Devon™ Magnetic Needle Counters (with BoxLocks)	11
131610	Curity™ Disposable Scalpels #10, 10 per box	16	31142154	Devon™ Magnetic Needle Counters (with BoxLocks)	11
131611	Curity™ Disposable Scalpels #11, 10 per box	16	31142170	Devon™ Foam Strip Needle Counters (with BoxLocks)	10
131615	Curity™ Disposable Scalpels #15, 10 per box	16	31142188	Devon™ Foam Strip Needle Counters (with BoxLocks)	10
139090	Curity™ Disposable Scalpels #10, bulk packaged	16	31142196	Devon™ Magnetic Needle Counters (no BoxLocks)	11
139091	Curity™ Disposable Scalpels #11, bulk packaged	16	31142204	Devon™ Foam Strip Needle Counters (with BoxLocks)	10
139093	Curity™ Disposable Scalpels #15, bulk packaged	16	31142212	Devon™ Foam Strip Needle Counters (with BoxLocks)	10
166024	Curity™ Frazier Instruments	4	31142238	Devon™ Magnetic Needle Counters (with BoxLocks)	11
166025	Curity™ Frazier Instruments	4	31142295	Devon™ Foam Strip Needle Counters (with BoxLocks)	10
166026	Curity™ Frazier Instruments	4	31142303	Devon™ Foam Strip Needle Counters (with BoxLocks)	10
166030	Curity™ Sigmoid Instruments and Sets	4	31142311	Devon™ Foam Block Needle Counters (with BoxLocks)	10
166031	Curity™ Sigmoid Instruments and Sets	4	31142345	Devon™ Foam Block Needle Counters (with BoxLocks)	10
166035	Curity™ Sigmoid Instruments and Sets	4	31142352	Devon™ Foam Block Needle Counters (with BoxLocks)	10
31140117	Devon™ Light Handle (Rigid)	9	31142360	Devon™ Foam Block Needle Counters (with BoxLocks)	10
31140125	Devon™ Light Handle (Rigid)	9	31142378	Devon™ Foam Block Needle Counters (with BoxLocks)	10
31140133	Devon™ Light Handle (Rigid)	9	31142386	Devon™ Foam Block Needle Counters (with BoxLocks)	10
31140208	Devon™ Light Gloves	9	31142394	Devon™ Foam Block Needle Counters (with BoxLocks)	10
31140216	Devon™ Light Gloves	9	31142402	Devon™ Foam Block Needle Counters (with BoxLocks)	10
31140240	Devon™ Light Sleeve (Semi-Rigid)	9	31142410	Devon™ Magnetic Needle Counters (no BoxLocks)	11
31140257	Devon™ Light Gloves	9	31142428	Devon™ Foam Block Needle Counters (with BoxLocks)	10
31140273	Devon™ Light Sleeve (Semi-Rigid)	9	31142451	Devon™ Foam Strip Needle Counters (with BoxLocks)	10
31140281	Devon™ Disposable Drape Clamps Sterile drape clamp	14	31142469	Devon™ Foam Block Needle Counters (with BoxLocks)	10
31140299	Devon™ Disposable Drape Clamps Non-sterile drape clamp	14	31142485	Devon™ Foam Block Needle Counters (with BoxLocks)	10
31140307	Devon™ Disposable Drape Clamps Sterile drape clamp	14	31142493	Devon™ Foam Block Needle Counters (with BoxLocks)	10
31140315	Devon™ Disposable Drape Clamps Non-sterile tube occluding clamp	14	31142501	Devon™ Foam Block Needle Counters (with BoxLocks)	10
31140331	Devon™ Reusable Drape Clamps Non-sterile	14	31142527	Devon™ Fog-Out	14
31140356	Devon™ Magentic Drapes with Neutral Zones 3¾" x 9" x 1" tray	19	31142535	Devon™ Instrument & Cautery Tip Cleaners	14
31140372	Devon™ Magentic Drapes with Neutral Zones 4" x 8½" x 2"	19	31142543	Devon™ Scope Adjustment Bag	16
31140380	Devon™ Magentic Drapes with Neutral Zones 3¾" x 9" x 1" tray	19	31142550	Devon™ Instrument Pouch	16
31140398	Devon™ Magentic Drapes with Neutral Zones 3¾" x 9" x 1" tray	19	31142576	Devon™ Blade Shield Scalpel Holder	20
31140422	Devon™ Magnetic Needle Finder	21	31142659	Devon™ Tube holder, removable fasten strap	15
31140430	Devon™ Magentic Drapes with Neutral Zones 16" x 20"	19	31142667	Devon™ Tube holder, removable fasten strap	15
31140505	Devon™ Magentic Drapes with Neutral Zones 10" x 16"	19	31142675	Devon™ Tube holder, fixed fasten strap	15
31140539	Devon™ Magentic Drapes with Neutral Zones 16" x 20"	19	31142691	Devon™ Suture Bag	16
31140547	Devon™ Hands Free Transfer Magentic Drapes 16" x 20", with stays	19	31142717	Devon™ Instrument & Cautery Tip Cleaners	14
31140562	Devon™ Hands Free Transfer Magentic Drapes 16" x 20"	19	31142790	Devon™ Instrument & Cautery Tip Cleaners	14
31140570	Devon™ Hands Free Transfer Magentic Drapes 10" x 16"	19	31142808	Devon™ Fog-Out	14
31140588	Devon™ Reusable Magentic Drapes 10" x 16"	20	31144119	Sponge Counting Sytems Devon™ Bag-It dispenser box, clear back	18
31140596	Devon™ Hands Free Transfer Magentic Drapes 10" x 16", with stays	19	31144127	Sponge Counting Sytems Devon Bag-It dispenser box, opaque back	18
31140612	Devon™ Magentic Drapes with Neutral Zones 10" x 16"	19	31144168	Devon™ rack for I.V. pole, single	18
31142113	Devon™ Magnetic Needle Counters (with BoxLocks)	11	31144176	Devon™ rack for wall, single	18
31142139	Devon™ Magnetic Needle Counters (with BoxLocks)	11	31144184	Devon™ rack for I.V. pole, triple	18

## Index

Code	Product	Page	Code	Product	Page
31144192	Devon™ rack for wall, triple	18	41461	Curity™ Sterile Wet Skin Scrub Components	24
31144226	Devon™ Safe-T-Strainer Instrument Handling Products	21	41514	Curity™ Sterile Dry Skin Scrub Trays	25
31144234	Devon™ Safe-T-Strainer Instrument Handling Products	21	41516	Curity™ Sterile Dry Skin Scrub Trays	25
31144283	Devon™ Safe-T-Strainer Instrument Handling Products	21	41522	Curity™ Sterile Dry Skin Scrub Trays	25
31144291	Devon™ Safe-T-Strainer Instrument Handling Products	21	41523	Curity™ Sterile Dry Skin Scrub Trays	25
31145660	Devon™ Mini Surgical Vessel Loops Yellow	15	41529	Curity™ Sterile Wet Skin Scrub Packs with PVP-I Solutions	24
31145678	Devon™ Maxi Surgical Vessel Loops Blue	15	41530	Curity™ Sterile Wet Skin Scrub Packs with PVP-I Solutions	24
31145686	Devon™ Maxi Surgical Vessel Loops White	15	41531	Curity™ Sterile Wet Skin Scrub Packs with PVP-I Solutions	24
31145694	Devon™ Maxi Surgical Vessel Loops Yellow	15	41533	Curity™ Sterile Vaginal Prep, Wet Skin Scrub Pack	25
31145702	Devon™ Mini Surgical Vessel Loops Blue	15	41535	Curity™ Sterile Wet Skin Scrub Packs with PVP-I Solutions	24
31145710	Devon™ Mini Surgical Vessel Loops Red	15	41536	Curity™ Sterile Vaginal Prep, Wet Skin Scrub Tray	25
31145728	Devon™ Mini Surgical Vessel Loops White	15	41550	Curity™ Sterile Wet Skin Scrub Components	24
31145744	Devon™ Maxi Surgical Vessel Loops Red	15	41556	Curity™ Sterile Wet Skin Scrub Components	24
31145751	Devon™ Utility Labels	12	41560	Curity™ Sterile Wet Skin Scrub Packs with One Step Prep PVP-I Gel	24
31145769	Devon™ Utility Labels Bulk packaged	12	41401	Curity™ Sterile Wet Skin Scrub Components	24
31145777	Devon™ Surgical Markers Regular tip, flexible ruler, 9 labels	11	41411	Curity™ Sterile Wet Skin Scrub Components	24
31145785	Devon™ Surgical Markers Regular tip, flexible ruler	11	41561	Curity™ Sterile Wet Skin Scrub Components	24
31145793	Devon™ Surgical Markers Fine tip, flexible ruler, 9 labels	11	41561	Curity™ Sterile Wet Skin Scrub Packs with One Step Prep PVP-I Gel	24
31145819	Devon™ Surgical Markers Regular tip	11	41589	Curity™ Sterile Wet Skin Scrub Trays with PVP-I Solutions	23
31145827	Devon™ Surgical Markers Fine tip, flexible ruler	11	41590	Curity™ Sterile Wet Skin Scrub Trays with PVP-I Solutions	23
31145843	Devon™ Surgical Markers with Ruler Cap Fine tip	12	41591	Curity™ Sterile Wet Skin Scrub Trays with PVP-I Solutions	23
31145868	Devon™ Surgical Markers with Ruler Cap Dual tip, 9 labels	12	5096	Sponge Counting Systems	18
31145884	Sterile Devon™ Utility Markers	12	5097	Curity™ rack with I.V. pole mounting bracket	18
31145892	Devon™ Surgical Markers with Ruler Cap Fine tip, flexible ruler	12	5098	Curity™ rack with scale and I.V. pole mounting bracket	18
31145900	Devon™ Surgical Markers with Ruler Cap Dual tip	12	51000177	Devon™ Utility Markers Non-sterile	12
31145918	Devon™ Surgical Markers with Ruler Cap Regular tip, 9 labels	12	54901200	Kwik-Sorb™ Solidification Products 1,200 ml single dose bottle	7, 18
31145926	Devon™ Surgical Markers with Ruler Cap Regular tip	12	54901500	Kwik-Sorb™ Solidification Products 1,500 ml single dose bottle	7, 18
31145934	Devon™ Surgical Markers with Ruler Cap Dual tip	12	54902000	Kwik-Sorb™ Solidification Products 2,000 ml single dose bottle	7, 18
31145942	Devon™ Surgical Markers with Ruler Cap Fine tip	12	54903000	Kwik-Sorb™ Solidification Products 3,000 ml single dose bottle	7, 18
31145959	Devon™ Surgical Markers with Ruler Cap Regular tip	12	54916000	Kwik-Sorb™ Solidification Products 16,000 ml multi dose bottle	7, 18
31146130	Solidification Products 16,000 ml multi dose bottle	7, 18	54930000	Kwik-Sorb™ Solidification Products 3 lb. bucket	18
31146148	Solidification Products 2,000 ml single dose bottle	7, 18	6022	Curity™ Lap Sponges — Sterile Tray Pack 18" x 18"	13
31146163	Solidification Products 1,500 ml single dose bottle	7, 18	6034	Curity™ Lap Sponges — Sterile Tray Pack 12" x 12"	13
31160560	Devon™ Surgical Markers with Ruler Cap Dual tip, flexible ruler	12	6045	Curity™ Lap Sponges — Sterile Tray Pack 4" x 18"	13
31161006	Solidification Products 1,200 ml single dose bottle	7, 18	605039	Merlin™ One Step Prep Applicator 39 ml applicator PVP/IPA	23
31163473	Solidification Products 3,000 ml single dose bottle	7, 18	6056	Curity™ Lap Sponges — Sterile Tray Pack 8" x 36"	13
31177366	Devon™ Surgical Markers Fine tip	11	606139	Merlin™ One Step Prep Applicator 39 ml applicator alcohol-free	23
4031	Kerlix™ Lap Sponges — Sterile Tray Pack 17" x 26"	13	6522	Curity™ Lap Sponges — Sterile Soft Pouch Pack 18" x 18"	13
4034	Kerlix™ Lap Sponges — Sterile Tray Pack 16" x 17"	13	6534	Curity™ Lap Sponges — Sterile Soft Pouch Pack 12" x 12"	13
4037	Kerlix™ Lap Sponges — Sterile Tray Pack 8" x 32"	13	6545	Curity™ Lap Sponges — Sterile Soft Pouch Pack 4" x 18"	13
4048	Kerlix™ Lap Sponges — Sterile Tray Pack 3" x 17"	13	6556	Curity™ Lap Sponges — Sterile Soft Pouch Pack 8" x 36"	13
4069	Kerlix™ Lap Sponges — Sterile Tray Pack 8" x 11"	13	6822	Curity™ Lap Sponges — Non-Sterile Bulk Pack 18" x 18"	13

## Index

Code	Product	Page	Code	Product	Page
6834	Curity™ Lap Sponges — Non-Sterile Bulk Pack 12" x 12"	13	8888290619	Argyle™ Conductive Bubble Tubing	6
7402	Curity™ Sterile Dry Skin Scrub Trays	25	8888301507	Argyle™ Connecting Tubes with Molded Connectors	5
7405	Curity™ Sterile Dry Skin Scrub Trays	25	8888301515	Argyle™ Connecting Tubes with Molded Connectors	5
77701	Blue OR towel, 1 pack	13	8888301523	Argyle™ Suction Tubing with Molded Connectors	5
77702	Blue OR towel, 2 pack	13	8888301523	Argyle™ Connecting Tubes with Molded Connectors	5
77704	Blue OR towel, 4 pack	13	8888301531	Argyle™ Connecting Tubes with Molded Connectors	5
77706	Blue OR towel, 6 pack	13	8888301549	Argyle™ Connecting Tubes with Molded Connectors	5
77708	Blue OR towel, 8 pack	13	8888301606	Argyle™ Connecting Tubes with Molded Connectors	5
8888270207	Argyle™ Non-Conductive Bubble Tubing, Sterile	6	8888301614	Argyle Suction Tubing with Molded Connectors	5
8888270306	Argyle™ Non-Conductive Bubble Tubing, Clean - Bulk Packaged, Vinyl	6	8888301614	Argyle™ Connecting Tubes with Molded Connectors	5
8888270603	Argyle™ Non-Conductive Bubble Tubing, Sterile - Individually	6	8888301622	Argyle™ Connecting Tubes with Molded Connectors	5
8888270702	Argyle™ Non-Conductive Bubble Tubing, Clean - Bulk	6	8888301630	Argyle™ Connecting Tubes with Molded Connectors	5
8888271007	Argyle™ Non-Conductive Bubble Tubing, Sterile	6	8888301705	Argyle™ Suction Tubing with Molded Connectors	5
8888271106	Argyle™ Non-Conductive Bubble Tubing, Clean - Bulk Packaged, Vinyl	6	8888301705	Argyle™ Connecting Tubes with Molded Connectors	5
8888271403	Argyle™ Non-Conductive Bubble Tubing, Clean - Bulk	6	8888301713	Argyle™ Connecting Tubes with Molded Connectors	5
8888271411	Argyle™ Non-Conductive Bubble Tubing, Sterile - Individually	6	8888301721	Argyle™ Connecting Tubes with Molded Connectors	5
8888271502	Argyle™ Non-Conductive Bubble Tubing, Sterile - Individually	6	8888301739	Argyle™ Connecting Tubes with Molded Connectors	5
8888271510	Argyle™ Non-Conductive Bubble Tubing, Clean - Bulk	6	8888501007	Argyle™ Flexible Yankauer Suction Instruments	3
8888272005	Argyle™ Non-Conductive Bubble Tubing, Clean - Bulk	6	8888501015	Argyle™ Flexible Yankauer Suction Instruments	3
8888272013	Argyle™ Non-Conductive Bubble Tubing, Sterile - Individually	6	8888501023	Argyle™ Flexible Yankauer Suction Instruments	3
8888275008	Argyle™ Non-Conductive Bubble Tubing, Sterile - Individually	6	8888501031	Argyle™ Flexible Yankauer Suction Instruments	3
8888280214	Argyle™ Non-Conductive Bubble Tubing	6	8888501411	Argyle™ Flexible Yankauer Suction Instruments	3
8888280412	Argyle™ Non-Conductive Bubble Tubing	6	8888501429	Argyle™ Flexible Yankauer Suction Instruments	3
8888280610	Argyle™ Non-Conductive Bubble Tubing	6	8888501619	Argyle™ Flexible Yankauer Suction Instruments	3
8888284505	Argyle™ Connecting Tubes with Funnel/Funnel Integral Connectors	5	8888501635	Argyle™ Sump Tip Flexible Yankauer Suction Instruments	3
8888284513	Argyle™ Connecting Tubes with Funnel/Funnel Integral Connectors	5	8888501742	Argyle™ Flexible Yankauer Suction Instruments	3
8888284521	Argyle™ Connecting Tubes with Funnel/Tapered Integral Connectors	5	8888501775	Argyle™ Sump Tip Flexible Yankauer Suction Instruments	3
8888284539	Argyle™ Connecting Tubes with Funnel/Funnel Integral Connectors	5	8888502005	Argyle™ Flexible Yankauer Suction Instruments	3
8888284547	Argyle™ Connecting Tubes with Funnel/Funnel Integral Connectors	5	8888504001	Argyle™ Flexible Yankauer Suction Instruments	3
8888284554	Argyle™ Connecting Tubes with Funnel/Funnel Integral Connectors	5	8888505016	Argyle™ Yankauer Suction Instruments, Rigid Type	2
8888284604	Argyle™ Connecting Tubes with Funnel/Funnel Integral Connectors	5	8888505024	Argyle™ Yankauer Suction Instruments, Rigid Type	2
8888284612	Argyle™ Connecting Tubes with Funnel/Funnel Integral Connectors	5	8888505032	Argyle™ Rigid Yankauer Suction Instruments	2
8888284620	Argyle™ Connecting Tubes with Funnel/Funnel Integral Connectors	5	8888505040	Argyle™ Rigid Yankauer Suction Instruments	2
8888284638	Argyle™ Connecting Tubes with Funnel/Funnel Integral Connectors	5	8888505115	Argyle™ Yankauer Suction Instruments, Rigid Type	2
8888284703	Argyle™ Connecting Tubes with Funnel/Funnel Integral Connectors	5	8888505123	Argyle™ Yankauer Suction Instruments, Rigid Type	2
8888284711	Argyle™ Connecting Tubes with Funnel/Tapered Integral Connectors	5	8888505164	Argyle™ Yankauer Selec-Trol, Rigid Type	2
8888284737	Argyle™ Connecting Tubes with Funnel/Funnel Integral Connectors	5	8888505172	Argyle™ Yankauer Selec-Trol, Rigid Type	2
8888284745	Argyle™ Connecting Tubes with Funnel/Funnel Integral Connectors	5	8888509703	Argyle™ Poole and Suction Sets	4
8888290213	Argyle™ Conductive Bubble Tubing	6	8888509737	Argyle™ Poole and Suction Sets	4
8888290411	Argyle™ Conductive Bubble Tubing	6	8888509745	Argyle™ Poole and Suction Sets	4





COVIDIEN, COVIDIEN with logo and Covidien logo are US and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company. © Covidien 2015

H62000R xxM 05/15

15 HAMPSHIRE STREET  
MANSFIELD, MA  
02048

800-962-9888  
508-261-8000

[WWW.COVIDIEN.COM](http://WWW.COVIDIEN.COM)