

## Aluminum Transport Chair

HCPCS: E1038



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F

- Weighs only 19 lbs
- Aluminum frame is lightweight and strong
- Composite 8" wheels are lightweight and maintenance free
- Both 17" and 19" come with swing-away footrests
- Comes with seat belt for added safety
- Padded armrests provide additional comfort
- Nylon upholstery is comfortable, lightweight and easy-to-clean Back folds down for storage and transport
- Aluminum casted back-release hinge allows back to fold down for easy storage and transport

### PRODUCTS

### Description

ATC17-BK	Transport Chair 17" Black, 1/cs
ATC17-BL	Transport Chair 17" Blue 1/cs
ATC17-RD	Transport Chair, Alum. 17" Red, 1/ea
ATC19-BK	Transport Chair, Blk, 19" Black, 1/cs
ATC19-BL	Transport Chair, Alum. 19in Blue, 1/cs
ATC19-RD	Transport Chair, Alum. 19" Red, 1/cv

### SPECIFICATIONS

### Description

Carton Shipping Weight	24 lbs
Overall Height	39.5"
Overall Length w/ Riggings	33"
Overall Width	20" (open)   9" (closed)   22" (open)   9" (closed)
Seat Dimensions	17" (W) x 16" (D)   19" (W) x 16" (D)
Seat To Floor Height	19"
Warranty	Limited Lifetime
Weight	19 lbs
Weight Capacity	300 lbs
Width Between Arms	15"   17"

### ACCESSORIES

### Description

STDS1040S	Cup Holder, Swivel, 1 BX, 6/CS, BX
-----------	------------------------------------

### PARTS

### Description

STDS871	Hand Grip, Sentra, 1/ea, EA
STDS1026	Front Bumper for Ft. Rst/Leg Rest, 1/ea, EA
STDS1051	Dust Cap, fork stem, 1/ea, EA
STDS1011	Wheel Bearing, 8mm x 22mm 1/ea, EA

PARTS	Description
STDS8M4601	ForkStemBearing,SSP2,Cruiser,1/ea, EA
STDS001000	Axle Cap 8mm, 1ea, EA
STDS1002-1	Axle for 8" caster, 1/ea, EA
STDS3J3900	Front Caster,8"TR39,1/ea, EA
STDS5D1070	SeatRailGuide5D,1ea, EA
STDS8M4603	BrakeTipSleeve,Cruiser,SSP, 1/ea, EA
STDS5D1069	Hanger Pins for ATC, SSP, Cruiser,4/pk, BAG
STDS3J091R	Footrest Rlse Lever,ATC,Black,Rt,1/ea, EA
STDS3J091L	Footrest Rlse Lever,ATC,Lt.1/e, EA
STDS3J0917	Back Upholstery,17" ATC 3J, 1/ea, EA
STDS3J0919	Back Upholstery,19"ATC, 3J,5D,9H,1/ea, EA
STDS5D09K7	Push Handle,Black,ATC 5D,1ea, EA
STDS5D09B7	Push Handle,Blue,ATC 5D,1ea, EA
STDS5D09R7	Push Handle,Red,ATC 5D,1ea, EA
STDS5D09S7	Push Handle,Silver,ATC 5D,1ea, EA
STDS5D092L	Brake Assembly,Left,ATC 5D,1ea, EA