

Compact Compressor Nebulizer

HCPCS: E0570 A7003-1,E0570+A7003



FIGURE A



FIGURE B

- All nebulizers compressors come with a dishwasher safe reusable JetNeb Plus
- The MQ5900 also comes with a disposable neb kit and the MQ5900P with a pediatric aerosol mask
- Powerful Piston-Driven workhorse compressor
- Durable, lightweight, quiet compressor is perfect for a wide range of patients
- Compact & portable design with an integrated handle for transport

PRODUCTS

Description

MQ5800	CompactComprsrNeb,w/Handle,1/ea, 6/cv
MQ5900	ompressorNeb,DisposNeb&NebJet,1/ea,6/cv
MQ5900P	CmpctNeb,Disps&JetNeb,w/PedMask1/ea,6/c

SPECIFICATIONS

Description

Compressor Free Air Flow	8 lpm
Electrical Requirements	120 VAC 60Hz
Maximum Pressure	30 psig or Greater / 248 KPa
Operating Humidity Range	15%-90% Non-Condensing
Operating Noise Level	<63 dBA
Operating Pressure	10 psig or Greater
Operating Temperature Range	50-104°F, (10-40°C)
Power Consumption	110W Max
Storage/Transport Humidity Range	10%-90%
Storage/Transport Temperature Range	-13 to 158°F, (-25 to 70°C)
Unit Dimensions	7.5" (W) X 3.7" (H) X 5.7" (D)
Warranty	Limited 5 Year
Weight	3.0 lbs

ACCESSORIES

Description

18031	Nebulizer Carry Bag, 1/ea, EA
18254	Reusable Neb Kit 10/case, CS
NEB KIT 700	NebKit,IncPediMask,Tub ,50/cs, CS
MASK 001P	Pediatric Aerosol mask CS/50, CS
MASK 001A	Adult Aerosol mask CS/50, CS
NEB KIT 600	NebKit,Inc Adlt Mask,Tub 50/cs, CS
NEB KIT 500	Neb Kit,Inc Mouth pc&Tub50/cs, CS
18260	JetNeb Plus Breath Enhanced Neb,12/CS, CS
MQ92692	Filters for MQ5800, MQ5900 5/pkg, PKG