



**FLEXIBILITY.**

**VARIETY. RELIEF.**

**Medline Compression Systems**  
Choose your level of therapeutic compression

# More choices. More comfort.

## A CoFlex® TLC Two-Layer Compression System

Low-profile short stretch two-layer compression system with absorptive and odor controlling first layer, easy hand tear outer layer and educational packaging.

Item No.	Description	Compression	Pkg.
AND7800	CoFlex TLC—for therapeutic compression	30-40 mmHg	8 kts/cs
AND7800XL	CoFlex TLC XL—for larger and/or longer limbs	30-40 mmHg	8 kts/cs
AND7802	CoFlex TLC Lite—for lighter compression	20-30 mmHg	8 kts/cs

## B Fourflex® and Threeflex® Multilayer Compression Bandage Systems

Consistent, therapeutic compression with award-winning educational packaging.

Item No.	Description	Compression	Pkg.
MSC4400	Four-layer compression bandage system	30-40 mmHg	8 kts/cs
MSC4400XL	Extra long four-layer compression bandage system Bandages are 25-30% longer	30-40 mmHg	8 kts/cs
MSC4300	Three-layer, light-compression bandage system	20-30 mmHg	8 kts/cs

## C Unna-Z™ and Unna-Z Stretch Unna Boot Bandages with Zinc

Improved knit design and plastic core for easy, comfortable application. Unna-Z Stretch has additional elasticity for conformability and easier application.

Item No.	Description	Compression	Pkg.
NONUNNA13	3" x 10 yds. (7.6 cm x 9.1 m)	*NA	12/cs
NONUNNA14	4" x 10 yds. (10.2 cm x 9.1 m)	*NA	12/cs
NONUNNA3	3" x 10 yds. (7.6 cm x 9.1 m) with calamine	*NA	12/cs
NONUNNA4	4" x 10 yds. (10.2 cm x 9.1 m) with calamine	*NA	12/cs
NONUNNAS130	3" x 10 yds. stretched (7.6 cm x 9.1 m), elastic	*NA	12/cs
NONUNNAS140	4" x 10 yds. stretched (10.2 cm x 9.1 m), elastic	*NA	12/cs
NONUNNAS30	3" x 10 yds. stretched (7.6 cm x 9.1 m) with calamine, elastic	*NA	12/cs
NONUNNAS40	4" x 10 yds. stretched (10.2 cm x 9.1 m) with calamine, elastic	*NA	12/cs

\*Not Applicable - Unna Boots do not themselves provide active compression. To provide light to moderate compression, Unna Boots are traditionally covered with a cohesive or elastic bandage.

## D MEDIGRIP™ Elastic Tubular Bandage

Offers compression when doubled; may also be used for joint support in a single layer.

Item No.	Size	Width	Application	Format	Compression	Pkg.
MSC9500	A	1.75" (4.5 cm)	Very small feet and arms	11 yd. roll	**Varies	1/bx
MSC9501	B	2.5" (6.3 cm)	Small hands and limbs	11 yd. roll	**Varies	1/bx
MSC9502	C	2.625" (6.8 cm)	Adult hands, arms or legs	11 yd. roll	**Varies	1/bx
MSC9503	D	3" (7.5 cm)	Large arms or legs	11 yd. roll	**Varies	1/bx
MSC9504	E	3.5" (8.75 cm)	Legs or small thighs	11 yd. roll	**Varies	1/bx
MSC9505	F	4" (10 cm)	Large knees or thighs	11 yd. roll	**Varies	1/bx
MSC9506	G	4.75" (12 cm)	Large thighs	11 yd. roll	**Varies	1/bx
MSC9504YD	E	3.5" (8.75 cm)	Legs or small thighs	1 yd. pre-cut	**Varies	30/cs
MSC9505YD	F	4" wide (10 cm)	Large knees or thighs	1 yd. pre-cut	**Varies	30/cs
MSC9506YD	G	4.75" wide (12 cm)	Large thighs	1 yd. pre-cut	**Varies	30/cs
MSC9507	J	6.75" (17 cm)	Small trunks	11 yd. roll	**Varies	1/bx
MSC9508	K	8.25" (21 cm)	Medium trunks	11 yd. roll	**Varies	1/bx

\*\*Compression levels will depend on leg circumference and the size of Medigrip chosen.

### Unique Benefits of CoFlex TLC

Includes Stocking for Enhanced Patient Comfort

- » Less likely to catch or roll on sheets or clothing

Absorbent Coating

- » Absorbent material of the first layer manages moisture for the whole leg

Hand-Tear Cohesive Bandage

- » Faster and easier application of second layer

Patient Information Card

Odor Control



Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093

Medline United States  
1-800-MEDLINE (633-5463)  
medline.com | info@medline.com

Medline Canada  
1-800-396-6996  
medline.ca | canada@medline.com

Medline México  
01-800-831-0898  
medlinemexico.com | mexico@medline.com

FOLLOW US [f](#) [in](#) [t](#) [BLOG](#)