

Dasher 4



The Dasher 4 scooter is designed for the user on the go. It has a 250 lb. weight capacity and has a quick-disconnect system for effortless disassembly. The unit disassembles into three parts, the heaviest weighing just 26 lb. The Dasher 4 maneuvers with ease and has an adjustable tiller for added comfort during operation. The Dasher 4 also comes standard with solid tires for a smoother ride and swivel seating system.

- + 250lb Weight Capacity
- + Swivel Seating
- + Dismantles into Three Parts
- + Curved Tiller
- + Interchangeable Fender Colors Included



- Pewter
- Burgundy
- Blue

Model Number: GK8
Medicare Billing Code: K0800

HCPCS codes provided should not be considered as legal advice, and does not guarantee reimbursement.

Technical Specifications

Overall Dimensions (L x W x H)	in	39.5" x 22" x 35"
Number, Size of Tires	Front, Rear	2, 7" (2, 8")
Suggested User Weight ¹	lb	Max. 250
Weight of Heaviest Part ²	lb	33
Weight ³	lb	94
Each Battery Module Weight	lb	22
Battery Size	Type x PCS	12V12AH x 2
Charger	Amp	Off board (2A)
Output Power of Motor	HP	0.5
Maximum Speed	MPH	3.8
Maximum Climbing Angle/ Safe Climbing Angle ⁴	%(°)	10(6)
Range (Per Charge) ⁶	Mile	7.5
Turning Radius	in	37"
Ground Clearance	in	1.5"
Handle Bar	Type	T-Bar
Model #		GK8
Standard Seat	W x D x BH	SF181614B

Note: ¹ On level Road. ² Without basket and seating. ³ With basket, seat & batteries. ⁴ Base on user weight of 200LB. ⁵ Based on user weight of 165LB. ⁶ After the batteries and mechanical moving parts fully break in.