



# Comperm Sizing

## Comperm® LF and Comperm®-PD LF

Tubular bandages that maximize support and minimize constriction.

Part #	Description	Sizing Recommendation	D
83010000	Comperm LF Size A	Very Small Limbs or Feet	1
83020000	Comperm LF Size B	Small Limbs	2
83030000	Comperm LF Size C	Hands, Feet or Small Lower Legs	2
83040000	Comperm LF Size D	Lower Legs or Arms	3
83050000	Comperm LF Size E	Knees or Small Thighs	3
83060000	Comperm LF Size F	Large Knees or Thighs	4
83070000	Comperm LF Size G	Large Thighs	5
83080000	Comperm LF Size J	Small Trunks	7
83090000	Comperm LF Size K	Medium Trunks	8
83100000	Comperm LF Size L	Large Trunks	1
83120000	Comperm - PD LF Small	Large Arms and Small Legs	3
83130000	Comperm - PD LF Large	Large Knees and Calves	4

The above guidelines are suggested usages. Each patient will vary in size and shape and their optimal size will vary. Comperm bandages (Comperm - PD) are for use with patients who need or would like extra cushioning to protect the limb.

### Comperm Compression Guidelines using mmHg:

For the best fit, measure the widest circumference of the extremity.

	Low Compression	Medium Compression	High Compression
<b>Circumference</b>	<b>5 - 10 mmHg</b>	<b>10 - 20 mmHg</b>	<b>20 - 30 mmHg</b>
6" - 10" (15 - 25cm)	Sizes D or C	Sizes C or B	Sizes B or A



### Comperm LF Compression Bandage size chart.pdf

Added yesterday at 1:32 PM • 290.07 KB

➤ Open in new tab   ✕ Delete