

**ASTHMA MANAGEMENT**

**ASTHMAMD™ PEAK FLOW METER**

**Description:**

The AsthmaMD™ Peak Flow Meter is a portable, pocket-sized asthma management tool to help patients measure their lung performance and anticipate attacks.

- Measures peak expiratory flow with universal 60-810 LPM range
- Accurate: <+/- 10% or 10 l/m
- Meets all ATS standards for portable peak flow meters <sup>1</sup>
- Designed for children and adults
- Single-patient Use
- Durable, lightweight, portable, cleanable, and easy to use
- Does not contain natural rubber latex or DEHP <sup>2</sup>
- Colorful, integrated markers allow patients to easily monitor changes in air flow with large visible scale
- Ergonomically shaped mouthpiece to assist in creating a secure seal
- Compatible with a physician-developed AsthmaMD™ Peak Flow Meter Mobile app\*



United States

To order this product, please call:

**800-397-5899**

[Click for Prices](#)

**VitalityMedical.com**  
Add Vitality to Your Life

Catalog #	Description	Qty/Case
1802	AsthmaMD™ Peak Flow Meter	10

References: 1. 510(k): K101380 2. Data on File, Noble International, Inc.

\*The AsthmaMD™ Peak Flow Meter app is available as a free download for iPhone/iPad and Android users from a third party and has incorporated video tutorials.

DISCLAIMER: The AsthmaMD™ mobile app is a third-party app offered by AsthmaMD, Inc. Teleflex has no control over such third-party app, nor does Teleflex review, approve, monitor, endorse, or make any representations or warranties with respect to such third-party app. Phone also not included.

AsthmaMD is a trademark of Genia Medical, Inc.

© 2017 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available

in all regions; contact customer service to confirm availability in your region.